



#### Non-animal method for predicting skin sensitization (OECD TG442C)

# **ADRA kit**

ADRA (**A**mino **a**cid **D**erivative Reactivity Assay) is a *in chemico* test method, representing an alternative (nonanimal) method for the evaluation of skin sensitization compounds. The kit comprises two amino acid derivates, NAC (N-(2-(1- **n**aphthyl)**a**cetyl)-L-**c**ysteine) and NAL ( $\alpha$ -N-(2-(1-**n**aphthyl)**a**cetyl)-L-Iysine), which are composed out of a nucleophilic region and a detection molecule (naphthalene rings).

ADRA is proposed to address the molecular initiating event of skin sensitization AOP (Adverse Outcome Pathway) by quantifying the reactivity of the test compounds against the model synthetic amino acid derivatives. The reaction is monitored by HPLC-analysis after 24 hours incubation at 25°C and 281 nm, determining the relative residual concentration of NAC and NAL in the reaction liquid.



## Advantages of ADRA

- NAC and NAL are detected at a relatively long wavelength (281 nm).
  - Preventing co-elution of the test chemical and the nucleophilic reagent.
- ADRA avoids precipitation by using only 1% of reactants of the existing method.
- Variability of result values are very low.
  - •••• In ADRA, the reaction is stopped by the addition of a fixing solution (2.5%TFA) before analyzing Multiple test chemicals can be assayed in a short period of time.
    - ··· The test procedure is performed using a single 96-well plate and a multichannel pipette.

#### Kit Reagents

	Components	Amount
	NAC	2 vials $\times$ 10 ml
	NAL	2 vials $\times$ 10 ml
	NAC Buffer (pH 8.0) premixed	2 vials $\times$ 300 ml
	NAL Buffer (pH 10.2) premixed	2 vials $\times$ 300 ml
	0.01 mol/I EDTA Solution	2 vials $\times$ 1 ml

Code No.	Product Name	Package Size
296-80901	ADRA Kit	1 Kit

Our company sells licensed products by Fujifilm Corporation.

### Additional reagents required

Code No.	Product Name	Package Size
204-02743	Trifluoroacetic Acid	25 ml
015-08633	Acetonitrile	31
217-01031	Ultrapure Water *	11
016-00346	Acetone **	500 ml
043-07216	Dimethyl Sulfoxide **	500 ml

\* Use water with low metal content.

\*\* Not to be used, if the test chemical dissolves in water or acetonitrile.

Phenylacetaldehyde (CAS RN<sup>®</sup> 122-78-1) is required as a positive control.

#### Apparatus required

- 1) Electronic balance – precision scale:  $\pm 0.1$  mg
- 2) Three Micropipettes – 10 µl, 100 µl and 1000 µl
- 3) 12 channel pipette - 150 µl
- 4) HPLC system -light-shielding auto-sampler for 96-well plates and 0.3 ml/min liquid feeding
- 5) UV detector – photodiode array (PDA) detector or absorbance detector (281 nm)
- 6) HPLC column
- 7) pH meter - precision scale: ±0.01 pH
- Incubator 25 °C 8)
- 96-well plates 9)
- 10) 500 mL plastic bottles
- 11) Vortex mixer
- Plate seal\* 12)
  - \* Use the seal having high sealability and solvent-resistant performance.
- Plate shaker 13)
- 14) Plate centrifuge

Note: to avoid NAC/NAL dimerization upon metal ion contamination, all consumables (except for components for HPLC analysis) must be made of polypropylene or polyethylene.

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: https://labchem-wako.fujifilm.com/ / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

# FUJIFILM Wako Laboratory Chemicals site https://labchem-wako.fujifilm.com



FUJIFILM Wako Pure Chemical Corporation 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Tel: +81 6 6203 3741 Fax: +81 6 6203 1999 ffwk-cservise@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation FUJIFILM Wako Chemicals Europe GmbH 1600 Bellwood Road, Richmond, VA 23237, U.S.A. Toll-Free (U.S. only): +1 877 714 1920 Tel: +1 804 271 7677 Fax: +1 804 271 7791 wkuslabchem@fujifilm.com

Fuggerstr 12, 41468 Neuss, Germany Tel: +49 2131 311 0 Fax: +49 2131 311 100 labchem\_wkeu@fujifilm.com