



UV Curing Accelerator

DPNG, IPEMA

UV curing technology uses less solvent and cures faster than heat-based radical polymerization, so it is highly productive, energy-saving, and eco-friendly. A universal problem in UV curing technology is reaction inhibition due to oxygen in the air, and this effect is particularly pronounced for thin films with a large surface area per volume. DPNG (Diprenyl Glycerin Ether) and IPEMA (Isoprenyl Methacrylate) are additives that reduce oxygen inhibition and accelerate curing through different mechanisms in UV curing under air. Please use it for research and development of polymer synthesis.

DPNG (<u>Diprenyl Glycerin Ether</u>)

IPEMA (<u>I</u>so<u>pre</u>nyl <u>M</u>eth<u>a</u>crylate)


Feature

Useful as a UV curing accelerator

•DPNG: possible to absorb and decompose oxygen by itself

• IPEMA: crosslinking agent with two olefins with different reactivity

<UV curing in the presence of DPNG and IPEMA>




Ex.	DPNG	IPEMA	Initiator	Conversion (after 30 sec.)	
1	(O ₂ blocked)		1 wt%	74.5%	
2	1 wt%	_	1 wt%	46.2%*	
3	_	1 wt%	1 wt%	44.5%	
4	_	_	5 wt%	41.5%	
5	_	_	3 wt%	20.2%	
6	_	_	1 wt%	4.3%	

Monomer: Pentaerythritol triacrylate. Initiator: 1-hydroxycyclohexyl phenyl ketone. UV irradiation intensity: 73 mW/cm². UV source: high pressure mercury lamp.

^{*} Calculated from the average value of 30 ± 0.3 seconds.

<Antioxidant function of DPNG in polyurethane resin>



(left) blank, (light) addition of DPNG

<Curling suppression of resin by addition of IPEMA>





(left) blank, (light) addition of IPEMA

*These products are developed by Kuraray Co., Ltd. The photo was provided by the company.

Product List

Code No.	Product Name	Structure	CAS RN® Storage conditions	Package Size
021-19741	1,3-Bis[(3-methyl-2-buten-1-yl)oxy]-2-	ООН	2337348-25-9	100mL
023-19745	propanol 【DPNG】		room temperature	500mL
132-19301	3-Methyl-3-buten-1-yl Methacrylate		156291-88-2	100mL
134-19305	【IPEMA】		room temperature	500mL

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: https://labchem-wako.fujifilm.com/ / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site https://labchem-wako.fujifilm.com



FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Tel: +81 6 6203 3741 Fax: +81 6 6203 1999 ffwk-cservise@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road, Richmond, VA 23237, U.S.A. Toll-Free (U.S. only): +1 877 714 1920 Tel: +1 804 271 7677 Fax: +1 804 271 7791 wkuslabchem@fujifilm.com

Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong Tel: +852-2799-9019 Fax: +852-2799-9808 wkhk.info@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH Fuggerstr 12, 41468 Neuss, Germany Tel: +49 2131 311 0 Fax: +49 2131 311 100

labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited FUJIFILM Wako (Guangzhou) Trading Corporation

Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road, Guangzhou, 510095, China Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai) Tel: +86-10-6413-6388(Beijing) wkgz.info@fujifilm.com