(109 × 210mm Size)

**FUJIFILM**

**Wako**

Code No. 293-79801

## Orexin A ELISA Kit *Wako*

### [Introduction]

Orexin is a neuropeptide that controls eating and sleep/wakefulness. Orexin occurs as two types, orexin A and orexin B, produced by cleavage of the same precursor, and they act by binding to orexin receptors. There are two types of orexin receptors, OX1R and OX2R, and it has been reported that OX1R exhibits about 50-time higher affinity to orexin A compared to orexin B and OX2R exhibits equal affinity to orexin A and B.

Orexin is considered to have a function in regulating wakefulness and sleep because orexin A was decreased remarkably in the cerebrospinal fluid of patients with narcolepsy, a sleep disorder.

This product is an ELISA kit for simple measurement of orexin A in human cerebrospinal fluid as well as rat cerebrospinal fluid, serum and plasma. The kit does not cross-react with orexin B.

### [Kit performance]

Range of the calibration curve	4.69 to 300 pg/mL
Measurement samples	Human cerebrospinal fluid, rat cerebrospinal fluid, rat plasma, rat serum
Necessary amount of sample	25 $\mu$ L
Duration of measurement	About 20 hours
Within-run reproducibility	CV < 5%
Between-run reproducibility	CV < 16%

### [Components of the kit]

1) Antibody-coated 96-well plate	1 plate
2) Orexin A Standard	300 pg $\times$ 1 vial
3) Anti Orexin A, Biotin-conjugated	12 mL $\times$ 1 bottle
4) Peroxidase-conjugated Streptavidin Solution	12 mL $\times$ 1 bottle
5) TMB Solution	12 mL $\times$ 1 bottle
6) Buffer	20 mL $\times$ 1 bottle
7) Stop Solution	12 mL $\times$ 1 bottle
8) Wash Solution (20 $\times$ )	50 mL $\times$ 1 bottle
9) Plate Seal	4 sheets

### [Measurement principle]

In each well of the measurement plate (96 wells), rabbit polyclonal orexin A antibody is immobilized. A standard solution or sample is placed in this well to bind orexin A with the antibody. Then, biotin-conjugated anti orexin A mouse monoclonal antibody and peroxidase-conjugated streptavidin are allowed to react. Finally, the peroxidase activity in the well is measured to determine the concentration of orexin A in the sample.

**[Apparatuses and equipment to be used]**

- Micropipets and tips (25 to 1000  $\mu\text{L}$ ) (8- or 12-channel multichannel pipet is recommended)
- Microplate reader (equipment that can measure the absorbance up to 3.0 at 450 nm)
- Microplate agitator or shaker
- Glass test tubes
- Microplate washer (continuous dispenser, needle dispenser, and aspirator or vacuum pump are recommended for manual method)
- Volumetric cylinder (1000 mL)
- Distilled water or deionized water

**[Methods for preparation of reagents]**

- Standard solutions**

Add 1 mL of the Buffer to the container of the reference standard to dissolve the content and prepare 300 pg/mL standard solution. Pipet 0.2 mL of this standard solution, and dilute this with 0.2 mL of the Buffer to prepare 150 pg/mL standard solution. Repeat the same dilution procedure and prepare 75.0, 37.5, 18.8, 9.38, and 4.69 pg/mL standard solutions. Use the Buffer as it is as 0 pg/mL standard solution.
- Washing solution**


Dilute 50 mL (whole volume) of the concentrated washing solution with 950 mL of distilled water and use.
- Use other reagents as they are.**

**[Measurement procedure]**

1. Equilibrate the components of the kit at room temperature (20 to 30  $^{\circ}\text{C}$ ).
2. Dispense 350  $\mu\text{L}$  of the washing solution to each well, discard the solution by sucking with the aspirator or turning over the plate, and remove the solution completely by lightly patting against paper towel. Repeat the procedure two more times to wash the plate 3 times in total.
3. Dispense 100  $\mu\text{L}$  of the Buffer to each well and, then, add 25  $\mu\text{L}$  of the standard solution or sample. Perform the procedure from the start of dispensing the standard solution to the end of dispensing the sample as quickly as possible (within 30 minutes).
4. Seal the measurement plate with a plate seal and agitate at room temperature for 18-24 hours (at about 100 rpm).
5. Remove the solution from each well, and perform the washing procedure similarly to 2. three times in total.
6. Add 100  $\mu\text{L}$  of Anti Orexin A, Biotin-conjugated to each well.
7. Seal the measurement plate with a Plate Seal and agitate at room temperature for 1 hour (at about 100 rpm).
8. Separate a necessary amount of TMB Solution about 1 hour before use and return to room temperature under the condition protected from light.

9. Remove the solution from each well, and perform the washing procedure similarly to 2. three times in total.
10. Add 100  $\mu\text{L}$  of Peroxidase-conjugated Streptavidin Solution to each well.
11. Seal the measurement plate with a Plate Seal and agitate at room temperature for 30 minutes (at about 100 rpm).
12. Remove the solution from each well, and perform the washing procedure similarly to 2. four times in total.
13. Add 100  $\mu\text{L}$  of TMB Solution to each well, leave to stand under the condition protected from light, and allow to react at room temperature for 30 minutes.
14. Add 100  $\mu\text{L}$  of the Stop Solution to each well.
15. Determine the absorbance at 450 nm using an absorption spectrometer for microplates.
16. Prepare a standard curve from the measured values of the absorbance at the respective concentrations of orexin A standard solutions using a commercially available software and 5 (or 4)-parameter regression equation and determine the concentration of orexin A in the sample. If double logarithmic plotting paper is used, plot the concentration of the standard solution along the horizontal axis and the absorbance of the standard solution at each concentration along the vertical axis to prepare the standard curve, apply the absorbance of the sample to the standard curve, and read the concentration of orexin A.

**Example of calibration curve**

**[Precautions for use]**

1. For collecting cerebrospinal fluid samples, it is strongly advised to add Triton X-100 to prevent a decrease in the measured values due to adhesion of orexin A in the sample to the container. For addition, prepare phosphate-buffered saline solution containing 0.3% Triton X-100 and mix in the sample in the ratio of 1 to 100 as follows. As an alternative, mixing equal amounts of cerebrospinal fluid and Buffer give the same effect.

#### Examples of addition

Amount of sample	Amount of 0.3% Triton X-100/PBS added
1 mL or less	10 $\mu$ L
1-2 mL	20 $\mu$ L
2-3 mL	30 $\mu$ L

2. After collecting a blood sample, separate serum or plasma and measure promptly. After collecting a cerebrospinal fluid sample, remove cellular blood and cell components by centrifugation and perform the measurement promptly. If not measured promptly, subdivide appropriately and store in a freeze at  $-30^{\circ}\text{C}$  or lower. Collect plasma in blood sampling tubes containing EDTA-2Na (1 mg/mL). Do not repeat freezing and thawing samples.
3. Prepare the reagents just before use as a rule. In particular, use the reference standard promptly after preparation. If the kit is used in divided times, use a glass container (such as a vial) for storage of the reference standard after dissolution. In divided uses, dilute newly from the 300 pg/mL standard solution. The reference standard (300 pg/mL) after dissolution is stable for 1 week of refrigerated storage at  $4^{\circ}\text{C}$  and for about 4 weeks of frozen storage at  $-30^{\circ}\text{C}$  or lower.
4. Perform the procedure from the start to the end of dispensing of the standard solution as quickly as possible (within 30 minutes).
5. In the case of a sample at a concentration higher than 300 pg/mL, dilute the sample with the Buffer attached to this kit.
6. For dilution of the standard solution, make sure to use a new tip at each dilution step. Use a new tip for each sample when pipetting into the well.
7. Make sure to perform the measurement in duplicate.
8. Agitate the measurement plate by the microplate agitator to allow reaction (except for coloring reaction). Agitate slowly taking care not to splash the reaction solution on the Plate Seal (at about 100 rpm).
9. Do not use the kits of different lots in combination.
10. Be sure to prepare a standard curve for each measurement as the color development level of the enzyme substrate may be affected slightly by the reaction temperature or time, or degree of agitation of the measurement plate.
11. After stopping the reaction, determine the absorbance promptly.
12. The concentrated washing solution may form precipitation during its storage but this precipitation dissolves by dilution with no quality issue.
13. Take care not to expose each reagent to intense light during storage or use.
14. Use TMB Solution after returning to room temperature under the condition protected from light.

**[Measurement examples]**

**■ Recovery test**

	No addition	Orexin A (5 pg/mL)		Orexin A (20 pg/mL)		Orexin A (100 pg/mL)	
	Measured value (pg/mL)	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery
Rat Cerebrospinal fluid [1]	51.82	54.96	96.7%	63.74	88.8%	120.77	79.6%
Rat Cerebrospinal fluid [2]	47.69	55.55	105.4%	63.40	93.7%	109.95	74.5%
Rat Cerebrospinal fluid [3]	38.05	43.99	102.2%	59.12	101.8%	101.16	73.3%
Human Cerebrospinal fluid [1]	12.25	14.29	82.8%	23.25	72.1%	65.16	58.1%
Human Cerebrospinal fluid [2]	47.21	50.45	96.6%	55.71	82.9%	107.12	72.8%
Human Cerebrospinal fluid [3]	24.41	25.99	88.4%	35.07	79.0%	83.51	67.1%
Rat Serum [1]	Not detected	5.64	112.8%	19.94	99.7%	89.58	89.6%
Rat Serum [2]	6.42	12.26	107.4%	25.05	94.8%	93.60	88.0%
Rat Plasma [1]	1.45	6.81	105.6%	20.30	94.6%	91.13	89.8%
Rat Plasma [2]	2.28	7.20	98.9%	21.64	97.1%	90.51	88.5%

**■ Dilution test**

	x 1	x 2		x 4		x 8	
	Measured value (pg/mL)	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery
Rat Cerebrospinal fluid [1]	47.47	19.77	83.3%	9.47	79.8%	4.46	75.2%
Rat Cerebrospinal fluid [2]	126.06	60.37	95.8%	30.21	95.9%	14.92	94.7%
Rat Cerebrospinal fluid [3]	116.04	56.39	97.2%	27.18	93.7%	12.99	89.6%
Rat Cerebrospinal fluid [4]	173.06	94.49	109.2%	49.82	115.2%	26.72	123.5%
Human Cerebrospinal fluid [1]	21.48	8.71	81.1%	4.65	86.6%	2.71	100.9%
Human Cerebrospinal fluid [2]	18.37	7.28	79.3%	3.61	78.6%	2.44	106.3%
Human Cerebrospinal fluid [3]	16.50	6.36	77.1%	3.08	74.7%	1.53	74.2%

	x 1	x 1.5		x 2		x 3	
	Measured value (pg/mL)	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery	Measured value (pg/mL)	Recovery
Human Cerebrospinal fluid [4]	58.59	39.63	101.5%	29.84	101.9%	20.12	103.0%
Human Cerebrospinal fluid [5]	26.14	20.87	119.8%	15.76	120.6%	11.86	136.1%

**[Outline of procedure]**

- Fully return plates and reagents to room temperature (20 to 30 °C).
- Dilution of concentrated washing solution : Dilute 20 times with distilled water.
- Preparation of standard solutions: Dissolve Orexin A Standard thoroughly by adding 1 mL of the Buffer to the vial to prepare 300 pg/mL standard solution. Then, dilute the solution with the Buffer as follows to prepare 150, 75, 37.5, 18.8, 9.38, and 4.69 pg/mL standard solutions.

Concentration (pg/mL)	300	150	75	37.5	18.8	9.38	4.69	0
Standard solutions (mL)	0.2*	0.2*	0.2*	0.2*	0.2*	0.2*	0.2*	0.2*
Buffer (mL)	0.2	0.2	0.2	0.2	0.2	0.2	0.2	0.2

\* : Standard solution of one-level higher concentration

- Antibody-coated 96-well plate
  - ↓ Wash three times.
- Buffer, 100  $\mu$ L/well
- Sample or standard solution, 25  $\mu$ L/well
  - ↓ Cover the plate with a Plate Seal and agitate at room temperature for 18-24 hours (at about 100 rpm).
  - ↓ Wash three times.
- Anti Orexin A, Biotin-conjugated, 100  $\mu$ L/well
  - ↓ Cover the plate with a Plate Seal and agitate at room temperature for 1 hour (at about 100 rpm).
  - ↓ Separate TMB Solution, and return to room temperature under protection from light (1 hour before use).
  - ↓ Wash three times.
- Peroxidase-conjugated Streptavidin Solution, 100  $\mu$ L/well
  - ↓ Cover the plate with a Plate Seal and agitate at room temperature for 30 minutes (at about 100 rpm).
  - ↓ Wash four times.
- TMB Solution, 100  $\mu$ L/well
  - ↓ Leave to stand for 30 minutes at room temperature under protection from light.
- Stop Solution, 100  $\mu$ L/well
  - ↓ Determine the absorbance at 450 nm by the plate reader.

**[Storage condition]**

Store at 2 to 10 °C.

**[Expiration date]**

Indicated on the label.

**[Packaging]**

96 times

---

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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