



EV-Up™ EV Production Basal Medium for MSC, AF MSC EV Production Supplement, AF

- 1 More EVs cultured in EV-Up™ can be collected than in serum media.
- 2 EV-Up™ produces EVs with high activity.
- 3 Cell viability can be maintained during EV production.

EV-Up™ EV Production Basal Medium for MSC, AF and EV-Up™ MSC EV Production Supplement, AF are medium and supplement for the effective exosomes (EVs) production from mesenchymal stem cells (MSCs). EV-Up™ as a set composed of the medium and the supplement can be used as a complete medium. These products are serum-free and animal component-free, and applicable to various growth media.

Products information Medium and its supplement are intended to be used as a set.


Code No.	Product name	Grade	Pkg. size
053-09451	EV-Up™ EV Production Basal Medium for MSC, AF 	For cell culture	95mL
298-84001	EV-Up™ MSC EV Production Supplement, AF 		For 100mL

Applicable cells

The complete medium is applicable for MSCs derived from various tissue sources.

bone marrow, umbilical cord, adipose tissue etc.

Procedure




*1 Mixture of EV-Up™ EV Production Basal Medium for MSC, AF and EV-Up™ MSC EV Production Supplement, AF.

The collected EVs can be isolated by the PS affinity method*2 using MagCapture™ Exosome Isolation Kit PS (Code No. 293-77601).

*2 EVs are captured specifically by phosphatidylserine (PS)-binding proteins in presence of metal ions. The captured EVs can be released afterwards with high purity by adding chelating agents such as EDTA.


1. Cell Viability

After the expansion of human bone marrow-derived MSCs in serum containing media, the medium was transferred to EV-Up™ medium and cultured for five days to produce EVs. MSCs cultured in EV-Up™ produced EVs without affecting the MSC viability, comparable high survival rate to conventional DMEM + 10% EV depleted FBS was obtained.




2. Number of EV particles

EVs isolated from various media supernatant by the PS affinity method were analyzed with NTA. MSCs cultured in EV-Up™ medium released 2.6 times more EVs than MSCs cultured in DMEM + 10% EV-depleted FBS. The particle diameter of the EVs were almost the same.



3. Anti-fibrotic Effect

5x10⁷ particles/mL of EVs isolated from various media supernatant by the PS affinity method were added to normal human fetal lung-diploid fibroblasts cells (TIG3) that were stimulated with TGFβ. And the fibrotic marker (αSMA) gene expression was quantified by RT-PCR. Significantly, MSC EVs produced in EV-Up™ media decreased the gene expression of fibrotic markers such as αSMA.



Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site
<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
 Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
 ffwk-cservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
 1600 Bellwood Road, Richmond, VA 23237, U.S.A.
 Toll-Free (U.S. only): +1 877 714 1920
 Tel: +1 804 271 7677 Fax: +1 804 271 7791
 wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
 Fuggerstr 12, 41468 Neuss, Germany
 Tel: +49 2131 311 0 Fax: +49 2131 311 100
 labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited
 Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road,
 Tsuen Wan, N.T., Hong Kong
 Tel: +852-2799-9019 Fax: +852-2799-9808
 wkhk.info@fujifilm.com

FUJIFILM Wako (Guangzhou) Trading Corporation
 Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road,
 Guangzhou, 510095, China
 Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai)
 Tel: +86-10-6413-6388(Beijing)
 wkgz.info@fujifilm.com