For Genetic Research microRNA Extractor[®] Kit for Purified EV

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This kit is for RNA extraction, including micro RNA, from purified Extracellular Vesicle (EV) derived from cell culture supernatants, serum, and plasma. microRNAs can be easily and efficiently obtained from EVs by using a proteolytic enzyme, a chaotropic agent, and spin columns. There is no need to use toxic chemicals such as phenol or chloroform.

This product is intended for laboratory use only.

【Features】

- Optimal for RNA extraction from purified EV solution
- Excellent RNA recovery
- Phenol and chloroform are not required
- Extracted RNA can be used for real-time quantitative PCR, microarrays, and next-generation sequencing

【Kit contents (for 20 extractions)】

This kit consists of 9 components.

| | |
|-------------------|-----------------|
| Lysis Solution | 9 mL × 1 tube |
| Reducing Agent | 80 μL × 1 tube |
| Protease Solution | 200 μL × 1 tube |
| Enhancer | 200 μL × 1 tube |
| Wash Solution I | 5 mL × 1 tube |
| Wash Solution II | 5 mL × 2 tubes |
| Elution Buffer | 1 mL × 2 tubes |
| Spin Column | 20 tubes |
| Collection Tube | 20 tubes |

【Other materials required】

Reagent :

- (1) 2-Propanol (Code No. 168-21675)
- (2) Ethanol (99.5) (Code No. 054-07225)
- (3) 1-Butanol (Code No. 022-16035)

Equipment :

- (1) 1.5 mL Micro centrifugation tube
- (2) Micropipette
- (3) Pipette tip
- (4) Cooling micro high-speed centrifuge
- (5) Vortex mixer

【Precautions】

1. Tools

Microcentrifugation tubes, pipette tips, and other experimental tools should be autoclaved, or use commercially available

RNase-free products. Plastic gloves and masks should be worn during the experiments to avoid contamination of RNase.

2. Reagents

Since the procedure mainly extracts RNA, handles the reagents with great care, and keeps the working reagents as sterile as possible.

3. Handling of samples and disposable tools used in the experiment

Handle the purified EV solution as infectious samples. Waste fluid generated during the RNA extraction protocol, pipette tips, microcentrifuge tubes, and other disposable tools used in the experiments should also be disposed of appropriately as infectious waste in accordance with the guidelines of your institution.

【Protocol】

Preparation of purified EV samples and precautions

Prepare 200 μL of purified EV solution in advance using the MagCapture™ Exosome Isolation Kit PS Ver.2 (Code No. 290-84103) or by some other means^{*1}. This solution is used as the sample for extraction.

- *1 If the volume is less than 200 μL, make up the sample to 200 μL with the solvent used to prepare the purified EV solution.

Preparation of reagents

● Wash Solution I (+2-Pro)

Add 7.5 mL of 2-propanol to Wash Solution I and mix well.

● Wash Solution II (+Ethanol)

Add 15 mL of ethanol to Wash Solution II and mix well.

● Lysis Solution (+ reducing agent)

Immediately before use, add 8.0 μL of the reducing agent to 1.0 mL of the Lysis Solution^{*2} and mix with a vortex mixer^{*3}.

● Spin column^{*4}

Place the spin column into a collection tube in advance.

- *2 The components of the Lysis Solution may precipitate. If this occurs, warm the solution in a water bath of ~ 37°C, and make sure that the precipitate is completely dissolved before use.

- *3 The prepared solution can be stored refrigerated (2-10°C) for about one week. Since the required amount of Lysis Solution changes depending on the number of samples for extraction, calculate the required amount in advance, transfer the solution to microcentrifugation tubes or other containers, and then add the reducing agent.

- *4 Spin Column may have deposits of filling material and scratches on then, but there is no problem with its performance of it.

Protocol of microRNA extraction

- (1) Transfer 10 μL of Protease Solution to a sterile 1.5 mL mi-

crocentrifugation tube.

- (2) Add 200 μL of the prepared **sample for extraction** to the microcentrifugation tube prepared in step (1).
- (3) Add 200 μL of the prepared **Lysis Solution (+ reducing agent)**^{*5}, mix with a vortex mixer, and spin down using a tabletop centrifuge.
- (4) Incubate the solution in a block incubator at 37°C for 30 minutes and spin down using a tabletop centrifuge.
- (5) Add 200 μL of 1-butanol^{*6} and mix with a vortex mixer.
- (6) Add 10 μL of Enhancer and mix with a vortex mixer.
- (7) Add 300 μL of 2-propanol and mix with a vortex mixer.
- (8) Spin down the solution using a tabletop centrifuge for about 2-3 seconds^{*7}.
- (9) Mix the solution by pipetting^{*8} and load half of the solution on the spin column.
- (10) Centrifuge the spin column at 8,000 $\times g$ for 1 minute at room temperature, remove the spin column, and discard the solution (filtrate) in the collection tube.
- (11) Mix the remaining solution by pipetting^{*8} and load the solution on the spin column.
- (12) Repeat step (10).
- (13) Load 500 μL of the **Wash Solution I (+2-Pro)** on the spin column.
- (14) Repeat step (10).
- (15) Load 700 μL of the **Wash Solution II (+Ethanol)** on the spin column.
- (16) Repeat step (10).
- (17) Repeat steps (15) and (16).
- (18) Set the spin column on a sterile 1.5 mL tube^{*9} and centrifuge the column at 8,000 $\times g$ for 3 minutes at room temperature to remove the washing solution altogether.
- (19) Set the spin column on a new sterile 1.5 mL tube.
- (20) Carefully load 20 μL of Elution Buffer^{*10} on the center of the membrane of the spin column and let it stand for 5 minutes at room temperature.
- (21) Centrifuge the column and tube at 8,000 $\times g$ for 1 minute at room temperature; the solution eluted into the 1.5 mL tube can be used as the RNA sample^{*11}.

*5 Do not add the **Lysis Solution (+ reducing agent)** directly to the proteolytic enzyme. The order of addition must be as follows : **Protease Solution** \rightarrow **sample for extraction** \rightarrow **Lysis Solution (+ reducing agent)**.

*6 If 1-butanol is not available, 2-propanol can be used instead. However, if 2-propanol is used, the recovered RNA fraction will be contaminated by more protein; therefore, 1-butanol should be used whenever possible

*7 Do not spin down for a prolonged period; removing the solution at the top of the microcentrifugation tube is sufficient.


*8 As microRNA may precipitate, be sure to mix the solution by pipetting.

*9 The collection tube can also be used continuously, but a new sterile 1.5 mL microcentrifuge tube is recommended to avoid contaminating the RNA sample with alcohol.

*10 The recovery rate may increase slightly when more Elution Buffer is loaded, but the RNA concentration in the eluate will decrease. The amount of Elution Buffer to be loaded should be between 6 and 100 μL , depending on the application after extraction.

*11 Store the collected RNA sample at -80°C until use in subsequent analysis.

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