## LabAssay™ Phospholipid

### 1. Introduction

Phospholipids are a class of lipids, and a major component of all biological membranes, along with glycolipids, cholesterol and proteins. LabAssay™ Phospholipid is based on an enzymatic methods using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of phospholipids in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

### 2. Kit Contents

	Components	Use State	Amount
(1)	Buffer Solution	liquid	50 mL/3 bottles
(2)	Chromogen Substrate	freeze dry	For 50 mL/3 bottles
(3)	Standard Solution	liquid	5 mL/1 bottle




※ Assay in a microplate ; 500 tests. Assay in a test tube ; 50 tests.

### 3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

### 4. Assay Principle

Phospholipids (lecithin, sphingomyelin and lysolecithin) in a sample are hydrolyzed to choline in a reaction catalyzed by phospholipase D. Choline so formed is oxidized by choline oxidase in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalyzed by peroxidase (POD), producing a blue pigment. The amount of phospholipids contained in the sample is determined by measuring the absorbance of the blue color.

### 5. Materials and Apparatuses to be Prepared

- 96wells microplate (transparent type) · Micropipette · Incubator maintaining at 37°C \*
- Plate mixer\*
- Microplate reader with 600 nm wavelength filter (\* if the microplate reader is not equipped.)

#### (For Test Tube method)

- Test tube · Pipette · Incubator maintaining at 37°C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

### 6. Preparation of Reagents to be used

#### ① Chromogen Substrate :

Dissolve 1 vial of Chromogen Substrate with 50 mL of Buffer Solution. After reconstitution, the solution should be stored at 2°C ~ 10°C and used within one week.

#### ② Standard Solution (Microplate method)

Standard Solution is prepared by dilution of the provided Choline Chloride Standard.

No.	Standard Solution	Distilled or deionized water	Sample volume	Phospholipid Concentration
1	50 $\mu\text{L}$	150 $\mu\text{L}$	2 $\mu\text{L}$	75.0 mg/dL
2	100 $\mu\text{L}$	100 $\mu\text{L}$	2 $\mu\text{L}$	150 mg/dL
3	Undiluted	–	2 $\mu\text{L}$	300 mg/dL
4	Undiluted	–	4 $\mu\text{L}$	596.1 mg/dL <sup>*1</sup>

\*1 The test sample volume is usually 2  $\mu\text{L}$ , but 4  $\mu\text{L}$  is taken in this case. The phospholipids concentration must be corrected accordingly as indicated in the table above.

## 7. Preparation of Samples

- Samples should be used immediately after collecting.
- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.
- This method measures lecithin, sphingomyelin, lysolecithin but not cephalin.

## 8. Procedure

### (1) Assay in a Microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 $\mu\text{L}$	300 $\mu\text{L}$	300 $\mu\text{L}$
Sample	Serum 2 $\mu\text{L}$	Standard solution 2 $\mu\text{L}$	–
Mix well and incubate at 37°C for 5 min. Measure the absorbance at 600 nm <sup>*2</sup> of the test sample and standard solution with the blank solution as the control.			

\*2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

### (2) Assay in a Test Tube


Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	Serum 20 $\mu\text{L}$	Standard solution 20 $\mu\text{L}$	–
Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm. Spectrophotometer : 600 nm <sup>*2</sup>			

\*2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

## 9. Standard Curve

[assay in a microplate]

## 10. Performance

### (1) Sensitivity [assay in a Test Tube]

- The absorbance is below 0.14, when measuring purified water as a sample.
- The absorbance is 0.12 to 0.58, when measuring 300 mg/dL phospholipid as a sample.

### (2) Specificity

The phospholipid concentration is less than  $\pm 15\%$ , when measuring the known concentration of a control serum as a sample.

## 11. Usage Notes

### (1) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The reaction completes in about 3 min. It is enough to incubate at 37°C for 5 min. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

### (2) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.

### (3) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

## [Reference]

Takayama, M., Itoh, S., Nagasaki, T. and Tanimizu, I. : *Clin. Chim. Acta.*, 79, 93 (1977).

## LabAssay™ Phospholipid

[Storage] Store at 2°C - 10°C (Do not freeze)

[Expiration date] Indicated on the label.

[Package] For 500 tests

[Cat #] 295-94401

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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