**FUJIFILM****Wako**

Code No. 052-09301 (50 μL)

## Exosomes, from COLO201 cells, purified

This product is highly purified exosomes from COLO201 cells (ATCC® CCL-224™, American Type Culture Collection), which is derived from human colon cancer tissue, by a PS affinity method. The exosome surface markers, CD9, CD63 and CD81 are expressed on these exosomes.

They can be used as positive control on qualitative analysis with Western blotting, ELISA and Flow cytometry. And it is possible to correct the variations in measured values between assays since the measured absorbance of the sample can be converted to the concentration by making the standard curve using this product as a calibrator, when measuring CD63 signals on exosomes using the anti-CD63 antibody, which is a control detection antibody in PS Capture™ Exosome ELISA Kits (Code No.297-79201 and 298-80601).

### [Storage]

Store at -20°C

Store in working aliquots is recommended. For example, store 10×5 μL aliquots at -20°C and thaw it at room temperature before use. Ten freeze-thaw cycles at -20°C or -80°C did not affect the quality.

### [Concentration]

10 μg/mL

This concentration is a corrected protein concentration value by a CD63 signal so that variations among lots of the quality can be reduced when it is used as a positive control and a calibrator.

### [Precaution for use]

This product is derived from COLO201 cells (ATCC® CCL-224™, American Type Culture Collection). The cells were negative for human immunodeficiency virus, hepatitis B virus, human papillomavirus, Epstein-Barr virus and cytomegalovirus test and other adventitious agent were not tested, so handling it as a potentially infectious sample is recommended.

### [Procedure]

#### ⟨To be used as positive control of exosomes of human CD9, CD63 and CD81⟩

Protocol

Dilute this product appropriate concentration and use as positive control of exosomes for Western Blotting, ELISA and Flow cytometry on qualitative analysis.

#### ⟨To be used as a calibrator for PS Capture™ Exosome ELISA Kits⟩

This product offers to correct the variations in measured values between assays since the measured absorbance of the sample can be converted to the concentration by making the standard curve using this product as a calibrator, when measuring CD63 signals on exosomes using the anti-CD63 antibody, which is a control detection antibody in PS Capture™ Exosome ELISA Kits (Code No.297-79201 and 298-80601).

Protocol

Dilute this product to appropriate concentrations with Reaction Buffer in order to make its absorbance value at 450 nm (absorbance at 620-650nm is subtracted) for ELISA measurement to be within the optical density range of 0.1 to 3.0 O.D.

—1/4—

### -Usage example-

To make a standard curve of PS Capture Exosome ELISA Kit (Streptavidin HRP), produce two-fold dilution series of 7 steps with a starting concentration of 10 ng/mL.

- (1) Add 2 μL of this product to 198 μL of Reaction Buffer, which is provided in the kit, to prepare 100-fold dilution (100 ng/mL).
- (2) According to the following table, prepare subsequent 7-steps of two-fold serial dilutions with Reaction Buffer from 10 ng/mL to 0.15625 ng/mL.

Final Conc. (ng/mL)	Starting Conc. (ng/mL)	Volume to dilute (μL)	Diluent (μL)
10.0	100	50	450
5.0	10.0	250	250
2.50	5.0	250	250
1.25	2.50	250	250
0.625	1.25	250	250
0.3125	0.625	250	250
0.15625	0.3125	250	250
0.00 (Blank)	—		250

### [Related products]

Code No.	Description	Size
293-77601	MagCapture™ Exosome Isolation Kit PS	10 purifications
290-80301	PS Capture™ Exosome Isolation Resin Kit	1 kit (0.5 mL Slurry)
297-79201	PS Capture™ Exosome ELISA Kit (Anti Mouse IgG POD)	96 reactions
297-79701	PS Capture™ Exosome Flow Cytometry Kit	300 reactions
016-27061	Anti CD63, Monoclonal Antibody (3-13)	20 μL
014-27643	Anti CD63, Monoclonal Antibody (3-13), Fluorescein Conjugated	For 100 tests
017-27753	Anti CD63, Monoclonal Antibody (3-13), Red Fluorochrome (635) Conjugated	For 100 tests
019-27713	Anti CD63, Monoclonal Antibody (3-13), Biotin Conjugated	100 μL
014-27763	Anti CD9, Monoclonal Antibody (1K)	100 μL
011-27773	Anti CD81, Monoclonal Antibody (17B1)	100 μL
290-35591	Magnet Stand	1 each
295-71701	microRNA Extractor™ SP Kit	50 reactions

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshimachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : +81-6-6203-3741  
Facsimile : +81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : +1-804-271-7677  
Facsimile : +1-804-271-7791  
<http://www.wakousa.com>

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