For Genetic Research

**PS Capture™ Exosome ELISA Kit
(Anti Mouse IgG POD)****【1. Introduction】**

Extracellular vesicles such as exosome have attracted attention as a messenger of cell-to-cell communication and a biomarker of diseases since they include proteins, mRNAs, microRNAs, and DNAs on their surface or inside and they are stably present in body fluids such as blood, urine, saliva, spinal fluid, and breast milk after being secreted from cells.

The kit includes reagents for enzyme-linked immunosorbent assay (ELISA) available for a qualitative analysis of extracellular vesicles purified from cell culture supernatant or samples from body fluid as well as a quantitative analysis of extracellular vesicles in samples of cell culture supernatant. It can detect extracellular vesicles, which have any surface marker protein, with high sensitivity by using a mouse monoclonal antibody against any surface marker protein of extracellular vesicles as a primary detection antibody and HRP-conjugated anti mouse IgG antibody of the kit as a secondary detection antibody after extracellular vesicles are captured by a plate on which proteins that specifically bind with phosphatidylserine (PS) on the surface of extracellular vesicles are immobilized. As a control primary detection antibody, anti human CD63 mouse monoclonal antibody is included in the kit. By using this, human CD63 positive extracellular vesicles can be detected.

This kit can easily detect surface marker proteins of extracellular vesicles purified from a sample with MagCapture™ Exosome Isolation Kit PS (Code No. 293-77601) with 50 to 1,000 times higher sensitivity than Western blot. Also, extracellular vesicles in samples of cell culture supernatant can be quantitatively measured by using extracellular vesicles purified from cell culture supernatant with MagCapture™ Exosome Isolation Kit PS as a reference standard.


Fig. 1 Measurement principle

【2. Summary of the measuring method】

Diluted cell culture supernatant samples or samples of extracellular vesicles purified by MagCapture™ Exosome Isolation Kit PS and diluted are incubated with stirring at room temperature for 2 hours in the wells of the Exosome Capture 96 Well Plate (PS binding protein immobilized microplate). After washing, add mouse monoclonal antibody against any surface marker protein of extracellular vesicles as a primary detection antibody or Control Primary Antibody Anti CD63 included with the kit and incubate that with stirring at room temperature for 1 hour. After washing, add Secondary Antibody HRP-conjugated Anti mouse IgG and incubate that with stirring at room temperature for 1 hour. After washing again, it is reacted with TMB Solution (coloring reagent) at room temperature for 30 minutes. Add Stop Solution and measure absorbance at 450 nm (reference wavelength 620 nm). Compare measurement values of each sample.

When the kit is used for a quantitative measurement, prepare a standard curve by using extracellular vesicles purified from samples with MagCapture™ Exosome Isolation Kit PS as a reference standard (Protein concentration or particle numbers should be measured beforehand) and plotting absorbance values of the dilution series of the reference standard, to determine the concentration of the sample from the standard curve.

【3. Intended purpose】

- (1) A qualitative analysis of extracellular vesicles purified from cell culture supernatant or samples from body fluids

The kit enables us to conduct a highly sensitive qualitative analysis of any surface marker protein of extracellular vesicles purified from cell culture supernatant or samples from body fluids with MagCapture™ Exosome Isolation Kit PS by using a mouse monoclonal antibody against any surface marker protein of extracellular vesicles as a primary detection antibody.

- (2) A quantitative analysis of extracellular vesicles in cell culture supernatant
Extracellular vesicles positive for any marker protein in cell culture supernatant can be quantitatively analyzed by using extracellular vesicles purified from cell culture supernatant with MagCapture™ Exosome Isolation Kit PS as a reference standard and using a mouse monoclonal antibody against any surface marker protein of extracellular vesicles as a primary detection antibody.

Note : Although Control Primary Antibody Anti CD63 as a control primary detection antibody in the kit can detect human CD63, it cannot detect mouse, rat, and bovine CD63. When a surface marker protein other than human CD63 is required to be detected, use an appropriate mouse monoclonal antibody.

Note : Since Secondary Antibody HRP-conjugated Anti mouse IgG as a secondary detection antibody of the kit can strongly react non-specifically with mouse IgG in a sample and weakly react non-specifically with human IgG and rat IgG, a quantitative analysis of serum or plasma samples including these IgGs should be avoided.

【4. Reagent supplied】

Components	State	Amount
(A) Exosome Capture 96 Well Plate	Use after washing	8 well × 12 strips/1 plate
(B) Plate Seal	—	4 sheets
(C) Reaction/Washing Buffer (10 ×)	Concentrated Use after pretreatment*	50 mL × 2 vials
(D) Exosome Binding Enhancer (100 ×)	Concentrated Use after pretreatment*	10 mL × 1 vial
(E) Control Primary Antibody Anti CD63 (100 ×)	Concentrated Use after pretreatment*	120 μL × 1 vial
(F) Secondary Antibody HRP-conjugated Anti mouse IgG (100 ×)	Concentrated Use after pretreatment*	120 μL × 1 vial
(G) TMB Solution	Ready for use	12 mL × 1 vial
(H) Stop Solution	Ready for use	12 mL × 1 vial
(I) Instruction Manual	—	1 copy

Note : each reagent should be prepared in accordance with section 10. Preparation of reagents.

【5. Storage and the expiration date of the kit】

Store the kit at 2-10 °C (Freezing is strictly prohibited). The kit is stable until the expiration date under the condition (The expiration date is printed on the label). Do not use expired reagents. Please use opened reagents as soon as possible since they may be affected by storage conditions.

【6. The storage method of each reagent when the kit is separately used】**(A) Exosome Capture 96 Well Plate**

When the plate strip is separately used, the remained unused strips should be returned to a zip seal pack and stored at 2-10 °C. The strips are stable until the expiration date.

(C) Reaction/Washing Buffer (10 ×)

After it is returned to a room temperature and required volume of the solution is dispensed, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

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After it is returned to a room temperature and required volume of the solution is dispensed, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

(E) Control Primary Antibody Anti CD63 (100 ×)

After required volume of the solution is dispensed soon after taking out from a refrigerator, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

(F) Secondary Antibody HRP-conjugated Anti mouse IgG (100 ×)

After required volume of the solution is dispensed soon after taking out from a refrigerator, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

(G) TMB Solution

After required volume of the solution is dispensed soon after taking out from a refrigerator, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

(H) Stop Solution

After required volume of the solution is dispensed soon after taking out from a refrigerator, the remained solution should be stored at 2-10 °C with their caps tightened. The solution is stable until the expiration date.

【7. Reagents required other than the kit】

- Use as a check box
- Purified water (distilled water)
- Primary detection antibody against any surface marker protein of extracellular vesicles (Except antibody for detection of human CD63)
~Recommended products~
Anti CD9 Antibody (HI9a) (Novus Biologicals, LLC : NB100-77915)
Anti CD81 Antibody (M38) (Novus Biologicals, LLC : NBP1-44861)
- MagCapture™ Exosome Isolation Kit PS (Code No. 293-77601) (When a reference standard is prepared)

【8. Equipment and materials required other than the kit】

- Use as a check box
- Test tubes for dilution of sample and antibody
- Glassware for preparation of Reaction/Washing Buffer (Graduated cylinder)
- Micropipette
- Syringe-type repeating dispenser such as Eppendorf multipipette or ThermoFisher Finnpipette (If available)
- Disposable reservoir (As necessary)
- Paper towel to remove the remained solution in wells.
- A vortex-type mixer
- A shaker for 96 well-plate (500 rpm)
Recommended products : Code No. 623-05671 MS 3 digital
Recommended products : Code No. 627-05691 MS 3 basic
- An automatic washer for 96 well-plate (If available)
Recommended products : Code No. 510-22411 HydroFlex™ M8Ch2
- A 96 well-plate reader (Required wavelength is 450 nm ± 10 nm and 600-650 nm)
Recommended products : Code No. 518-84231 Infinite® F50R [450 nm, 620 nm]
- Software for data analysis

【9. Preparation of a reference standard of an extracellular vesicle (When a quantitative measurement is conducted)】

Purify extracellular vesicles from a sample with MagCapture™ Exosome Isolation Kit PS in accordance with the instruction manual. In regard to the purified extracellular vesicles, protein concentration shall be measured with such as BCA method or particle numbers shall be measured by Nanoparticle Tracking Analysis (by using such as NanoSight LM-10)^{※1}. Purified extracellular vesicles should be refrigerated^{※2}.

※1 Conduct each measurement in accordance with an appropriate protocol.

※2 Since responsiveness of extracellular vesicles tends to be weakened by freezing and thawing, refrigeration storage is recommended.

【10. Preparation of reagents】

Reagents described as “Ready for use” in 【4. Reagent supplied】 can be used as it is after it is returned to a room temperature. Reagents described as “Concentrated and use after pretreatment” shall be prepared as follows. Only required amount of each reagent shall be prepared at the time of use.

Preparation of the Reaction/Washing Buffer (1 ×)

After Reaction/Washing Buffer (10 ×)^{※3} is diluted 10-fold with purified water (distilled water), which has been returned to a room temperature, add 1/100 volume of Exosome Binding Enhancer (100 ×) to the diluted solution^{※4} and mix them well^{※5}.

- ※3 Since ingredients of Reaction/Washing Buffer (10 ×) may be deposited by refrigeration storage, after returning that to a room temperature, check no presence of deposition before use.
- ※4 Be sure to add Exosome Binding Enhancer (100 ×) since it is an essential ingredient for extracellular vesicles to bind with a plate.
- ※5 Since ingredients of the Reaction/Washing Buffer (1 ×) tend to be deposited after addition of Exosome Binding Enhancer (100 ×), please use it within 8 hours after preparation.

(e.g.) When solution for 96-well reactions is prepared
 After adding 900 mL of purified water (distilled water) to 100 mL of Reaction/Washing Buffer (10 ×) to make 1000 mL, add 10 mL of Exosome Binding Enhancer (100 ×) and mix them well.

Preparation of the reaction solution for anti CD63 antibody as a control detection antibody (When human CD63 is measured as a surface antigen)

Add 1/100 volume of Control Primary Antibody Anti CD63 (100 ×) to Reaction/Washing Buffer (1 ×) and mix them well.

(e.g.) When solution for 96-well reactions is prepared
 Add 100 μL of Control Primary Antibody Anti CD63 (100 ×) to 10 mL of Reaction/Washing Buffer (1 ×) and mix them well.

Preparation of the reaction solution for a primary detection antibody (When a surface antigen other than human CD63 is measured)

Dilute mouse monoclonal antibody against any surface protein of extracellular vesicles with Reaction/Washing Buffer (1 ×)^{※6} to an appropriate concentration. A rough standard of the antibody concentration in a reaction solution for a primary detection antibody is 200-500 ng/mL.

(e.g.) When the reaction solution for a primary detection antibody (250 ng/mL) is prepared for 96-well reactions from the primary detection antibody (1 mg/mL)
 To prepare 40-fold dilution, add 195 μL of Reaction/Washing Buffer (1 ×) to 5 μL of the primary detection antibody (1 mg/mL) and mix them. Then add 100 μL of the 40-fold dilution to 9.9 mL of Reaction/Washing Buffer (1 ×) and mix them well.

Preparation of the reaction solution for a secondary detection antibody

Add 1/100 volume of Secondary Antibody HRP-conjugated Anti mouse IgG (100 ×) to Reaction/Washing Buffer (1 ×)^{※6} and mix well.

(e.g.) When it is used for all 96 wells
 Add 100 μL of Secondary Antibody HRP-conjugated Anti-mouse IgG (100 ×) to 10 mL of Reaction/Washing Buffer (1 ×) and mix them well.

Preparation of the dilution series of a reference standard of extracellular vesicles for preparation of a standard curve (When a quantitative measurement is conducted)

Dilute a reference standard of extracellular vesicles prepared in 【9. Preparation of a reference standard of an extracellular vesicle】 to an appropriate concentration with Reaction / Washing Buffer (1 ×)^{※6} in order to make its absorbance value at 450 nm (absorbance at 620 nm is subtracted) for ELISA measurement to be within the range of 0.1-3.0.

(e.g.) Regarding a reference standard of extracellular vesicles (Protein concentration is 20 μg/mL) : when 7 steps of 2-fold serial dilution are prepared from 20 ng/mL of a reference standard of extracellular vesicles

- 1) Prepare a 10-fold dilution (2000 ng/mL) of a reference standard of extracellular vesicles by mixing 5 μL of the reference standard of extracellular vesicles (20 μg/mL) with 45 μL of Reaction/Washing Buffer (1 ×).
- 2) Mix a reference standard of extracellular vesicles with Reaction/Washing Buffer (1 ×) in accordance with the following table to prepare 7 steps of 2-fold serial dilution from 20 ng/mL.

Concentration of standard solution (ng/mL)	Volume of standard solution	Reaction/Washing Buffer (1 ×)
20.0	10-Fold Dilution (2000 ng/mL) : 5 μ L	495 μ L
10.0	Standard solution (20.0 ng/mL) : 250 μ L	250 μ L
5.00	Standard solution (10.0 ng/mL) : 250 μ L	250 μ L
2.50	Standard solution (5.0 ng/mL) : 250 μ L	250 μ L
1.25	Standard solution (2.5 ng/mL) : 250 μ L	250 μ L
0.625	Standard solution (1.25 ng/mL) : 250 μ L	250 μ L
0.313	Standard solution (0.625 ng/mL) : 250 μ L	250 μ L
0.00 (Blank)	—	250 μ L

【11. Dilution of samples】

Dilute samples with Reaction/Washing Buffer (1 ×)^{※6} using a tube, etc. so that absorbance at 450 nm for ELISA (value obtained by subtracting the absorbance at a secondary wavelength of 620 nm) falls in the range between 0.2 to 2.5. Also, when measuring undiluted samples, please add 1/100 volume of Exosome Binding Enhancer (100 ×) to them.

※6 Dilute a primary detection antibody, a secondary detection antibody, a reference standard of extracellular vesicles, and samples with Reaction/Washing Buffer (1 ×) containing Exosome Binding Enhancer prepared in 【6. The storage method of each reagent when the kit is separately used】. In other cases, there is also the possibility that using a TBS buffer containing BSA (Final concentration : 1 w/v %) and Exosome Binding Enhancer (Final concentration : 1 ×) may decrease background and increase responsiveness. Since ingredients of PBS-based buffer are deposited after addition of Exosome Binding Enhancer (100 ×), please do NOT use it.

【12. Operating method】

In advance of following each washing step, prepare the reagents to be dispensed next in accordance with 【6. The storage method of each reagent when the kit is separately used】.

- 1) Wash each well of Exosome Capture 96 Well Plate (strips to be used) with 300-350 μ L of Reaction/Washing Buffer (1 ×) 3 times^{※7}. Then, wipe off the remained fluid in the wells by turning the plate upside down and gently tapping it on a folded paper towel.
- 2) Dispense 100 μ L of sample dilution, standard dilution (When quantitative analysis is conducted), or Reaction/Washing Buffer (1 ×) as a blank in each well.
- 3) Label a plate seal^{※8} and react them by stirring at about 500 rpm using a microplate shaker at room temperature (20-25 °C) for 2 hours^{※9}.
- 4) After reaction, remove the reaction solution and wash each well with 300-350 μ L of Reaction/Washing Buffer (1 ×) 3 times. Then, wipe off the remained fluid in the wells by turning the plate upside down and gently tapping it on a folded paper towel.
- 5) Dispense 100 μ L of either anti CD63 antibody reaction solution or any primary detection antibody reaction solution in each well.
- 6) Label a plate seal^{※8} and react them by stirring at about 500 rpm using a microplate shaker at room temperature (20-25 °C) for 1 hour^{※9}.
- 7) After reaction, remove the reaction solution and wash each well with 300-350 μ L of Reaction/Washing Buffer (1 ×) 3 times. Then, wipe off the remained fluid in the wells by turning the plate upside down and gently tapping it on a folded paper towel.
- 8) Dispense 100 μ L of the secondary detection antibody reaction solution in each well.
- 9) Label a plate seal^{※8} and react them by stirring at about 500 rpm using a microplate shaker at room temperature (20 to 25 °C) for 1 hour^{※9}.

- 10) After reaction, remove the reaction solution and wash each well with 300 to 350 μL of Reaction/Washing Buffer ($1 \times$) 5 times. Then, wipe off remained fluid in the wells by turning the plate upside down and gently tapping it on a folded paper towel.
- 11) Dispense 100 μL of TMB Solution, which is returned to a room temperature, in each well and stir them with a microplate shaker for about 1 minute.
- 12) Label a plate seal^{※8} and react them by standing at room temperature (20-25 $^{\circ}\text{C}$) for 30 minutes.
- 13) Add 100 μL of Stop Solution, which is returned to a room temperature, in each well.
- 14) After the plate is stirred using a microplate shaker for about 5 seconds, measure the absorbance at 450 nm and 620 nm^{※10} with a 96-well microplate reader.

※7 Deposits may be observed on surface of a plate. It does not affect quality of a product.

※8 When a plate strip is separately used, cut a plate seal into a size of the strip for use.

※9 Since detection sensitivity becomes lower under the condition of standing and the error between wells tends to become larger, reaction by stirring at about 500 rpm using a rotary microplate shaker is recommended.

※10 When the absorbance at 620 nm cannot be measured, measure it at a wavelength in the range of 620-650 nm.

【13. Calculation】

In the following calculation, use the value of 450 nm by subtracting the absorbance value of 620 nm (600-650 nm) from the absorbance value of 450 nm.

< Qualitative measurement >

Calculate the value by subtracting the absorbance value of blank (450 nm) from the absorbance value of diluted samples (450 nm) and compare the values between samples.

< Quantitative measurement >

- 1) Calculate the value by subtracting the absorbance value of blank (450 nm) from the absorbance value of a diluted sample (450 nm)^{※11}. ... An absorbance value of a diluted sample
- 2) Calculate the value by subtracting the absorbance value of blank (450 nm) from the absorbance values of the dilution series of the reference standard (450 nm). ... An absorbance value of a reference standard
- 3) Plot the concentration of the reference standard on the horizontal axis and its absorbance value on the vertical axis to prepare a standard curve^{※12}.
- 4) Read a concentration value corresponding to an absorbance value of a diluted sample with the standard curve^{※13}. Multiply the dilution ratio of the sample and the read concentration together to determine the measurement.

※11 When an absorbance value of a diluted sample is deviated from the range of a standard absorbance curve, measure it again after diluting it with an appropriate scale.

※12 Prepare a standard curve for each measurement.

※13 Usage of a third-degree polynomial with 4 or 5 parameters is recommended for computer software processing.


【14. Reference data on performance of the kit】

■1 Qualitative analysis of extracellular vesicles purified from various cell culture supernatants

Add 1 ng of extracellular vesicles purified from a various cell culture supernatant to a well and surface marker levels were compared by a qualitative analysis with the kit using a primary detection antibody, which recognized cell surface markers, CD9, CD63, or CD81 each. Also, as reference comparative data, 150 ng of extracellular vesicles purified from a various cell

culture supernatant was electrophoresed and protein levels of each surface markers were detected by Western blot using a primary detection antibody, which recognized cell surface markers, CD9, CD63, or CD81 each. Then, qualitative analysis was conducted.

<Qualitative analysis of extracellular vesicles>

<Reference comparative data>

Purified Exosome : 150 ng


Fig. 2 Qualitative analysis of extracellular vesicles purified from various cell culture supernatants

■2 Qualitative analysis of extracellular vesicles purified from healthy human serum

40, 20, or 10 ng of extracellular vesicles purified from 6 healthy human serum samples were added to each well and qualitative analysis was conducted using a control primary detection antibody which recognized the cell surface marker CD63.

<Qualitative analysis of extracellular vesicles>



Fig. 3 Qualitative analysis of extracellular vesicles purified from healthy human serum

■3 Reference data : Standard curve of extracellular vesicles purified from cell culture supernatant


An example of a standard curve of extracellular vesicles purified from cell culture supernatant of COLO 201 cells (Detected by anti CD63 antibody)

Concentration (ng/mL)	Abs.(450 nm - 620 nm)			
	1	2	Mean	∠Blank
0	0.051	0.0522	0.0516	-
0.3125	0.0859	0.0882	0.0871	0.0355
0.625	0.132	0.1339	0.1330	0.0814
1.25	0.2198	0.2217	0.2208	0.1692
2.5	0.3812	0.382	0.3816	0.3300
5	0.7137	0.7307	0.7222	0.6706
10	1.4068	1.406	1.4064	1.3548
20	2.6552	2.6261	2.6407	2.5891

An example of a standard curve of extracellular vesicles purified from cell culture supernatant of K562 cells (Detected by anti CD63 antibody)

Concentration (ng/mL)	Abs.(450 nm - 620 nm)			
	1	2	Mean	∠Blank
0	0.0462	0.0528	0.0495	-
1.25	0.0870	0.0911	0.0891	0.0396
2.5	0.1360	0.1446	0.1403	0.0908
5	0.2355	0.2387	0.2371	0.1876
10	0.4320	0.4205	0.4263	0.3768
20	0.8192	0.7850	0.8021	0.7526
40	1.5618	1.5124	1.5371	1.4876
80	2.9187	2.8005	2.8596	2.8101

■4 Reference data : Detection sensitivity of extracellular vesicles purified from cell culture supernatant
Detection sensitivity of extracellular vesicles purified from each sample of cell culture supernatant by the assay (Detected by anti CD63 antibody)

	Detection limit (Blank +3.3SD)	Determination limit (Blank +10SD)
Purified exosomes of COLO201	0.11 ng/mL (11pg)	0.34 ng/mL (34pg)
Purified exosomes of K562	0.50 ng/mL (50pg)	1.34 ng/mL (134pg)

< Reference comparative data >

Detection sensitivity of extracellular vesicles purified from cell culture supernatant of COLO201 cells by Western blot (Detected by anti CD63 antibody)

Purified exosomes of COLO201


Fig. 4 Result of Western blot analysis
 sample : extracellular vesicles purified from cell culture supernatant of COLO201 cells with MagCapture™ Exosome Isolation Kit PS
 ★ : Detection limit

■5 Reference data : Precision of the quantitative analysis of cell culture supernatant samples
Intra-assay error (Detected by anti CD63 antibody)

A standard curve was prepared using extracellular vesicles purified from cell culture supernatant of COLO201 cells and K562 cells as a standard. Then concentrations of 4-step dilution samples of cell culture supernatant of K562 cells (1:100 to 1:800) and 4-step dilution samples of cell culture supernatant of COLO201 cells (1:200 to 1:1600) were measured (n = 6) and CV (%) of measurement values were calculated.

Cell culture sup. of COLO201 (Dilution ratio)	Assay value (ng/mL)								
	1	2	3	4	5	6	Mean	SD	CV(%)
1 : 1600	1.18	1.17	1.16	1.17	1.15	1.14	1.16	0.01	1.3
1 : 800	2.31	2.32	2.32	2.30	2.31	2.32	2.31	0.01	0.4
1 : 400	4.82	4.78	4.77	4.78	4.81	4.78	4.79	0.02	0.4
1 : 200	9.94	9.89	9.76	9.85	9.77	9.96	9.86	0.09	0.9

Cell culture sup. of K562 (Dilution ratio)	Assay value (ng/mL)								
	1	2	3	4	5	6	Mean	SD	CV(%)
1 : 800	3.93	3.88	3.88	3.83	3.78	3.86	3.86	0.05	1.3
1 : 400	7.71	7.82	7.64	7.64	7.48	7.46	7.62	0.13	1.8
1 : 200	15.47	14.89	15.32	15.32	14.85	14.69	15.09	0.32	2.1
1 : 100	32.03	31.77	31.21	31.59	30.31	30.03	31.16	0.82	2.6

Inter-assay error (Detected by anti CD63 antibody)

A standard curve was prepared using extracellular vesicles purified from cell culture supernatant of COLO201 cells and K562 cells as a standard. Then concentrations of 4-step dilution samples of cell culture supernatant of K562 cells (1:100 to 1:800) and 4-step dilution samples of cell culture supernatant of COLO201 cells (1:200 to 1:1600) were measured by 3 different assays and CV (%) of measurement values were calculated.

Cell culture sup. of COLO201 (Dilution ratio)	Assay value (ng/mL)					
	1	2	3	Mean	SD	CV(%)
1 : 1600	1.09	1.07	1.12	1.09	0.02	2.1
1 : 800	2.13	2.13	2.19	2.15	0.03	1.6
1 : 400	4.26	4.39	4.42	4.36	0.09	2.0
1 : 200	8.65	8.80	8.96	8.81	0.16	1.8


Cell culture sup. of K562 (Dilution ratio)	Assay value (ng/mL)					
	1	2	3	Mean	SD	CV(%)
1 : 800	3.10	3.30	4.09	3.50	0.52	15.0
1 : 400	6.53	6.76	7.45	6.91	0.48	6.9
1 : 200	13.44	14.53	15.27	14.41	0.92	6.4
1 : 100	27.09	29.19	31.34	29.21	2.12	7.3

■ 6 Reference data : Dilution linearity of cell culture supernatant samples

A standard curve was prepared using extracellular vesicles purified from cell culture supernatant of COLO201 cells and K562 cells as a standard. Then dilution linearity of 5-step dilution samples of cell culture supernatant of COLO201 cells (1:100 to 1:1600) and 5-step dilution samples of cell culture supernatant of K562 cells (1:50 to 1:800) was evaluated.


Dilution linearity of cell culture supernatant samples of COLO201 cells (Detected by anti CD63 antibody)

Dilution		Assay value (ng/mL)	Expected value (ng/mL)	% of expected
Ratio	Factor (x)			
1 : 1600	0.000625	0.89	0.91	98.4
1 : 800	0.00125	1.82	1.72	105.6
1 : 400	0.0025	3.44	3.52	97.8
1 : 200	0.005	7.04	6.78	103.9
1 : 100	0.01	13.6	-	-

Dilution linearity of cell culture supernatant samples of K562 cells (Detected by anti CD63 antibody)

Dilution		Assay value (ng/mL)	Expected value (ng/mL)	% of expected
Ratio	Factor (x)			
1 : 800	0.00125	3.02	3.07	98.2
1 : 400	0.0025	6.15	6.19	99.3
1 : 200	0.005	12.4	12.7	97.5
1 : 100	0.01	25.4	26.1	97.2
1 : 50	0.02	52.2	-	-

■ 7 Reference data : Spike and Recovery test with cell culture supernatant samples
Spike and Recovery on 1:400 diluted cell culture supernatant samples of COLO201 cells

Three concentrations (2.5, 5, 10 ng/mL) of extracellular vesicles purified from cell culture supernatant of COLO201 cells were added to 1:400 diluted samples of cell culture supernatant of COLO201 cells and concentrations of samples were measured to calculate the recovery rate.

Diluted samples of cell culture supernatant of COLO201 (1:400)

Spiked value (ng/mL)	Assay value (ng/mL)	Recovery value (ng/mL)	Recovery rate (%)
-	4.23	-	-
2.5	7.03	2.81	112.2
5.0	9.19	4.96	99.2
10.0	14.3	10.1	100.5

Spike and Recovery on 1:200 diluted cell culture supernatant samples of K562 cells

Three concentrations (10, 20, 40 ng/mL) of extracellular vesicles purified from cell culture supernatant of K562 cells were added to 1:200 diluted samples of cell culture supernatant of K562 cells and concentrations of samples were measured to calculate the recovery rate.

Diluted samples of cell culture supernatant of K562 (1:200)

Spiked value (ng/mL)	Assay value (ng/mL)	Recovery value (ng/mL)	Recovery rate (%)
-	14.3	-	-
10.0	24.6	10.3	103.0
20.0	33.8	19.5	97.5
40.0	53.4	39.1	97.8

【15. Precautions】

- Wear gloves, glasses, and laboratory coats when handling the assay materials.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. In case of contact with these, wash skin/eyes thoroughly with water and seek medical attention when necessary.
- In treating assay samples of human or animal origin, be careful for possible biohazards. This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Used samples and tips should be rinsed in 1% formalin, 2% glutaraldehyde or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by autoclaving before disposal. Dispose consumable and unused component in accordance with applicable regional/national regulatory requirements.
- In order to avoid dryness of wells, contamination of foreign substances, bias of temperature, and evaporation of dispensed reagents during incubation, never forget to label a plate seal.
- ELISA can be easily affected by the laboratory environment. Follow the operating method and temperature in the laboratory and in the incubator should be strictly maintained at 20-25 °C. In addition, avoid airstream velocity over 0.4 m/sec including wind from an air conditioner and humidity less than 30%.
- Measurements should be conducted considering cross contamination between wells.

【16. Related Products】

Code No.	Description	Size
299-77603	MagCapture™ Exosome Isolation Kit PS	2 purifications
293-77601		10 purifications
298-80601	PS Capture™ Exosome ELISA Kit (Streptavidin HRP)	96 reactions
297-79701	PS Capture™ Exosome Flow Cytometry Kit	300 reactions
290-80301	PS Capture™ Exosome Isolation Resin Kit	1 Kit (0.5 mL Slurry)
016-27061	Anti CD63, Monoclonal Antibody (3-13)	20 µL
012-27063		100 µL
018-27641	Anti CD63, Monoclonal Antibody (3-13), Fluorescein Conjugated	25 tests
014-27643		100 tests
013-27711	Anti CD63, Monoclonal Antibody (3-13), Biotin Conjugated	20 µL
019-27713		100 µL
011-27751	Anti CD63, Monoclonal Antibody (3-13), Red Fluorochrome (635) Conjugated	25 tests
017-27753		100 tests
018-27761	Anti CD9, Monoclonal Antibody (1K)	20 µL
014-27763		100 µL
015-27771	Anti CD81, Monoclonal Antibody (17B1)	20 µL
011-27773		100 µL
058-09261	EV-Save™ Extracellular Vesicle Blocking Reagent	1 mL

【17. References】

- 1) W. Nakai, *et al.*, *Sci. Rep.*, **6**, 33935 (2016).
- 2) A. Kraft, W81 Award. (2017).
- 3) S. Osawa, *et al.*, *Biochem. Biophys. Res. Commun.*, **488**, 232-238 (2017).
- 4) S. Nagashima, *et al.*, *J. Virol.*, **91**(22), e00822-17 (2017).
- 5) T. Yoshida, *et al.*, *Curr. Protoc. Cell Biol.*, **77**(1), 3.45.1-3.45.18 (2017).
- 6) S. Saito, *et al.*, *Sci. Rep.*, **8**, 3997 (2018).
- 7) H. Kobayashi, *et al.*, *Nagoya J. Med. Sci.*, **80**(2), 141-153 (2018).
- 8) Y. Obata, *et al.*, *JCI Insight*, **3**(8), e99680 (2018).
- 9) R.C. Middleton, *et al.*, *J. Extracell. Vesicles*, **7**(1), 1456888 (2018).
- 10) H. Kawahara, and R. Hanayama, *Biol. Pharm. Bull.*, **41**, 1119-1125 (2018).
- 11) M. Santiana, *et al.*, *Cell Host & Microbe*, **24**, 208-220 (2018).
- 12) T. Antes, *et al.*, *J. Nanobiotechnol.*, **16**, 61 (2018).
- 13) H. Kimura, *et al.*, *Clinical Cancer Research*, **25**(6), 1936-1947 (2019).
- 14) G. Gorji-Bahri, *et al.*, *Current Gene Therapy*, **18**, 336-350 (2018).
- 15) H. Xu, *et al.*, *Anal. Chem.*, **90**(22), 13451-13458 (2018).
- 16) H. Kuroda, *et al.*, *Mol Pharm.*, **16**, 292-304 (2019).
- 17) J. Juan, *et al.*, *Stem. Cell Res. Ther.*, **10**(1), 95 (2019).
- 18) L. Zhi, *et al.*, *Cell Death Discovery*, **5**, 79 (2019).
- 19) G.K. Patel *et al.*, *Sci Rep.*, **9**(1), 5335 (2019).
- 20) Y.L. Tai *et al.*, *J. Biomed. Sci.* **26**(1), 35 (2019).
- 21) Y. Sugihara *et al.*, *PLoS One.*, **14**(5), e0217394 (2019).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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