LabAssay™ Phospholipid (Choline Oxidase•DAOS method)

For Research Use Only. Not for use in diagnostic procedures.

[Introduction]

Phospholipids are a class of lipids, and a major component of all biological membranes, along with glycolipids, cholesterol and proteins. LabAssay™ Phospholipid is based on an enzymatic methods using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of phospholipids in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

[Kit contents]

	Kit contents	State	Amount
(1)	Buffer Solution (50 mmol/L Good' s buffer, pH7.5)	liquid	50 mL/3 bottles
(2)	Chromogen Substrate After reconstitution Phospholipase D 0.47 units/mL Choline Oxidase 2.0 units/mL Peroxidase 4.2 units/mL <i>N</i> -Ethyl- <i>N</i> -(2-hydroxy -3- sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) 0.77 mmol/L 4-Aminoantipyrine 0.24 units/mL Ascorbate Oxidase 3.9 units/mL	freeze dry	For 50 mL/3 bottles
(3)	Standard Solution Choline Chloride 54 mg/dL (Corresponding to 300 mg/dL Phospholipid)	liquid	5 mL/1 bottles


※Assay in a microplate;500 tests. Assay in a test tube; 50 tests.

[Storage and expiration]


When the complete kit is stored at 2 °C ~ 10 °C, the kit is stable until the expiration date shown on the label on the container.

[Assay principle]

Phospholipids (lecithin, sphingomyelin and lysolecithin) in a sample are hydrolyzed to choline in a reaction catalyzed by phospholipase D. Choline so formed is oxidized by choline oxidase in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalyzed by peroxidase(POD), producing a blue pigment. The amount of phospholipids contained in the sample is determined by measuring the absorbance of the blue color.

LabAssay™ Phospholipid (Choline Oxidase-DAOS method)

[Materials and apparatuses to be prepared]

- 96wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37 °C *
- Plate mixer *
- Microplate reader with 600 nm wavelength filter (* if the microplate reader is not equipped.)

(For test tube method)

- Test tube ·Pipette ·Incubator maintaining at 37 °C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

[Preparation of reagents to be used]

①Chromogen Substrate :

Dissolve 1 vial of Chromogen Substrate with 50 mL of Buffer Solution. After reconstitution, the solution should be stored at 2 °C~10 °C and used within one week.

②Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard

No.	Phospholipid Standard Solution	Distilled or deionized water	Sample volume	Phospholipid Concentration
1	50 µL	150 µL	2 µL	75.0 mg/dL
2	100 µL	100 µL	2 µL	150 mg/dL
3	Undiluted	—	2 µL	300 mg/dL
4	Undiluted	—	4 µL	596.1 mg/dL* ¹

* 1 The test sample volume is usually 2 µL, but 4 µL is taken in this case. The phospholipids concentration must be corrected accordingly as indicated in the table above.

[Procedure]

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	300 µL	300 µL	300 µL
Sample	Serum 2 µL	Standard solution 2 µL	—
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance at 600 nm* ² of the test sample and standard solution with the blank solution as the control.		

* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Chromogen reagent	3.0 mL	3.0 mL	3.0 mL
Sample	Serum 20 µL	Standard solution 20 µL	—
	Mix well and incubate at 37 °C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm. Spectrophotometer : 600 nm*2		

* 2 In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

[Standard curve]

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.14, when measuring purified water as a sample.
- The absorbance is 0.12 to 0.58, when measuring 300 mg/dL phospholipid as a sample.

(2) Specificity

The phospholipid concentration is less than $\pm 15\%$, when measuring the known concentration of a control serum as a sample.

[Usage Notes]

(1) Sample

- Samples should be used immediately after collecting.
- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.
- This method measures lecithin, sphingomyelin, lysolecithin but not cephalin.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.

LabAssay™ Phospholipid (Choline Oxidase·DAOS method)

- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The reaction completes in about 3 min. It is enough to incubate at 37 °C for 5 min. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.

(4) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

[Reference]

Takayama, M., Itoh, S., Nagasaki, T. and Tanimizu, I. : *Clin. Chim. Acta* .,79, 93 (1977) .

LabAssay™ Phospholipid (Choline Oxidase·DAOS method)

[Storage condition] Store the kit at 2 °C~10 °C

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