

Wako

# EZGlyco<sup>®</sup> O-Glycan Prep Kit

Revolutionary steamlined, rapid and reliable O-glycan sample preparation

#### Introduction

EZGlyco<sup>®</sup> O-Glycan Prep Kit provides a streamlined, rapid and reliable O-glycan sample preparation. The complete kit includes all reagents to release, enrich, label and purify O-glycans for HPLC and LC-MS analysis. The kit adopts a unique O-glycan releasing method using a high reactive amine and an organic superbase, which yields high recovery and less degradation (peeling) of released O-glycans.

Applications include characterization of O-glycosylation of recombinant glycoproteins in upstream R&D, quality control and process control. It is also suited for O-glycan biomarker discovery in biological samples such as mucins, serum and exosome, and O-glycan analysis of functional foods such as dairy products.



# Advantages of EZGlyco® O-Glycan Prep Kit

# O-glycan Sample Prep in 5 Hours

O-glycan sample preparation is completed in 5 hours with its unique O-glycan releasing reagents. Thus, HPLC/LC-MS analysis is achieved on the same day. The kit contains 2-AB labeling reagent for a sensitive detection with LC and LC-MS analysis.

# **High Recovery**

A combination of the unique O-glycan releasing reagent and O-glycan Capturing Beads enables high recovery of released O-glycans from various glycoproteins.

# Minimized "Peeling"

Glycan releasing reagents minimize decomposition of O-glycans (peeling), resulting in O-glycan characterization with higher accuracy.

## Easy and Safe Operation

With a simple protocol to follow, O-glycan sample preparation is completed without any special laboratory equipment.



# Workflow of EZGlyco® O-Glycan Prep Kit

For HPLC and LC-MS-ready sample preparation – The streamlined operation can be carried out in **5 hours** without any complicated processes and need for special laboratory instruments.



## **Comparision Data**

Comparison of 2-AB labeled bovine fetuin O-glycan analysis prepared with 2 methods: EZGlyco<sup>®</sup> O-Glycan Prep kit vs Hydrazinolysis



## **Reference data: Dose-dependency (Linearity)**

- Varying amount of bovine fetuin (0.5 200 μg) was subjected to O-glycan preparation using EZGlyco<sup>®</sup> O-Glycan Prep Kit (N=1).
- 2. 1 µL of the recovered solution containing O-glycans was analyzed using HILIC mode UHPLC.



EZGlyco<sup>®</sup> O-Glycan Prep Kit.

## Reference data: Reproducibility

 20 μg of bovine fetuin was subjected to O-glycan preparation using EZGlyco<sup>®</sup> O-Glycan Prep Kit (N=3 x 3 days).



#### Intra-assay variability

Peak ratio of Day1								
Peak#	entry1	entry2	entry3	3 CV				
1	6.3%	6.3% 6.3%		0.2%				
2	73.1%	73.1%	73.0%	0.1%				
3	17.6%	17.6%	17.6%	0.3%				
4	4 3.1% 3		3.1%	0.8%				
Total peak ratio of Day1								
	entry1	entry2	entry3	cv				
Total peak area	605,921	609,083	606,215	0.3%				

Average of peak ratio (N=3)							
Peak#	Day1	Day2 Day3		cv			
1	6.3%	6.2% 6.3%		0.4%			
2	73.1%	73.1% 72.8%		0.3%			
3	17.6% 17.5% 17.		17.7%	0.6%			
4	3.1%	3.1%	3.1% 3.2%				
Average of total peak area (N=3)							
	Day1	Day2	Day3	cv			
Total peak area	607,073	613,145	710,713	9.0%			

EZGIyco<sup>®</sup> O-Glycan Prep Kit demonstrates glycan profile CV < 5% and total recovery CV< 10%.

### **Performance Test**

2-AB labeled bovine fetuin O-glycan analysis prepared with different sample preparation kits.

- 1. 200 μg of bovine fetuin was subjected to O-glycan preparation using a) EZGlyco<sup>®</sup> O-Glycan Prep Kit, b) Product A (alkaline β-elimination), and c) Product B (alkaline β-elimination).
- 2. Recovered O-glycans were analyzed using HILIC mode UHPLC.



## Application #1

O-glycan profiles in standard glycoproteins

- Each sample was dissolved in 10 μL of pure water and subjected to O-glycan preparation using EZGlyco<sup>®</sup> O-Glycan Prep Kit.
- 2. 1 µL of the recovered solution containing O-glycans was analyzed using HILIC mode UHPLC

O-glycans were detected in all samples with a low peeling ratio (less than 10%).

#### Application #2

O-glycan profiles in human serum

- 20 μL of human serum was dried using a centrifugal evaporator, the dried human serum was dissolved in 10 μL of pure water and subjected to O-glycan preparation using EZGlyco® O-Glycan Prep Kit.
- 2. 1 µL of the recovered solution containing O-glycans was analyzed using HILIC mode UHPLC.

Low sample amount of human serum is sufficient for O-glycan analysis with EZGlyco<sup>®</sup> O-Glycan Prep Kit.

## **Ordering Information**

	Wako Code Product name		Package size						
	635-46299 EZGlyco® O-Glycan Prep Kit		10 tests						
Kit	t components:								
1.	Protocol		5.	Filter Column					
2.	Glycan Release	Reagent A	6.	2-Aminobezamide					
3.	Glycan Release	Reagent B	7.	Reducing Reagent					
4.	Glycan Capturir	ig Beads	8.	Cleanup Column					
Required Equipment, Labware, and Reagents:									
•	Acetic acid (AcC reagent grade*	DH), • 1	• F	ilter Column	•	Heating block for use at			
			-Aminobezamide		37°C and 50°C				
•	Acetonitrile (AC	N),	• F	Reducing Reagent	•	Vortex mixer			
	reagent grade*	• 1	.5-mL microcentrifuge	•	Microcentrifuge (used at				
<ul> <li>Methyl alcohol ( reagent grade*</li> </ul>		MeOH), • P		Pipette and tips for		500 and 3,000 x g)			

\*These products are also available at FUJIFILM Wako. For more information please visit www.e-reagent.com

1000, 200, and 20 µL

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit our online catalog to search for other products from FUJIFILM Wako; https://labchem-wako.fujifilm.com / This leaflet may contain products that can not be exported to your country due to regulations. / Bulk quote requests are welcome. Please contact us.

#### References •

1. Kameyama A., et al., *A practical method of liberating O-linked glycans from glycoproteins using hydroxylamine and an organic superbase,* Biochem Biophys Res Commun. 513(1), 186-192 (2019)

#### Manufacturer •

#### Sumitomo Bakelite Co., Ltd.

5-8, Higashi-Shinagawa 2-Chome, Shinagawa-ku, Tokyo, 140-0002, Japan

Telephone:+81-3-5462-4831Fax:+81-3-5462-4835

E-mail: s-bio@sumibe.co.jp

#### **Distribution & Inquiries**

FUJIFILM Wako Chemicals Europe GmbH Fuggerstr. 12 • 41468 Neuss, Germany

Telephone: +49-2131-311-271 Fax: +49-2131-311-110

E-mail: labchem\_wkeu@fujifilm.com

#### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan

Telephone: +81-6-6203-3741 Fax: +81-6-6203-1999

E-mail: ffwk-cservice@fujifilm.com



## https://labchem-wako.fujifilm.com

# FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road Richmond • VA 23237, USA

Telephone: +1-804-714-1920 Fax: +1-804-271-7791

E-mail: wkuslabchem@fujifilm.com