



# **Tissue Optical Clearing Reagent**

Whole-brain/ body imaging with single-cell resolution

# CUBIC



# **CUBIC Trial Kit Summary**

Dr. Hiroki R. Ueda *et al.* developed a simple, efficient, and scalable brain / body clearing method and imaging / computational analysis pipeline, CUBIC (clear, unobstructed brain/body imaging cocktails and computational analysis). CUBIC offers a fluorescent-preserving, high-performance and device-free tissue clearing method based on hydrophilic Solutions (ScaleCUBIC Solutions).

It enables reproducible whole-organ and whole-body clearing and rapid 3D imaging with LSFM. CUBIC also provides processing and analysis of 3D images for extracting biological information.

Therefore, CUBIC presents a platform of whole-organ/body imaging with single-cell resolution and image informatics, which enables a wide range of users to perform experiments targeting cellular and organ layers with multiple samples.

# **Benefits of CUBIC**

- Easy-to-use
- •No special equipment required
- •Follow up 2D H&E

- High transparency
- Applicable for many organs
- Compatible with IF, FP and other fluorescent labels

# **CUBIC Applications**

\*\*Typically, 7:3 mixture of Mounting Solution 1 and 2 (RI ~1.49) can be used.

	Fixation	Wash	Preprocessing	Delipidation		Wash	Permeabilization	Clearing	Imaging
	Perfusion fixation + post fixation 4%PFA/PBS(-)	PBS	50 % Scale CUBIC-1 Solution	Scale CUBIC-1 Solution	Scal e CUBIC-1 Solution	PBS	50 % Scale CUBIC-2 Solution	Scal e CUBIC-2 Solution	Mounting**
Processing temperature	4 ℃	RT	RT	37 ℃	37 ℃	RT - 37 ℃	RT	RT - 37 ℃	RT
Processing time	1 day	1 day	1 day	2 days	6 - 8 days*	1 day	1 day	2 days	



Figure 1. Transmittance images of mouse whole organs before and after clearing with ScaleCUBIC Solutions.





Figure 2. Transmittance images of mouse 2 mm slices before and after clearing with ScaleCUBIC Solutions.

# **CUBIC Application for fluorescent protein expression**


	Fixation	Wash	Preprocessing	Delipidation		Wash	Fluorescent labels	Wash	Permeabilization
·	Perfusion fixation + post fixation 4%PFA/PBS(-)	PBS	50 % Scale CUBIC-1 Solution	Scale CUBIC-1 Solution	Scale CUBIC-1 Solution	PBS	PI(5 μg / mL)	PB(10 mM)	50 % Scale CUBIC-2 Solution
Processing temperature	4 ℃	RT	RT	37 ℃	37 ℃	RT - 37 ℃	RT	RT	RT
Processing time	1 day	1 day	1 day	2 days	6 - 8 days*	1 day	2 days	> 1 day	1 day
	Clearing	Imaging							

Clearing	Imaging		
Sca/e CUBIC-2 Solution	Mounting:**		
RT - 37 ℃	RT		
2 days			

# An enlarged view of brain



An enlarged view of heart






Figure 3. LSFM imaging of CAG-EGFP Tg Mouse (8 w Male) whole brain, whole heart and dissected liver.




Figure 4. Tissue clearing of various human organs with CUBIC Tissue clearing and 3D imaging of deparaffinized samples with CUBIC




Figure 5. Tissue clearing and 3D imaging of pathological specimens with CUBIC.

### [References]

- 1) E. A. Susaki.et al.: Cell, 157(3),726 (2014).
- 2) K. Tainaka.et al.: Cell, 159(4), 911 (2014).
- 3) E. A. Susaki.et al.: *Nature Protocols*, **10**,1709 (2015).
- 4) S. Nojima.et al.: Scientific Reports ,7,9269 (2017).

Product name	Wako cat. No.	Package Size	Storage	Grade
CUBIC Trial Kit	290-80801	1kit	Keep at RT	Tissue Optical Clearing Reagent

- Listed products are intended for laboratory research use only, and not to be used for drug, food or human use.
- This brochure may contain products that cannot be exported to your country due to regulations.
- Bulk quote requests for some products are welcome. Please contact us.

# **FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan

Tel: 81-6-6203-3741 Fax: 81-6-6203-1999

# **FUJIFILM Wako Chemicals U.S.A. Corporation**

Toll-Free (U.S. only): 1-877-714-1920 **Head Office** (Richmond, VA)

Tel: 1-804-714-1920 Fax: 1-804-271-7791

# **FUJIFILM Wako Chemicals Europe GmbH**

### **European Office:**

Fuggerstraβe 12, D-41468, Neuss, Germany

Tel: 49-2131-311-0 Fax: 49-2131-311100