

Antifade Mountant for fluorescent proteins /oil-immersion objectives

High-resolution images with

Super Clear Mount

Super Clear Mount is designed to minimize spherical aberrations, enabling high-resolution imaging with high-NA oil-immersion objective lenses. Super Clear Mount is particularly powerful in super-resolution microscopy. Super Clear Mount is useful not only for thick brain tissues, but also for thin samples for cell biology and histology; Super Clear Mount should work better than conventional mounting media with high-NA oilimmersion objective lenses. Super Clear Mount offers the following benefits:



• High-refractive index: ~1.52

- High-resolution imaging
- Compatible with CLSM and SRM
- Suppress photobleaching







Figure 1. High-resolution confocal images of YFP-H line mouse 220 um brain with Super Clear Mount.

Figure 3. HEK293T cells expressing membrane EGFP were stained with DAPI (nuclei) and MitoTracker (mitochondria) and imaged using 3D SR-SIM.



Figure 4. Super-resolution (Airyscan) images of the medulla of the bsh-Gal4 UAS-CD8GFP flies, mounted Super Clear Mount and imaged with an oil-immersion objetive.

Resolution and brightness is maintained throughout depths in Super Clear Mount

Super Clear Mount refractive index (RI) of \sim 1.52, which is optimized to that of cover glass, immersion oil. and objectives.

With Super Clear Mount, the fluorescence level is maintained throughout all depths.

Super Clear Mount can be used to obtain super – resolution images of fluorescently labeled cells or tissue samples on microscope slides or coverslips.



Figure 5. Point-spread-function (PSF) analysis of fluorescent microspheres (diameter, 100 nm) in Super Clear Mount.

Super Clear Mount Summary



Figure 6. Full width at half maximum (FWHM) for lateral and axial directions.



Figure 7. Confocal axial scans of rhodamine dye solutions from the interface with the glass coverslip.

Different medias $100 \times / 1.40 \text{ Oil}$		Super Clear Mount
immersion coverslip sample Non high-resolution images with different medias	Curing	Non-curing and mounted sample stock for 1yr >.
	Refractive index (RI)	1.52
	Recommend sample size	0 ~ 200 μm
ton ingli resolution inlages with different inclus.	Photobleaching protection	+++
Super Clear Mount 100 × / 1.40 Oil RI 1.52 1.52 1.52	Stock temperature	RT
	+++ : best	

High-resolution images with Super Clear Mount.

Product name	Wako cat. No.	Package Size	Storage Condition	Grade
Super Clear Mount	190-18661	10ml	Keep at RT	Cell Biology

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use.

This brochure may contain products that cannot be exported to your country due to regulations.

Bulk quote requests for some products are welcome. Please contact us.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan Tel: 81-6-6203-3741 Fax: 81-6-6203-1999

FUJIFILM Wako Chemicals U.S.A. Corporation

Toll-Free (U.S. only): 1-877-714-1920 Head Office (Richmond, VA) Tel: 1-804-714-1920 Fax: 1-804-271-7791

FUJIFILM Wako Chemicals Europe GmbH

European Office: Fuggerstraβe 12, D-41468, Neuss, Germany Tel: 49-2131-311-0 Fax: 49-2131-311100