

Wako

QuickGene **NUCLEIC ACID** ISOLATION SYSTEMS

PATENTED POROUS MEMBRANE FOR HIGH-PURITY ISOLATION

- + Quick and easy operation
- + Semi and fully automated processes
- + Intact DNA isolation from varying volume and diverse samples origins








INTRODUCTION

The **QuickGene series** uses a patented porous membrane to realize a **low pressure** nucleic acid extraction with high yield and purity.

The product line combines semi to fully automated systems adapted to customers' needs. This enables quick and easy isolation of nucleic acid from various samples including whole blood/tissues/cells/plants/virus and others, as well as diverse sample volumes.

The **QuickGene** innovative isolation method contributes to the reduction of the users' burden and supports the high throughput in routine analytical laboratory work.




Adsorption of nucleic acids

Owing to their hydrophilic properties, nucleic acids are adsorbed onto the membrane, while protein and lipids, which are comparatively hydrophobic, tend to seep out of the membrane.


CORE TECHNOLOGY

The nucleic acid adsorption medium used in QuickGene series is a patented porous membrane developed through application of advanced polymer membrane production technology. Compared to conventional isolation glass fiber membranes, the QuickGene material is ten times thinner, thus providing rapid and reliable isolation at low pressure, preventing damage of the final product.




DNA/RNA ISOLATION PROCESS

QUICK & EASY PROCESSING




AUTOMATED PROCESS DEPENDING ON PROTOCOLS

Completely automated DNA/RNA extraction process from primary sample to final DNA storage, ensuring no cross-contamination and application errors. The **QuickGene-Auto12S** and **QuickGene-Auto24S** benchtop device system provides short operation time due to the simultaneous processing of multiple samples using pre-filled reagent strips.




SEMI-AUTOMATED PROCESS

The compact and handy design of the **QuickGene-Mini8L** and **QuickGene-Mini480** enables gentle nucleic acid isolation without centrifugation or removal and transfer of liquid – saving valuable time. The 2-in-1 multifunction characteristics of the devices allow the exchange of adapters and holder sets as well as enable to process 1 to 48 samples per run, volume ranging from 200 µl to 2 ml.



SELECTION GUIDE



Application guide is available under: https://www.wako-chemicals.de/files/quickgene_application_guide.pdf

QuickGene-AutoS RNA cultured cell kit

AS-RC

| INSTRU | JMENTS | | Bar BuckCare | |
|-----------|-------------------|----------|--------------|--------------|
| CODE | PRODUCT | Ŀ | 0 | |
| 631-46159 | QuickGene-Auto12S | KIKURADO | | |
| 638-46169 | QuickGene-Auto24S | | 90 | |
| 636-31579 | QuickGene-Mini480 | | | |
| 628-05969 | QuickGene-Mini8L | | | |
| | | | | Timming Them |

7

KITS & CONSUMABLES

| CODE | PRODUCT | | | |
|-----------|----------------------------|--|--|--|
| 635-46179 | AS-DB AutoS DNA Blood Kit | | | |
| 632-46189 | AS-DT AutoS DNA Tissue Kit | | | |
| 630-23549 | DB-S DNA Whole Blood Kit S | | | |
| 633-23539 | DB-L DNA Whole Blood Kit L | | | |
| 637-23559 | DT-S DNA Tissue Kit | | | |
| 632-29589 | DT-L DNA Tissue Kit | | | |
| 631-23579 | RB-S RNA Blood Cell Kit S | | | |
| 638-23609 | RT-S2 RNA Tissue Kit SII | | | |
| 634-23569 | Plasmid PL-SII | | | |
| 635-23619 | SP DNA Whole Blood Kit | | | |
| 632-23629 | SP DNA tissue Kit | | | |
| 638-23589 | RNA Cultured Cell Kit S | | | |
| 635-23599 | RNA cultured cell HC Kit | | | |



MKURABO

For further information on our products or to place an order, please contact us.

FUJIFILM Wako Chemicals Europe GmbH Fuggerstr. 12 41468 Neuss, Germany

> Phone + 49 2131 311 271 Fax + 49 2131 311 110

labchem_wkeu@fujifilm.com www.wako-chemicals.de

