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Prepared: September 2019 (1st version)

Reagent for Research Use

Cell **Suspension** and Preservation Solution

Cellstor-S

Instructions for Use

Note:

- (1) This product is a reagent for research use and cannot be administered to humans.
- (2) This product is not intended for medical use in humans or animals.

Cellstor-S is a cell suspension and preservation solution containing dextran 40 and trehalose hydrate for refrigerated and room-temperature storage. More than 90% of the cells were confirmed to be viable when human adipose-derived mesenchymal stromal cells were suspended using the product and stored for 24 hours (storage conditions: cell concentration, 5×10^5 cell/mL; storage temperature, 5°C or 25°C). When human adipose-derived mesenchymal stromal cells are suspended in the product, the cells remain in suspension for at least 1 hour. 1) Since this product is formulated with dextran 40, it is not suitable for washing cells when performing centrifugation. In such cases, using "Cellstor-W" solution, which does not contain dextran 40, is recommended.

Composition and Description

1. Composition

The ingredients and their quantities in a bag (250 mL) are as follows:

	- '
Ingredients	250 mL
Dextran 40	12.5 g
Trehalose Hydrate	8.29 g
Calcium Chloride Hydrate	0.05 g
Potassium Chloride	0.075 g
Sodium Chloride	1.5 g
Sodium L-Lactate	0.775 g
pH Adjuster	Appropriate Amount
Water for Injection	Appropriate Amount

Electrolyte concentration (mEq/L)

Na ⁺	K ⁺	Ca ²⁺	Cl-	L-Lactate ⁻
130	4	3	109	28

This product is a clear, colorless to pale yellow, sterile aqueous solution. The pH and osmotic pressure ratio are as follows:

рН	6.0-8.5
Osmotic Pressure Ratio (relative to normal saline)	1.0-1.6

Usage

- 1. Remove the protective seal from the rubber stopper.
- 2. Disinfect the outer surface of the bag by spraying with ethanol
- · The sterility of the outer surface of the product bag (including the surface of the rubber stopper inside the protective seal) is not guaranteed.
- Disinfect all outer surfaces with ethanol when transferring the product to the cleanliness control area and sterile working area, etc. Do not sterilize by autoclaving, irradiation with γ-rays or electron beams, EOG, or hydrogen peroxide, as these may affect
- The product name is printed on the bag. In order to facilitate identification, use the attached label sticker as necessary.

- 3. To extract the solution, carefully insert a needle vertically into the center of the circle indicated on the rubber stopper, and withdraw the solution.
- · If the stopper is not punctured vertically, it may cause contamination of the sample with minute fragments of the rubber stopper, or leakage of the contents. Do not repeatedly puncture the rubber stopper with a needle at the same point.

*Precautions for Use

- 1. Perform all procedures aseptically.
- 2. Open the outer bag immediately before use.
- 3. The duration of storage varies depending on the type and state of the cells and the storage conditions. Perform a confirmation test on the cells to be used.
- 4. The product bag is not designed for cell storage. When storing cells, use a container suitable for cell storage.
- 5. If an oil-based ink pen is used on the surface of the product bag, organic solvents such as xylene may be eluted into the product. Therefore, do not write on the bag with an oil-based ink pen.
- 6. Cells suspended in the product are less likely to settle. When performing centrifugation, centrifuge after diluting cells 5- to 15-fold with Cellstor-W, sold separately, or culture medium.
- 7. When this product is administered extravascularly to animals, edema may occur due to the colloid osmotic pressure of dextran 40. (For extravascular administration, using Cellstor-W, sold separately, is recommended).
- 8. When this product is administered to rats in an animal study, dextran 40 may cause allergic symptoms such as edema.²⁾

Precautions for Handling

- 1. Do not freeze the product. Store the product at room temperature, as freezing it may damage the bag containing the product. In order to prevent precipitation of dextran, store the bag horizontally (not vertically) in an area that does not experience marked changes in temperature.
- 2. Avoid strong physical shock or contact with sharp objects as this may cause leakage.
- 3. Do not use the product in the following instances:
- If the bag is damaged
- · If foreign matter is observed in the bag
- · If insoluble dextran (scales or precipitate) is observed
- 4. When discarding the product, follow the disposal methods prescribed by the relevant local authority.

Storage and Expiration Date

Storage: Store at room temperature.

Expiration date: 3 years

Packaging

 $250 \text{ mL} \times 1 \text{ bag (soft bag)}$

*[References]

- 1) Fujita Y, et al. Regen Ther. 2020;14;95-102.
- 2) Delitheos AK, et al. Int Archs Allergy appl Immun. 1976;50:436-445.

Use Example (Storage method of human adipose-derived mesenchymal stromal cells)1)

- 1) Centrifuge the cells and remove the supernatant. (Wash the cells with Cellstor-W, sold separately, as necessary.)
- 2) Add the product and resuspend the cells.
- 3) Refrigerate or store the sample at room temperature.

Manufactured by:
Otsuka Otsuka Pharmaceutical Factory, Inc.
115 Kuguhara, Tateiwa, Muya-cho, Naruto, Tokushima, Japan