

Lysyl Endopeptidase[®], recombinant, Biopharmaceutical Analysis Grade (rLys-C)

This product is a recombinant Lysyl Endopeptidase[®] expressed in *E. coli*. Lysyl endopeptidase[®] is a serine protease that cleaves peptide bonds at the carboxy-terminus of lysine residues with high specificity.

Taking advantage of its excellent specificity, Lysyl Endopeptidase[®] is used for peptide fragmentation and peptide mapping for analysis of the primary structure of proteins.

This product is checked residual DNA and Host cell protein for biopharmaceutical analysis.

Features

- Expressed in *E. coli*
- Checked residual DNA / Host cell protein
- Missed cleavage rate is lower than competitor's Lys-C



Specification

Requirement	Specification
Appearance	Lyophilisate
Activity	≥2.0AU/mg
DNA residual test	≤10ng/mg
HCP assay	≤1.0μg/mg
Endotoxin testing	<20EU/mg
Electrophoresis test (SDS-PAGE)	to pass test

Cat. No.	Product Name	Grade	Pkg. Size	Storage
124-06871	Lysyl Endopeptidase [®] , recombinant, Biopharmaceutical Analysis Grade	for Cell Culture	20 μg/vial	Keep at -20°C.

BSA Digestion Analysis

We incubated BSA with each protease for 1 hour and 18 hours.

After incubations, we analyzed the missed cleavage rate and number of peptides.

	Our Product	Competitor's Product
Specificity (Missed cleavage rate after 1 hour's incubation)	0%	0%
Specificity (Missed cleavage rate after 18 hours' incubation)	10%	20%
Activity (Number of peptides after 18 hours' incubation)	41	35

➡ The data suggested that our rLys-C have higher activity and specificity in long time reaction than the competitor's product.

Procedure

1. Sample preparation

- ① Dissolve or dilute the protein sample to be digested with 25 mmol/L Tris-HCl, 1 mmol/L EDTA, pH 8.5~9.0.
- ② Add Disulfide threitol (DTT) or β -mercaptoethanol to the solubilized protein at the final concentration of 5 mmol/L.
- ③ Incubate for 30 min at room temperature.
- ④ Add iodoacetamide to the solubilized protein at the final concentration of 10 mmol/L.
- ⑤ Incubate in dark for 30 min at room temperature.

2. rLys-C preparation

- ① Dissolve lyophilized powder of rLys-C with 100 μ l 12.5 mmol/L Tris-HCl, pH 8.5~9.0.

3. Sample digestion

- ① Add the prepared rLys-C to sample solution according to the enzyme : protein mass ratio of 1:20~1:100, and incubate at 35°C~37°C for 2~18 hours.
- ② The final concentration of 0.5~1.0v/v% trifluoroacetic acid is added to stop the reaction.

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: <https://labchem-wako.fujifilm.com/> / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site
<https://labchem-wako.fujifilm.com>



FUJIFILM Wako Pure Chemical Corporation
1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan
Tel: +81 6 6203 3741 Fax: +81 6 6203 1999
ffwk-cservice@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation
1600 Bellwood Road, Richmond, VA 23237, U.S.A.
Toll-Free (U.S. only): +1 877 714 1920
Tel: +1 804 271 7677 Fax: +1 804 271 7791
wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH
Fuggerstr 12, 41468 Neuss, Germany
Tel: +49 2131 311 0 Fax: +49 2131 311 100
labchem_wkeu@fujifilm.com