## Mature BDNF ELISA Kit *Wako*, High Sensitive

### 1. Introduction

Brain-derived neurotrophic factor (BDNF) is known to mediate neurogenesis, neuroprotection, synapses formation, and other biological processes to play key roles in the brain. Therefore, BDNF is expected to serve as a clinical marker for depression and other psychological diseases. In addition to nerve-related diseases, it is reportedly associated with cardiac diseases such as heart failure, thereby becoming a research target in a broad range of areas. BDNF has a precursor, proBDNF, which undergoes processing to convert it to mature BDNF (mBDNF). As such, proBDNF and mBDNF have been reported to have distinct mechanisms of action. Mature BDNF ELISA Kit *Wako*, High Sensitive is an ELISA kit capable of specifically measuring mBDNF concentrations in a sample.

With the use of a luminescence detection system, this product is more sensitive than our other products: approximately 35 times more sensitive compared with our standard product (code No. 296-83201 Mature BDNF ELISA Kit *Wako*). In addition, because it has lower cross-reactivity to proBDNF than conventional products, specific and highly sensitive measurement of mBDNF is possible.

### 2. Performance of the Kit

Range of standard curve concentrations	0.116-50 pg/mL
Analyte	mature BDNF*
Analyte samples	Mouse serum, plasma, and brain lysate Rat serum and plasma Human serum, plasma, and saliva
Required volume	13 $\mu$ L (as diluted 4-fold)
Assay time	Approx. 4 hours
Detection method	Luminescent detection

\* This kit has 1.30% cross-reactivity to recombinant human proBDNF, and 0.32% cross-reactivity to recombinant mouse proBDNF.

### 3. Materials Supplied

Component	State	Volume
Antibody-coated Plate	Use after washing	1 plate 96-wells (8 $\times$ 12 wells)
Mature BDNF Standard	Use after reconstitution	1 vial
Buffer	Ready-to-use	60 mL/1 vial
Biotin-conjugated Antibody	Use after reconstitution	1 vial

Peroxidase-conjugated Streptavidin Solution	Use after dilution	100 $\mu$ L/1 vial
Luminescent Reagent 1	Ready-to-use	6 mL/1 vial
Luminescent Reagent 2	Ready-to-use	6 mL/1 vial
Wash solution (10 $\times$ )	Use after dilution	100 mL/1 vial
Plate Seal	Ready-to-use	4 sheets
Manual	—	1 book

### 4. Measurement Principle

The microplate is coated with anti-Mature BDNF monoclonal antibody. In each well, the standard solution or sample and Biotin-conjugated anti BDNF polyclonal antibody are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is measured to determine Mature BDNF in the sample.

### 5. Equipment or Materials Required but not Supplied

- Purified water (distilled water)
- Test tubes for dilution of standard solution and samples
- Glass utensils for dilution of Wash Solution (e.g. volumetric cylinders and beakers)
- Pipettes with disposable tips (one capable of accurately pipetting 10  $\mu$ L of liquid, and one capable of accurately pipetting 200-500  $\mu$ L of liquid)
- Continuous dispensing pipettes for 100  $\mu$ L of liquid
- Water absorbents such as paper towels (used to remove the liquid remaining on the washed microplate)
- Stirrer (vortex type)
- Microplate shaker (approximately 600 to 800 rpm)
- Automatic washer (preferred) or washing bottles for 96-well microplates
- 96-well microplate reader for luminescence detection
- Software for data analysis

### 6. Reagents Preparation

Be sure to return kit reagents to room temperature (20-25°C) before use (approximately 2 hours).

The reagents described as "Ready-to-use" in the 3. Materials Supplied section shall be used as they are after being returned to room temperature. For the reagents described as "Use after reconstitution" or "Use after dilution," prepare only as much reagent as needed for measurements as instructed below.

#### 6-1. Standard Solution

Dissolve the Mature BDNF Standard in the specified volume (see Appendix)<sup>\*1</sup> of purified water added to prepare a standard stock solution (10.0 ng/mL). Thereafter, prepare standard solutions at various concentrations by mixing purified water and the standard stock solution using the kit component buffer, previously returned to room temperature, as follows :

\*1 Since the volume of purified water to be added varies among different lots, please confirm the volumes specified in Appendix. The following is an example formulation.

Standard solution concentration (pg/mL)	Volume of standard solution	Volume of Buffer
500.0	stock solution (10.0 ng/mL) : 20 $\mu$ L	380 $\mu$ L

↓

50.0	Standard solution (500.0 pg/mL) : 50 $\mu$ L	450 $\mu$ L
18.2	Standard solution (50.0 pg/mL) : 100 $\mu$ L	175 $\mu$ L
6.61	Standard solution (18.2 pg/mL) : 100 $\mu$ L	175 $\mu$ L
2.40	Standard solution (6.61 pg/mL) : 100 $\mu$ L	175 $\mu$ L
0.874	Standard solution (2.40 pg/mL) : 100 $\mu$ L	175 $\mu$ L
0.318	Standard solution (0.874 pg/mL) : 100 $\mu$ L	175 $\mu$ L
0.116	Standard solution (0.318 pg/mL) : 100 $\mu$ L	175 $\mu$ L
0 (Blank)	-	175 $\mu$ L

#### 6-2. Biotin-conjugated Antibody

Dissolve in 100  $\mu$ L of purified water, and dilute 100-fold with the Buffer.

#### 6-3. Peroxidase-conjugated Streptavidin Solution

Dilute 100-fold with the Buffer.

#### 6-4. Luminescent Reagent 1 and Luminescent Reagent 2

Mix Luminescent Reagent 1 and Luminescent Reagent 2 in a 1 : 1 ratio (vol./vol.) and stored protected from light 15-30 minutes before use.

e.g. a mixture of Luminescent Reagent 1 (6 mL) and Luminescent Reagent 2 (6 mL) (when using all 96 wells)

#### 6-5. Wash Solution (10 $\times$ )

Dilute 10-fold with purified water (distilled water) before use. e.g., 100 mL of Wash Solution (10 $\times$ ) + 900 mL of purified water (distilled water) (when using all 96 wells)

○The other reagents are supplied Ready-to-use.

### 7. Reagent Stability and Storage

#### 7-1. Antibody-coated Plate

Unused antibody-immobilized strips should be returned to the attached zip seal pack and stored at 2-10 $^{\circ}$ C as they are. It is stable until the expiration date.

#### 7-2. Mature BDNF Standard

Store the prepared standard stock solution (10.0 ng/mL) at 2-10 $^{\circ}$ C, and use within 2 weeks. The standard solutions that have been diluted should be used immediately. Do not store. While in a freeze-dried state, it is stable until the expiration date.

#### 7-3. Buffer

When using a portion of the solution, transfer a volume slightly larger than required to another container, close the lid tightly, and store the remaining portion at 2-10 $^{\circ}$ C (do not return to room temperature). It is stable until the expiration date.

#### 7-4. Biotin-conjugated Antibody

Store the Biotin-conjugated Antibody in solution in purified water at 2-10 $^{\circ}$ C, and use within 2 weeks. The solution that has been diluted with the buffer should be used immediately. Do not store.

While in a freeze-dried state, this product is stable until the expiration date.

#### 7-5. Peroxidase-conjugated Streptavidin Solution

When using the kit in several divided portions, take out the solution from the refrigerator, dilute, immediately close the lid tightly, and store the remaining stock solution at 2-10 $^{\circ}$ C (do not return to room temperature). This product is stable until the expiration date. Discard the unused portions of diluted solutions.

#### 7-6. Luminescent Reagent 1 and Luminescent Reagent 2

When using the kit in several divided portions, take out the solution from the refrigerator, prepare, immediately close the lid tightly, and store the remaining stock solution at 2-10 $^{\circ}$ C (do not return to room temperature). This product is stable until the expiration date.

#### 7-7. Wash Solution (10 $\times$ )

When storing the Wash Solution (10 $\times$ ), close the lid tightly, and store at 2-10 $^{\circ}$ C. This product is stable until the expiration date. Discard the unused portions of diluted Wash Solutions.

### 8. Method of Sample Preparation

After diluting the sample with Buffer, allow to stand at room temperature for 15 minutes before measurement. It is recommended that you optimize the dilution ratio to fall within the range of standard curve concentrations.

Sample	Dilution ratio
Mouse plasma and serum	4-fold
Mouse brain lysate Lysis buffer : RIPA buffer + Protease inhibitor	2-320-fold
Rat plasma	30-fold
Rat serum	50-fold
Human plasma	10-80-fold
Human serum	10-800-fold
Human saliva	2-8-fold

### 9. Assay Procedures

- (1) Remove the microplate protectant liquid, fill each well with a previously prepared Wash Solution, and wash four 4 times (\*①). Thereafter, invert the microplate on a paper towel or the like, and remove the liquid remaining in the wells by gentle tapping.
- (2) Dispense 50  $\mu$ L of the standard solution at each concentration to each well for standard measurement.
- (3) Dispense 50  $\mu$ L of the sample, previously diluted with Buffer in the optimum ratio, to each sample measurement well.
- (4) Stir using a microplate shaker or the like (\*②).
- (5) Apply a Plate Seal (\*③), and allow the microplate to stand

at room temperature (20-25°C) for 2 hours.

- (6) After completion of the reaction, discard the reaction liquid, fill each well with Wash Solution, and wash four 4 times (\*①).

Thereafter, invert the microplate on a paper towel or the like, and remove the liquid remaining in the wells by gentle tapping.

- (7) Dispense 50  $\mu$ L of Biotin-conjugated Antibody Solution to each well.

Stir using a microplate shaker or the like (\*②).

- (8) Apply a new Plate Seal (\*③), and allow to stand at room temperature (20-25°C) for 1 hour.

- (9) After completion of the reaction, perform the washing procedure described in (6).

- (10) Dispense 50  $\mu$ L of Peroxidase-conjugated Streptavidin Solution to each well.

Stir using a microplate shaker or the like (\*②).

- (11) Apply a Plate Seal (\*③), and allow the microplate to stand at room temperature (20-25°C) for 30 minutes.

- (12) After completion of the reaction, perform the washing procedure described in (6).

- (13) Dispense 50  $\mu$ L of prepared luminescent reagent mixture to each well.

Stir for 1 minute using a microplate shaker.


- (14) After stirring, determine luminescence intensity using a 96-well microplate reader for luminescence detection. It is recommended that measurements be taken over a period of 10 to 20 minutes after addition of luminescent reagent.

- (15) Generate a standard curve with the standard solution concentration (pg/mL) plotted on the X-axis, and luminescence intensity on the Y-axis. Read the concentration (pg/mL) corresponding to the luminescence intensity of the diluted sample. Multiply the concentration reading by the sample dilution ratio (e.g., 4-fold for mouse plasma and serum), and use the product as the assay value.

\*In computing using computer software, it is recommended to use a tertiary polynomial expression or four or five parameters.

For (\*①), (\*②), and (\*③), refer to the measurement procedure outline.

## 10. Standard Curve (Example)

## 11. Precautions

1. This kit should be used under the guidance of a person who has been trained in ELISA techniques or a trainer.
2. For manual measurements, this kit should be used by an operator who can stably obtain reproducible pipetting results.
3. Wear protective gloves, goggles, and a protective gown during preparatory work and operation of this kit.
4. Avoid contact between the skin and the reagents. In the event any reagent in this kit accidentally comes into contact with wounds, the eye, mouth, skin, and other parts of the body, take first-aid measures, such as by immediately washing off the reagent, and seeking medical attention as required.
5. Do not eat or smoke in the places where this kit is being used.
6. Samples should be carefully handled as they are potentially infectious. This kit contains animal-derived components.
7. Used samples, used consumables, etc. should be immersed in 1% formalin, 2% glutaraldehyde, or  $\geq$  0.1% sodium hypochlorite solution for 1 hour or more. Alternatively, such items may be autoclaved and discarded. Used consumables and unused chemicals should be discarded in accordance with the rules of your facility and local laws and regulations.
8. Reagents from different lots should not be used when preparing the mixture. During the standing reaction in each step, be sure to affix a Plate Seal to stop the well from drying out, prevent the entry of foreign matter, prevent uneven temperature distributions, and stop evaporation of dispensed reagents.
9. ELISA techniques are susceptible to the measurement environment. Room temperatures during measurements and the places used for standing reactions should be strictly controlled at 20-25 °C (on the experimental bench or in the incubator). And avoid measurements in an air flow (including air movement caused by an air-conditioner) and at low humidity. No heat sources (computers, incubators, etc.) should be placed close to the microplate to be reacted.

## 12. Assay Procedure Summary

Before performing any measurements, read through the operating manual, and confirm the sample conditions, measurement conditions, and measurement methods.

Bring the plate and reagents to room temperature (20-25°C) about 2 hours before use.

Preparation of Wash Solution (10×) : Dilute 10-fold with purified water, previously returned to room temperature.

Preparation of standard solutions (example) : Dissolve the Mature BDNF Standard in the specified volume (see Appendix)\* of purified water added to prepare a standard stock solution (10.0 ng/mL). Thereafter, prepare standard solutions at various concentrations using Buffer, previously returned to room temperature. For the volume of purified water to be added, refer to the Appendix.

\*Since the volume of purified water to be added varies among different lots, please confirm the volumes specified in Appendix. The following is an example formulation.

Concentration (pg/mL)	500.0	50.0	18.2	6.61	2.40	0.874	0.318	0.116	0
Standard solution (μL)	20	50	100*	100*	100*	100*	100	100	-
Buffer (μL)	380	450	175	175	175	175	175	175	175

\*Standard solution at a 1-step higher concentration.

**Antibody-coated Plate**

↓ Wash 4 times (\*①)

**Sample or standard solution** **50 μL/well**

↓ Stir the solution (\*②), react at room temperature (20-25°C) for 2 hours, and allow to stand (\*③).

\*Preparation of Biotin-conjugated Antibody (dissolve in 100 μL of purified water, and dilute 100-fold with Buffer, previously adjusted to room temperature)

↓ Wash 4 times (\*①)

**Biotin-conjugated Antibody** **50 μL/well**

↓ Stir the solution (\*②), react at room temperature (20-25°C) for 1 hour, and allow to stand (\*③).

\*Preparation of Peroxidase-conjugated Streptavidin Solution (dilute 100-fold with Buffer, previously adjusted to room temperature)

↓ Wash 4 times (\*①)

**Peroxidase-conjugated Streptavidin Solution** **50 μL/well**

↓ Stir the solution (\*②), react at room temperature (20-25°C) for 30 minute, and allow to stand (\*③).

\*Preparation of luminescent reagent (mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1:1 ratio (vol./vol.).

↓ Wash 4 times (\*①)

**Luminescent Reagent** **50 μL/well**

↓ Stir at room temperature (20-25°C) for 1 minute.

**Measurement of luminescence intensity (measure between 10 and 20 minutes after sample preparation)**

(\*①) In each washing operation, dispense the Wash Solution (1×) into wells, gently agitate the filled plate on the palm for about 10 seconds, and then empty the wells. After washing the wells 4 times consecutively, reverse the plate, and tap it against paper towel to remove the washing solution completely. After removal of the Wash Solution, immediately dispense the next

solution with care not to dry the wells. Use of a pipet set at the liquid volume of 300 μL may be appropriate for dispensing the Wash Solution into each well.

(\*②) Three repeats of agitation at 600 to 800 rpm for 10 seconds may be appropriate.

(\*③) After stirring is complete, cover with a Plate Seal. Remove the liner from the Plate Seal, and apply its adhesive side to the plate for affixation. Do not re-use any Plate Seal.

**[Storage]** Store at 2-10°C.

**[Expiration date]** Indicated on the label.

**[Package]** For 96 assays.

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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