LBIS™ Mouse/Rat High Molecular Weight Adiponectin ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Mouse/Rat High Molecular Weight Adiponectin ELISA Kit is a sandwich ELISA system for quantitative measurement of mouse or rat High Molecular Weight Adiponectin. This is for research use only.

2. Introduction

Adiponectin is one of cytokines secreted from adipocyte (adipocytokines, adipokines), and controls lipid metabolism and insulin sensitivity. Adiponectin is an important substance showing anti-diabetic, anti-atherosclerotic, and anti-inflammatory actions. In blood, adiponectin monomer molecules are associated to form trimer, hexamer or 12 - 18-mers. The monomer form of adiponectin is found only in adipocytes, and homomultimer, i.e. trimer, hexamer and 12 - 18-mers formed from adiponectin monomer are found in plasma. Trimer (low molecular weight complexes/LMW) is formed by interaction of non-covalent bonds of triple helix area and by hydrophobic interaction between globular Cq1 domain, and hexamer (middle molecular weight complexes/MMW) or larger complexes (higher molecular weight/HMW) are formed through disulfide bonds between cysteines at 39 of LMW-complexes.

Adiponectin is believed to influence on cell growth, angiogenesis, and tissue remodeling by isolating various growth factors through binding with distinct affinity. The affinity is dependent upon the type of complexes, LMW, MMW and HMW.

The blood levels of HMW were reported to reflect BMI, sex, effect of body weight decrease, glucose tolerance, liver insulin sensitivity, metabolic syndrome or DM2 more than total adiponectin levels. This means that assay for HMW is expected to be more useful for analysis of metabolic syndrome and DM2. FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse/Rat High Molecular Weight Adiponectin ELISA Kit is specific to HMW.

3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse/Rat High Molecular Weight Adiponectin ELISA Kit, standards or diluted samples are incubated in monoclonal anti-adiponectin antibody-coated wells to capture HMW adiponectin. After 2 hours incubation and washing, HRP (horse radish peroxidase)-labeled anti-adiponectin antibody is added, and incubated for 90 minutes. After washing, HRP-complex remaining in wells is reacted with a TMB Solution for 30 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to HMW adiponectin concentration. The standard curve is prepared by plotting absorbance against standard HMW adiponectin concentrations. HMW adiponectin concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

- Assay range

The assay range of the kit is 3.13 ng/mL - 200 ng/mL.

- Specificity

The kit uses two monoclonal antibodies specific to mouse/rat HMW adiponectin.

| Sample | Cross reaction | Sample | Cross reaction |
|-----------------------------|----------------|-----------------------|----------------|
| Mouse adiponectin (HMW) | 100% | Rat adiponectin (HMW) | 100% |
| Mouse adiponectin (Hexamer) | < 5% | Tested at 1000 ng/mL. | |

No cross reaction at 1000 ng/mL : Mouse adiponectin (Trimer), Mouse adiponectin (Monomer), Mouse MCH, Mouse TNF- α , Mouse INF- γ , Mouse Insulin, Mouse Leptin, Rat adiponectin (Monomer), Rat TNF- α , Rat INF- γ , Rat Insulin, Rat Leptin

- Precision of assay

Within assay variation (2 samples, 5 replicates assay), mean CV was within 10%.

- Reproducibility

Between assay variation (3 samples, 4 days, duplicate assay), mean CV was within 10%

- Recovery test

HMW Adiponectin was added in 3 concentrations to 2 serum samples and was assayed.

The recoveries were 94.4% - 105%

- Dilution test
2 serum samples were serially diluted by 3 steps.
The dilution curves showed excellent linearity. ($R^2=0.999$)

5. Reference Assay Data

| Strain | Week | Gender | numbers | Results (ng/mL) | | Remarks |
|----------|------|--------|---------|-----------------|-----|---|
| | | | | Mean | SD | |
| BALB/c | 6 w | ♂ | 10 | 2369 | 743 | Serum, ad libitum feeding ad libitum feeding |
| ICR | 6 w | ♂ | 10 | 2119 | 802 | Serum, ad libitum feeding Serum |
| CD (rat) | 8 w | ♂ | 9 | 3220 | 670 | Serum, ad libitum feeding Serum |

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for HMW Adiponectin levels independently.

6. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person.
- Wear gloves and laboratory coats when handling assay materials.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the Stop Solution because it is sulfuric acid. The Stop Solution and the TMB Solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic Stop Solution and TMB Solution containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

7. Reagents supplied

| Components | Use Status | Amount |
|---|-----------------------------------|----------------------|
| (A) Antibody-coated Plate | Use after washing. | 96 wells/1 plate |
| (B) Adiponectin Standard | Concentrated. Use after dilution. | 200 μ L/1 bottle |
| (C) Buffer | Ready for use. | 60 mL/1 bottle |
| (D) Peroxidase-conjugated Antibody Solution | Concentrated. Use after dilution. | 100 μ L/1 bottle |
| (F) TMB Solution | Ready for use. | 12 mL/1 bottle |
| (H) Stop Solution | Ready for use. | 12 mL/1 bottle |
| (I) Wash Solution (10 \times) | Concentrated. Use after dilution. | 100 mL/1 bottle |
| Plate Seal | – | 3 sheets |

[Storage and Stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Adiponectin Standard (2000 ng/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C) Buffer] & [(F) TMB Solution]

If not opened, store at 2°C - 10°C. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Peroxidase-conjugated Antibody Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10×)]

If not opened, store at 2°C - 10°C. Dispose any unused diluted buffer.

8. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water). Test tubes for preparation of standard solution series.

Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle). Pipettes (disposable tip type). One should be able to deliver 5 μL - 10 μL precisely, and another for 50 μL - 500 μL. Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 50 μL. Paper towel to remove washing buffer remaining in wells. A vortex-type mixer. A shaker for 96 well-plate (600 rpm - 1200 rpm). An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. A 96 well-plate reader (450 nm ± 10 nm, 620 nm : 600 nm - 650 nm) Software for data analysis.

9. Preparation of Samples

This kit is intended to measure HMW Adiponectin in mouse or rat serum, plasma or culture medium. The necessary sample volume for the standard procedure is 10 μL.

Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results.

Hemolytic and hyperlipemic serum samples are not suitable.

***To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 80 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter. Sample dilution should be carried out with the (C) Buffer using small test tube such as PP, PE or glass, before assay. You can choose dilution rate 25 - 50 × if necessary.

| Example of dilution : | Rate | 50× | (25×) |
|-----------------------|------|-----|-------|
| Sample (μL) | | 10 | (10) |
| Buffer solution (μL) | | 490 | (240) |

Storage and stability

Adiponectin in samples will be inactivated if stored at 2°C - 10°C. If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freezing and thawing.

10. Preparation of Reagents

◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Adiponectin Standard (2000 ng/mL)]

Make a serial dilution of original standard solution to prepare each standard solution. An example is shown below.

| Volume of standard solution | (C) Buffer | Concentration (ng/mL) |
|------------------------------|------------|-----------------------|
| Original solution : 50 μL | 450 μL | 200 |
| 200 ng/mL solution : 200 μL | 200 μL | 100 |
| 100 ng/mL solution : 200 μL | 200 μL | 50 |
| 50 ng/mL solution : 200 μL | 200 μL | 25 |
| 25 ng/mL solution : 200 μL | 200 μL | 12.5 |
| 12.5 ng/mL solution : 200 μL | 200 μL | 6.25 |
| 6.25 ng/mL solution : 200 μL | 200 μL | 3.13 |
| 0 (Blank) | 200 μL | 0 |

[(D) Peroxidase-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

[(I) Wash Solution (10×)]

Dilute 1 volume of the concentrated Wash Solution (10×) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated Wash Solution (10×) and 900 mL of deionized water (or distilled water).

11. Assay Procedure

Remove the cover sheet of the 96 well-plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the well with washing buffer and discard 3 times (*①), then strike the plate upside-down onto several layers of paper towels to remove residual buffer in the wells.
- (2) Pipette 50 μ L of diluted samples to the designated sample wells.
- (3) Pipette 50 μ L of standard solution to the wells designated for standards.
- (4) Shake the plate gently on a plate shaker (*②)
- (5) Stick a Plate Seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette 50 μ L of (D) Peroxidase-conjugated Antibody Solution to all wells, and shake as step (4).
- (8) Stick a Plate Seal (*③) on the plate and incubate the plate for 90 minutes at 20°C - 25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette 50 μ L of (F) TMB Solution to wells, and shake as step (4).
- (11) Stick a plate seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (12) Add 50 μ L of the (H) Stop Solution to all wells and shake as step (4).
- (13) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm*) immediately using a plate reader.


*Refer to 15. Summary of Assay Procedure for notes of *①, *② and *③.

12. Technical Tips

- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells \times 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

13. Calculations

- (1) Prepare a standard curve using semi-logarithmic or two-way logarithmic section paper by plotting absorbance* (Y-axis) against HMW adiponectin concentration (ng/mL) on X-axis. Physiological or pathological situation of animals should be judged comprehensively taking other examination results into consideration.
 - (2) Using the standard curve, read the HMW adiponectin concentration of a sample at its absorbance*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.
- * We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.
- (3) Assay results from this product and results from other commercially available mouse/rat (total) adiponectin ELISA kits cannot be compared because of the difference of the assay systems.

HMW adiponectin assay standard curve (an example)
Absorbance may change due to assay environment.

14. Trouble Shooting


- Low absorbance in all wells
Possible explanations :
 - 1) The standard or samples might not be added.
 - 2) Reagents necessary for coloration such as Peroxidase-conjugated Antibody Solution, or TMB Solution might not be added.
 - 3) Wrong reagents related to coloration might have been added. Wrong dilution of Peroxidase-conjugated Antibody Solution.
 - 4) Contamination of enzyme inhibitor (s).
 - 5) Influence of the temperature under which the kits had been stored.
 - 6) Excessive hard washing of the well plate.
 - 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- The OD of blank is higher than that of the lowest standard concentration (3.13 ng/mL)
Possible explanations :

Improper or inadequate washing. (Change washing repetition from 3 times to 4 - 6 times after the reaction with Peroxidase-conjugated Antibody Solution.)
- High coefficient of variation (CV)
Possible explanation :
 - 1) Improper or inadequate washing.
 - 2) Improper mixing of standard or samples.
 - 3) Pipetting at irregular intervals.
- Q-1 : Can I divide the plate to use it for the other testing?
A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon
- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?
A-2 : When we manufacture 96 well-plate, we put protective solution in wells.

15. Summary of Assay Procedure : Use as a check box

* First, read this instruction manual carefully and start your assay after confirmation of details.

- Bring the well-plate and all reagents to **20°C - 25°C for 2 hours**.
- Wash Solution (10×) concentrate must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.
- Adiponectin Standard dilution example : Mix the original standard solution by pipetting so as not to make bubbles, and dilute using the (C) Buffer (of 20°C - 25°C).

- Antibody-coated Plate
- ↓ Washing 3 times (*①) *⑤
- Diluted Samples / Standards 50 µL
- ↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))
- Dilute (D) Peroxidase-conjugated Antibody Solution to **1 : 100** by using (C) Buffer and use. (Dilute reagents during the first reaction.)
- ↓ Washing 3 times (*①) *⑤
- Peroxidase-conjugated Antibody Solution 50 µL
- ↓ Shaking (*②), Incubation for 90 minutes at 20°C - 25°C. (Standing (*③))
- ↓ Washing 3 times (*①) *⑤
- TMB Solution 50 µL
(After dispense, the color turns to blue depending on the concentration.)
- ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- Stop Solution 50 µL
(After dispense, the color turns to yellow depending on the concentration.)
- ↓ Shaking (*②) (Immediately shake.)
- Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately. (Ref. wave cancels the dirt in the back of plate.)

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution.

Standard of plate-washing pressure : 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

*② Guideline of shaking : 600 rpm - 1200 rpm for 3 times for 10 seconds.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.


*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet Example

| | Strip 1&2 | Strip 3&4 | Strip 5&6 | Strip 7&8 | Strip 9&10 | Strip 11&12 |
|---|------------|-----------|-----------|-----------|------------|-------------|
| A | 200 ng/mL | Sample 1 | Sample 9 | Sample 17 | Sample 25 | Sample 33 |
| B | 100 ng/mL | Sample 2 | Sample 10 | Sample 18 | Sample 26 | Sample 34 |
| C | 50 ng/mL | Sample 3 | Sample 11 | Sample 19 | Sample 27 | Sample 35 |
| D | 25 ng/mL | Sample 4 | Sample 12 | Sample 20 | Sample 28 | Sample 36 |
| E | 12.5 ng/mL | Sample 5 | Sample 13 | Sample 21 | Sample 29 | Sample 37 |
| F | 6.25 ng/mL | Sample 6 | Sample 14 | Sample 22 | Sample 30 | Sample 38 |
| G | 3.13 ng/mL | Sample 7 | Sample 15 | Sample 23 | Sample 31 | Sample 39 |
| H | 0 (Blank) | Sample 8 | Sample 16 | Sample 24 | Sample 32 | Sample 40 |

Assay Worksheet

16. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Mouse/Rat High Molecular Weight Adiponectin ELISA Kit

[Storage] Store the kit at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the label.

[Package] For 96 tests

[Cat #] 297-90701

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : +81-6-6203-3741
Facsimile : +81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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1000ng/mL 添加で交差性が認められなかったもの: Mouse Adiponectin (Trimer), Mouse Adiponectin (Monomer), Mouse MCH, Mouse TNF- α , Mouse INF- γ , Mouse Insulin, Mouse Leptin, Rat Adiponectin (Monomer), Rat TNF- α , Rat INF- γ , Rat Insulin, Rat Leptin

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