LabAssay™ Triglyceride

Please, read this instruction carefully before use.

1. Introduction

Lipids in serum consist of triglycerides, cholesterol, phospholipids, free fatty acids and slight amounts of fat-soluble components such as fat-soluble vitamins and carotenes. Triglycerides, as major components of very low density lipoprotein (VLDL) and chylomicrons, play an important role in metabolism as energy sources and transporters of dietary fat.

LabAssay™ Triglyceride is based on an enzymatic method using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of triglycerides in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

2. Kit Contents

| | Components | State | Amount |
|-----|---------------------|------------|---------------------|
| (1) | Buffer Solution | liquid | 105 mL/1 bottle |
| (2) | Chromogen Substrate | freeze dry | For 105 mL/1 bottle |
| (3) | Standard Solution | liquid | 4 mL/1 bottle |

※ Assay in a microplate ; 350 tests. Assay in a test tube ; 35 tests.


3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

4. Assay Principle

Triglycerides in a sample are hydrolyzed to glycerol and free fatty acids in a reaction catalyzed by lipoprotein lipase (LPL). Glycerol is converted to glycerol-3-phosphate by glycerolkinase (GK) in the presence of ATP.

Glycerol-3-phosphate so formed is oxidized by glycerol-3-phosphate oxidase (GPO) in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalyzed by peroxidase (POD), producing a blue pigment. The amount of triglycerides contained in the sample is determined by measuring the absorbance of the blue color.

5. Materials and Apparatuses to be Prepared

- 96wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37°C *
- Plate mixer *
- Microplate reader with 600 nm wavelength filter (* if the microplate reader is not equipped.)

(For Test Tube method)

- Test tube ·Pipette ·Incubator maintaining at 37°C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

6. Preparation of Reagents to be used

① Chromogen Reagent :

Prepare Chromogen Reagent by Dissolving 1 bottle of Chromogen Substrate (for 105 mL) to 105 mL of Buffer Solution. After reconstitution, the solution should be stored at 2°C - 10°C and used within 10 days.

② Standard Solution (Microplate method)

Standard Solution is prepared by dilution of the provided Standard Solution.

| No. | Standard Solution | Distilled or deionized water | Sample volume | Triglyceride Concentration |
|-----|-------------------|------------------------------|---------------|----------------------------|
| 1 | 100 μ L | 200 μ L | 2 μ L | 100 mg/dL |
| 2 | 200 μ L | 100 μ L | 2 μ L | 200 mg/dL |
| 3 | undiluted | – | 2 μ L | 300 mg/dL |
| 4 | undiluted | – | 4 μ L | 596 mg/dL ^{*1} |
| 5 | undiluted | – | 6 μ L | 888 mg/dL ^{*1} |

^{*1} The test sample volume is usually 2 μ L, but 4 or 6 μ L are taken in this case.

The triglyceride concentration must be corrected accordingly as indicated in the table above.

7. Preparation of Samples

- Samples should be used immediately after collecting.
- Hemolysis causes an increase of the triglycerides value.
- Ascorbic acid causes a slightly negative effect on the assay.
- Bilirubin may not significantly affect the assay.

8. Procedure

(1) Assay in a Microplate

Perform the assay in the wells according to the following table scheme.

| | Test | Standard | Blank |
|-------------------|--|--------------------------------|-------------|
| Chromogen reagent | 300 μ L | 300 μ L | 300 μ L |
| Sample | Serum 2 μ L | Standard solution 2 μ L | – |
| | Mix well and incubate at 37°C for 5 min. Measure the absorbance at 600 nm ^{*2} of the test sample and standard solution with the blank solution as the control. | | |

^{*2} In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

(2) Assay in a Test Tube


Perform the assay in a test tube according to the following table scheme.

| | Test | Standard | Blank |
|-------------------|---|---------------------------------|--------|
| Chromogen reagent | 3.0 mL | 3.0 mL | 3.0 mL |
| Sample | Serum 20 μ L | Standard solution 20 μ L | – |
| | Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm. Spectrophotometer : 600 nm ^{*3} | | |

^{*3} In two wavelength assay, measure using Main wavelength 600 nm/Sub wavelength 700 nm.

9. Standard Curve

[assay in a microplate]

10. Performance

(1) Sensitivity [assay in a Test Tube]

- The absorbance is below 0.10, when measuring purified water as a sample.
- The absorbance is 0.09 to 0.25, when measuring 300 mg/dL triglycerides as a sample.

(2) Specificity

The triglycerides concentration is less than $\pm 12\%$, when measuring the known concentration of control serum as a sample.

11. Usage Notes

(1) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purpose other than those described herein.
- The vial is sealed at reduced pressure. Slowly remove the stopper in order not to blow the powder in the vial.
- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction induced by Chromogen Reagent is completed in about 5 min at 37°C. The reaction is also completed within 15 min at room temperature of over 20°C.
- When the incubation is continued at 37°C, the absorbance gradually decreases. Therefore, the incubation should be done within 20 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This product calculates the concentration of triolein, which exists most within an organism. The concentration of glycerol is 31.2 mg/ dL (corresponding to 300 mg/dL triolein).
- This kit is for research use only. Not for diagnostic use.

(2) Safety Precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(3) Waste

- The waste should be processed appropriately according to regulations.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.
- Triglyceride Standard Solution contains 0.05% sodium azide. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagents.

[Reference]

1. Spayd, R. W., Bruschi, B., *et al.*: *Clin Chem.*, 24, 1343 (1978).

LabAssay™ Triglyceride

| | |
|-------------------|------------------------|
| [Storage] | Store at 2°C - 10°C |
| [Expiration date] | Indicated on the label |
| [Package] | For 350 tests |
| [Cat #] | 291-94501 |

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