



# FUJIFILM Cellular Dynamics

FUJIFILM Cellular Dynamics, Inc., Madison, WI USA

**Background:** The FDA Modernization Act 2.0 underscores the urgent need for non-animal-based new approach methodologies (NAMs) for chemical safety assessment. Human induced pluripotent stem cell-derived cardiomyocytes (iPSC-CMs) are emerging as critical tools in in vitro cardiotoxicity studies, with iPSC-derived cardiac fibroblasts and endothelial cells enabling the development of advanced co-culture and 3D cardiac models. However, serum-free, off-the-shelf media that support both acute and chronic toxicity testing across a broad pharmacological range remain scarce.

**Results:** MEA results demonstrated consistent field potential durations (FPDs), with long-term culture in iCSFM showing slight FPD shortening, while enhancing assay reliability for albumin-binding compounds. A panel of hERG blockers (E-4031, dofetilide, ondansetron) showed dose-dependent effects on cardiac action potential morphology. Metabolic assays revealed decreased maximal respiration with cardiotoxic molecules, while cardiac contractility decreased significantly with chronic idarubicin exposure and increased contraction amplitudes with (S)-BayK-8644. The incorporation of 3D co-culture models showed improved biological relevance, capturing positive inotropic responses and immune response toxicology.

## Media, Methods, and Assay Workflows




**iCell Cardiomyocytes<sup>2</sup> in iCMM on 0.1% Gelatin w/ CalBryte 520 [2  $\mu$ M]**




**iCell Cardiomyocytes<sup>2</sup> in iCSFM on 0.1% Gelatin w/ CalBryte 520 [2  $\mu$ M]**

**iCell Cardiomyocytes<sup>2</sup> in iCTAM on 0.1% Gelatin w/ CalBryte 520 [2  $\mu$ M]**

## Serum-free Medium for Acute and Chronic Toxicity Assessments



Condition	Amplitude (ΔF/F₀ dye)
ICM dye	~2900
ICM	~1050
ICSM	~1050
C-Tox	~500



## Cardiac Co-culture: 3D Models and Culture Media

Lot #	Heart Rate (BPM)
A	~51.5
B	~49.5
C	~49.0
D	~50.0
E	~50.5
F	~50.5
G	~49.0
Whole	~50.0

↓ iCell Macrophages 2.0  
↓ + Bay K 8641 [100 nM]

**No** **No**

**Yes** **No** **✓**

**✓** **✓**

0.500 s

- A. iCell CardioSpheres can be generated combining the 3 cells times in the specified ratio using complete co-culture medium (iCMM + iCell Cardiac Co-culture Supplements Kit components) following this protocol.
- B. Uniform sphere formation from 5-10K cells per well in various ultra-low attachment (ULA) plate types.
- C. Spheroid diameter varies depending on number of cells and ULA plate type, but ranges from 350-450  $\mu\text{m}$ .
- D. Preparation of cardiospheres from different lots and combinations of cells show similar functionality. In this example, the Beats Per Min (BPM) of cells in 96w plate on DIV 14 averaged  $\sim 50$  bpm.
- E. Dosing with compounds like isoproterenol results in a positive inotropic response ( $\uparrow$  peak amplitude).
- F. iCell CardioSpheres typically are cultured in iCMM + iCell Cardiac Co-Culture Supplements, but we have demonstrated that supplemented iCSFM also results in functional 3D spheroids. Data is from DIV 14.
- G. Increasing the model complexity via addition of iCell Macrophages 2.0 not only cleans up cellular debris, but also the cells modify baseline and dosed  $\text{Ca}^{2+}$  handling properties.

iCell Cardiomyocytes<sup>2</sup> provide an in vitro test system that recapitulates the metabolism and physiology of native human cardiomyocytes. Complementary cell types including iCell Cardiac Fibroblasts and iCell Endothelial Cells are essential for making more complex and biologically relevant cell models. The work presented here highlights the utility and flexibility of using human iPSC-derived cell types in 3D as a promising in vitro model for measuring compound effects on human heart tissues in high throughput format for drug discovery studies.