

# Development of cardiac contractile function evaluation system by measuring impedance of human iPS cell-derived cardiomyocytes using a multipoint electrode array

多点電極アレイを用いたヒトiPS細胞由来心筋細胞のインピーダンス測定による心収縮機能評価系の開発

OReiko Hara<sup>1</sup>, Hayato Miyoshi<sup>1</sup> 1. FUJIFILM Corporation

○原 玲子<sup>1</sup>, 三好 隼人<sup>1</sup> 1. 富士フイルム株式会社

## Background

### The importance of cardiac contraction evaluation in drug development

- Side effects on the cardiovascular system are one of the main reasons for terminating drug development. Therefore, it is crucial to accurately understand the effects of compounds on the cardiovascular system at an early stage.
- Potential side effects of drugs on cardiac systolic function may increase the risk of adverse cardiovascular events or heart failure.
- However, currently, there is no established in vitro evaluation methods for cardiac systolic function that can be utilized during the early stages of drug development.

### Trends in cardiac contraction evaluation using iPS-derived cardiomyocytes

- There are multiple evaluation methods available for assessing cardiac contractile function using human iPS cell-derived cardiomyocytes (iPSC-CMs) that involve both maturation techniques and detection methods. In this study, we opted for electrical stimulation due to its high throughput and practicality.
- We previously presented the findings of electrical stimulation and impedance measurements at JSOT2020. In this study, we aimed to validate the reproducibility of those results, expand the scope of drug evaluation data, and provide additional biological data such as calcium concentration measurements.

## Purpose of research

We examined whether it is possible to evaluate the effect of drugs on cardiac contractile function by measuring the impedance of iPSC-CMs using a multi-electrode array (MEA) with electrical stimulation.

## Materials and Methods


### Measuring equipment

- ✓ MEA : Maestro Pro (Axion BioSystems)
- ✓ MEA plate : Cytoview (Axion BioSystems)


### Reagents

- ✓ Cell : iCell® Cardiomyocytes<sup>2</sup> (FUJIFILM Cellular Dynamics, Inc.)
- ✓ Coating agent : Fibronectin
- ✓ Media : iCell® CM Plating Medium , iCell® CM Maintenance Medium

### Protocol




### Analysis method




### Test compounds and Mechanisms

| Compounds           | Inotropy | Mechanisms                                       |
|---------------------|----------|--|
| Isoproterenol       | positive | β adrenergic receptor agonist                    |
| Digoxin             | positive | Na <sup>+</sup> /K <sup>+</sup> ATPase inhibitor |
| Foskolin            | positive | Adenylyl cyclase activator                       |
| Omeacantiv mecarbil | positive | Myosin activator                                 |
| Propranolol         | negative | β adrenergic receptor antagonist                 |
| Mexiletine          | negative | Na <sup>+</sup> channel blocker                  |
| (S)-(-)-Bay K 8644  | positive | L-type Ca <sup>2+</sup> channel activator        |
| Pimobendan          | positive | Calcium sensitizer with PDE-3                    |
| Ryanodine           | negative | Ryanodine receptor inhibitor                     |
| Milrinone           | positive | PDE-3 inhibitor                                  |
| Verapamil           | negative | L-type Ca <sup>2+</sup> channel blocker          |



## Result


### Study1 : FFR (Force Frequency Relationship)



Electrical stimulation +  
Electrical stimulation -


Electrical stimulation improves FFR.

### Study3 : Change of Contractility by compounds Comparison with and without electrical stimulation



It was confirmed that the response increased when electrical stimulation was applied for five of these compounds.


### Study2 : Conduction velocity



Electrical stimulation +  
Electrical stimulation -

Electrical stimulation improves  
Conduction velocity.

### Study4 : Change of Contractility by compounds with electrical stimulation comparison to negative control DMSO




The responses of 11 compounds with different mechanisms of action on the contractility of iPSC-CMs were evaluated.

One-way ANOVA Dunnett's test was employed, because all compounds were compared to the same DMSO background data.

### Study5 : Effect of electrical stimulation on biology of iPSC-CMs

#### 1. Gene expression

RNA extraction : RNeasy Mini Kit (QIAGEN)  
qPCR : QuantiTect SYBR Green RT-PCR Kit (QIAGEN)  
Primers : Obtained from TAKARA BIO  
Measurements : CFX384 (Bio-Rad)



- Electrical stimulation did not alter the gene expressions of iPSC-CMs.
- The expression levels of ATP2A2 and PDE3B were found to be comparable to those of the human heart.
- The expression levels of TNNI3, RYR2, and PDE3A were lower than those of the human heart.

#### 2. Ca<sup>2+</sup> concentration


Calcium staining : EarlyTox Cardiotoxicity Kit (Molecular Devices)  
Measurements : FDSS $\mu$ CELL (Hamamatsu Photonics)



- The effect of electrical stimulation on the amplitude of intracellular calcium concentration was evaluated.
- Electrical stimulation did not increase the intensity of the intracellular calcium peak.

#### 3. cAMP concentration

The cAMP concentration was calculated by measuring the luminescence intensity using cAMP-Glo (Promega).



- Electrical stimulation led to an increase in intracellular cAMP concentration.
- Exposure to isoproterenol was found to significantly elevate the cAMP concentration compared to DMSO.
- The intracellular concentration of cAMP remained unaffected by propranolol.

## Discussion

We confirmed that electrical stimulation significantly enhances the maturation of iPSC-CMs, improving factors such as FFR, conduction velocity, cAMP concentration, and drug responses. This method effectively detected the effects of compounds with diverse pharmacological actions on cardiomyocyte contractility. However, we observed that certain compounds had less impact on iPSC-CMs than anticipated. We hypothesize that low gene expression levels, inadequate intracellular calcium concentrations, and insufficient cAMP concentrations could contribute to the diminished compound responses.

Further research is necessary to promote the maturation of iPSC-CMs concerning morphology (myofibril alignment, sarcomere structure, T-tubule organization, etc.), electrophysiology and calcium handling. These improvements will enable a more comprehensive assessment of drugs with a wider range of actions using this evaluation system.

## Conclusion

The MEA assay of iPSC-CMs with electrical stimulation is a promising method for promoting the maturation of iPSC-CMs and detecting the effects of various compounds on cardiomyocyte contractility.

**FUJIFILM**  
Value from Innovation

