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Cell Wash and Preservation Solution

**Reagent for
Research Use**

Cellstor-W

Instructions for Use

Note:

- (1) This product is a reagent for research use and cannot be administered to humans.
- (2) This product is not intended for medical use in humans or animals.

Cellstor-W is a cell wash and preservation solution containing trehalose hydrate for refrigerated and room-temperature storage. More than 90% of the cells were confirmed to be viable when human adipose-derived mesenchymal stromal cells were suspended using the product and stored for 24 hours (storage conditions: cell concentration, 5×10^5 cell/mL; storage temperature, 5°C or 25°C).¹⁾ Unlike Cellstor-S, which contains dextran 40 and inhibits cells settling, this product does not contain dextran 40 and is suitable for washing when performing centrifugation. When needing to store cells in a state of suspension for about 1 hour, using "Cellstor-S" containing dextran 40 is recommended.

Composition and Description

1. Composition

The ingredients and their quantities in a bag (250 mL) are as follows:

Ingredients	250 mL
Trehalose Hydrate	8.29 g
Calcium Chloride Hydrate	0.05 g
Potassium Chloride	0.075 g
Sodium Chloride	1.5 g
Sodium L-Lactate	0.775 g
pH Adjuster	Appropriate Amount
Water for Injection	Appropriate Amount

Electrolyte concentration (mEq/L)

Na ⁺	K ⁺	Ca ²⁺	Cl ⁻	L-Lactate ⁻
130	4	3	109	28

2. Description

This product is a clear, colorless to pale yellow, sterile aqueous solution. The pH and osmotic pressure ratio are as follows:

pH	6.5-8.5
Osmotic Pressure Ratio (relative to normal saline)	1.0-1.6

Usage

1. Remove the protective seal from the rubber stopper.
2. Disinfect the outer surface of the bag by spraying with ethanol disinfectant.
 - The sterility of the outer surface of the product bag (including the surface of the rubber stopper inside the protective seal) is not guaranteed.Disinfect all outer surfaces with ethanol when transferring the product to the cleanliness control area and sterile working area, etc. Do not sterilize by autoclaving, irradiation with γ -rays or electron beams, EOG, or hydrogen peroxide, as these may affect the product.
 - The product name is printed on the bag. In order to facilitate identification, use the attached label sticker as necessary.

3. To extract the solution, carefully insert a needle vertically into the center of the circle indicated on the rubber stopper, and withdraw the solution.
 - If the stopper is not punctured vertically, it may cause contamination of the sample with minute fragments of the rubber stopper, or leakage of the contents. Do not repeatedly puncture the rubber stopper with a needle at the same point.

Precautions for Use

1. Perform all procedures aseptically.
2. Open the outer bag immediately before use.
3. The duration of storage varies depending on the type and state of the cells and the storage conditions. Perform a confirmation test on the cells to be used.
4. The product bag is not designed for cell storage. When storing cells, use a container suitable for cell storage.
5. If an oil-based ink pen is used on the surface of the product bag, organic solvents such as xylene may be eluted into the product. Therefore, do not write on the bag with an oil-based ink pen.

Precautions for Handling

1. Do not freeze the product. Store the product at room temperature, as freezing it may damage the bag containing the product.
2. Avoid strong physical shock or contact with sharp objects as this may cause leakage.
3. Do not use the product in the following instances:
 - If the bag is damaged
 - If foreign matter is observed in the bag
4. When discarding the product, follow the disposal methods prescribed by the relevant local authority.

Storage and Expiration Date

Storage: Store at room temperature.

Expiration date: 3 years

Packaging

250 mL \times 1 bag (soft bag)

*[References]

- 1) Fujita Y, et al. Regen Ther. 2020;14:95-102.

Use Example (Storage method of human adipose-derived mesenchymal stromal cells)¹⁾

- 1) Centrifuge the cells and remove the supernatant. (Wash the cells with Cellstor-W as necessary.)
- 2) Add the product and resuspend the cells.
- 3) Refrigerate or store the sample at room temperature.

Manufactured by:
Otsuka Pharmaceutical Factory, Inc.
115 Kuguhara, Tateiwa, Muya-cho, Naruto, Tokushima, Japan