# RINGER'S SOLUTION (1/4 STRENGTH)

## DILUENT

### 1 INTENDED USE

Ringer's Solution 1/4 strength is a diluent used to prepare dispersion solutions for the microbiological control of water, milk, dairy and other food products. It is also widely used to prepare serial tenfold dilutions.

The typical composition corresponds to that defined in numerous standards ; many of those listed in the Bibliography section.

### 2 PRINCIPLES

Ringer's Solution 1/4 strength is an isotonic solution that prevents osmotic shock from occurring, when bacteria are removed from their usual environment.

### 3 TYPICAL COMPOSITION

The composition can be adjusted in order to obtain optimal performance.

For 1 liter of solution :

- Sodium chloride .....	2,250 g
- Potassium chloride.....	0,105 g
- Calcium chloride.....	0,120 g
- Sodium hydrogenocarbonate .....	0,050 g

pH of the complete ready-to-use media at 25 °C : 7,0 ± 0,2.

### 4 PREPARATION

- Place one tablet (BR001) in 500 mL of distilled or demineralized water.
- Stir slowly until complete dissolution.
- Dispense in tubes (9 mL) or in vials (90 or 225 mL).
- Sterilize in an autoclave at 121°C for 15 minutes.
- Cool to room temperature.

✓ **Reconstitution :**  
1 tablet in 500 mL

✓ **Sterilization :**  
15 min at 121 °C

### 5 INSTRUCTIONS FOR USE

#### Preparation of dispersion solutions :

- Aseptically add 10 or 25 g of the product to analyze to a tared flask containing 90 or 225 mL of the medium in order to prepare a 1/10 suspension (or a solution).
- Homogenize thoroughly.

#### Preparation of tenfold dilutions :

- Add 1 mL of homogenized sample or stock suspension to a tube containing 9 mL of diluent.
- Homogenize thoroughly.
- Repeat this operation until the desired dilution is obtained.

## 6 QUALITY CONTROL

---

**Dehydrated media** : white tablets.

**Prepared medium** : colorless solution, limpid.

Typical culture response after maintaining diluted for 45-60 minutes at 20-25 °C, followed by inoculation onto Casein Soy Agar :

Microorganisms		Growth
<i>Escherichia coli</i>	WDCM 00013	± 30 % colonies / T <sub>0</sub>
<i>Staphylococcus aureus</i>	WDCM 00034	± 30 % colonies / T <sub>0</sub>
<i>Enterococcus faecium</i>	WDCM 00178	± 30 % colonies / T <sub>0</sub>
<i>Enterococcus hirae</i>	WDCM 00089	± 30 % colonies / T <sub>0</sub>
<i>Enterococcus faecalis</i>	WDCM 00176	± 30 % colonies / T <sub>0</sub>

## 7 STORAGE / SHELF LIFE

---

**Tablets** : 2-30 °C.

The expiration date is indicated on the label.

**Prepared solution in vials or tubes (\*)** : 180 days at 2-25 °C.

(\*) Benchmark value determined under standard preparation conditions, following manufacturer's instructions.

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