

Characterization of human iPS cell-derived intestinal epithelial like cells (F-hiSIECTM) and their utility as pharmacological models

OSeiichi Mochizuki¹, Yuki Imakura¹, Yasushi Morohashi¹, Nao Yamazaki¹, Takahiro Iwao², Tamihide Matsunaga², Kentaro Nakamura¹ 1 Bio Science & Engineering Laboratory. FUJIFILM Corporation

2 Dept. Clinical Pharmacy, Graduate School of Pharmaceutical Sciences. Nagoya City Univ.

0. Introduction

We have developed small intestinal epithelial cells (F-hiSIEC[™]) differentiated from human iPS cells. F-hiSIEC[™] has properties similar to those of the human small intestine and have been shown to be useful as a surrogate in vitro model of the human small intestine for drug absorption, immune and inflammatory evaluation, intestinal toxicity evaluation etc. ¹⁾. However, the prediction of drug metabolism in the gastrointestinal tract (Fg value) deviates from that in humans, and improvement in prediction accuracy is needed. In this study, we examined whether a novel method of F-hiSIEC[™] (New medium & novel culture method) could improve the accuracy of Fg prediction. 1) Imakura et al., Biochem. Biophys. Res. Commun.

Characteristics of F-hiSIEC[™]

- Forms barrier function on trans wells
- Expression of transporters and metabolic enzymes is similar to that of humans
- Multiple cell type markers are expressed, including absorptive epithelial cells, goblet cells, and M cells
- Co-culture with **enterobacteria** and culture of **norovirus** is possible

Purpose of this study



1. Method

Overview of differentiation protocol

Cryopreservation

3. Novel method for drug metabolism evaluation

692 (2024) 149356.



ayu														Day.
End	loderm	al cel	ls	Intest	inal s	tem c	ells		Inte	stinal	epith	elial o	cells	
How	v to us	se F-ł	niSIE	 С ^{тм}	 M	 С: Ме	dium C	hange		Usa	ble te	erm f	or a t	test
DAY	0	1	2	3	4	5	6	7	8	9	10	11	12	13
ex. ′	Seeding	MC		МС		MC		MC		MC		MC		
ex. 2	Thu Seeding	Fri MC	Sat	Sun	Mon MC	Tue	Wed MC	Thu	Fri MC	Sat	Sun	Mon MC	Tue	Wed

F-hiSIEC[™] differentiate into mature small intestinal cells by culturing for 9-11 days after cell seeding.

2. Examples of F-hiSIEC[™] applications

Fatty acid receptor expression and response



The novel method increased the expression of pharmacokinetic-related genes. Barrier function was maintained in the novel method as in the conventional method.

4. Fg estimation of CYP3A4 substrates

vell) 09 vell

Permeation

Fg estimate for Midazolam

It was shown that fatty acid receptor-mediated immune responses can be assessed.



TNF α +**IFN** γ +Acetate **TNF** α +**IFN** γ +**Propionate TNF** α +**IFN** γ +**Butyrate**

On day 8 of cell seeding, inflammatory cytokines (100 ng/mL TNFa, 50 ng/mL IFN γ) were added to the basal side of the cell culture insert and short-chain fatty acids to the apical side. The gene expressions and TEER values were evaluated 24 hours after addition. **p* < 0.05, ***p* < 0.01.



Decreased TEER and increased inflammatory cytokine gene expression due to inflammatory conditions were suppressed by the addition of short-chain fatty acids.

Presence of various intestinal epithelial cells







The estimated Fg values of Midazolam in the novel method ranged from 0.62 to 0.70, closer to human Fg values than in the conventional method.



		Estimated Fg Values ("								
Substrate	In vivo human	Without SF ((SF=1.0)	With S	SF					
	Fg value ³⁾⁴⁾	Conventional	Novel	Conventional (SF=11.7)	Novel (SF=1.6)					
Felodipine	0.35	0.88 (+153%)	0.42 (+19%)	0.40 (+14%)	0.31 (-12%)					
Midazolam	0.55	0.89 (+62%)	0.70 (+27%)	0.41 (-25%)	0.59 (+7%)					
Verapamil	0.73	0.98 (+34%)	0.86 (+18%)	0.78 (+7%)	0.80 (+10%)					
Sildenafil	0.78	Not tested	0.90 (+15%)	-	0.85 (+9%)					
rmspe	-	-	-	0.17	0.09					

The presence of goblet cells (MUC2), endocrine cells (REG4), Paneth cells (LYZ), M cells (GP2), and tufted cells (DCLK1) was demonstrated.



Five days after cell seeding, the cells were cultured on cell culture inserts, and Raji cells (lymphoblast-like cells) were co-cultured on the basal side of the cell culture inserts. On day 10, fluorescent particles (0.2 µm) were added to the apical side and incubated for 2h. The number of fluorescent particles passed through to the basal side was counted by flow cytometry.

Co-culturing with Raji cells suggested that uptake by M cells could be assessed.

Parentheses represent difference between estimated Fg value and *in vivo* human Fg value. These differences and root mean squared percentage error (rmspe) were calculated by the following equation (2) and (3), respectively.

Estimated Fg =
$$\frac{A \text{ to B transport}}{(A \text{ to B transport}) + (Metabolite production}) \times SF} \dots (1)$$

Difference =
$$\frac{(\text{Estimated Fg - In vivo human Fg})}{In vivo human Fg} \dots (2) \quad \text{rmpse} = \sqrt{\frac{\sum (Difference)^2}{N}} \dots (3)$$

3) Verma et al., J Med Chem 53: 1098-1108 (2010) 4) Gertz et al., Curr Drug Metab 9: 785-795 (2011)

The estimated Fg values showed good correlation with the measured values. The SF for absolute value prediction was also smaller for the new method, and the accuracy of the predicted values after correction by SF was higher.

Conclusion : F-hiSIEC[™] has characteristics similar to the human small intestine and is a useful pharmacological model. The novel culture method improved the metabolic function of F-hiSIEC[™] and the accuracy of Fg prediction for CYP3A substrate drugs.