
A novel affinity-based method for the isolation of highly purified extracellular vesicles

Wataru Nakai¹, Takeshi Yoshida^{1,2}, Diego Diez³, Yuji Miyatake^{1,2}, Takahiro Nishibu⁴, Ryo Ukekawa⁴, Naoko Imawaka⁴, Ken Naruse⁴, Yoshifusa Sadamura⁴ & Rikinari Hanayama^{1,2,5} 1 Laboratory of Immune Network, WPI Immunology Frontier Research Center (IFReC), Osaka University, Japan 2 Department of Immunology, Kanazawa University Graduate School of Medical Sciences, Japan 3 Quantitative Immunology Research Unit, WPI Immunology Frontier Research Center (IFReC), Osaka University, Japan 4 Life Science Research Laboratories, Wako Pure Chemical Industries Ltd, Japan 5 PRESTO, Japan Science and Technology Agency (JST), Japan

MagCapture[™] Exosome Isolation Kit PS



Method	Tim4-affinity method (MagCapture™ Exosome Isolation Kit PS)	Ultracentrifu- gation	Polymeric precipitation	Antibody- based affinity purification	
EV's Purity		••			
State of vesicles	Intact	Intact	Intact	Not Intact	s
Operability	Easy and Stable	Easy	Easy and Fast	Easy and Stable	
Recovery amount					c F


			_
Code No.	Product Name	Package Size	
299-77603	MagCapture™ Exosome Isolation Kit PS	2 tests	
293-77601	MagCapture™ Exosome Isolation Kit PS	10 tests	

MagCpature™ Exosome Isolation Kit PS



Sample Type: cell culture supernatant, serum, plasma, urine, etc.

◆ MagCpature[™] Exosome Isolation Kit PS can purify any EVs which expose phosphatidylserine on the outer surface of their lipid bilayer. It has been confirmed that human, mouse, and bovine EVs can be purified by this isolation kit.




sEVs from K562 culture supernatant, human serum, or human urine were purified by the Tim4-affinity method or conventional methods (Ultracentrifugation, Polymeric precipitation). Purified sEVs were analyzed by Silver stain or Western blotting.

Tim4-affinity is an ideal method to isolate sEVs containing exosome with much higher purity than that by using other conventional methods.

Tim4-affinity method is suitable for isolating sEVs from biofluids including serum and urine.

Particle analysis of sEVs purified by the Tim4-affinity methods



sEVs from K562 cells were examined by transmission electron microscope (TEM) and nanoparticle tracking analysis (NTA) using NanoSight

Tim4-affinity

Size Distribution Mean 106 ± 4.1 nm

Microarray analysis of miRNA

miRNA microarray analysis of sEVs from COLO201 cells was performed

The appearance of sEVs isolated by the Tim4-affinity method matched the typical saucer-like shape as previous reported*. The mean size of sEVs purified by the conventional methods was larger than that by the Tim4-affinity method due to increased populations of aggregated or fused sEVs larger than 200nm. Tim4-affinity method could isolate sEVs with higher quality than that by using conventional methods.

*Raposp, G. et al. (1996)

The Uptake of sEV Ultracentrifugation Ultracentrifugation (5μg, 2.8x10¹⁰ particles sEVs) (5μg sEVs) (2.8x10¹⁰ particles sEVs) sEVs from K562 cells were Phase labeled by the fluorescence dye and the uptake of these sEVs by HeLa cells was examined. Hoechst33342 (5μg, 1.8x10¹⁰ particles) 5μg, 2.8x10¹⁰ particles) ug, 1.6x10¹⁰ particles 5μg, 3.4x10¹⁰ particles (5.6μg, 1.8x10¹⁰ particles) (4.1μg, 2.8x10¹⁰ particles) sEVs purified by the Tim4-affinity method were more efficiently

incorporated into the HeLa cells than that by the ultracentrifugation.


Tim4-affinity Tim4-affinity miRNA microarray analysis of sEVs from COLO201 cells revealed high correlation of 993 miRNA profiles between Tim4-UC Polymer 284 affinity method and ultracentrifugation method.

Tim4-affinity method

Conclusion

- Tim4-affinity method can isolate high purity and quality extracellular vesicles from cell culture supernatant and biofluid.
- Tim4-affinity method can purify intact extracellular vesicles.
- sEVs purified by the Tim4-affinity method were efficiently taken up by the recipient cells.
- High correlation of microRNA profiles between Tim4-affinity method and ultracentrifugation was

PS Capture[™] Exosome ELISA Kit



- 1. sEV binds to the Tim4 immobilized on the 96-well plate
- 2. Detect the sEV with sEV surface antigen antibody and POD-labeled Anti-Mouse Antibody

Code No.	Product Name	Package Size
297-79201	PS Capture™ Exosome ELISA Kit (Anti-mouse IgG POD)	96 Reactions

Linearity-of-Dilution Assessments

sEV standard: purified sEVs from K562 or COLO201 Culture supernatnat (sEVs were purified by the MagCapture™ Exosome Isolation Kit PS)

Sample: K562 or COLO201 cell culture supernatant.

Culture Supernatant.

Concentration (E6 particles/ml


by the 3D-Gene (TORAY)


Primary Antibody: Anti-human CD63 Antibody Cell culture sup. of K562 Ratio Factor (×) ng/mL ng/mL 0.125 30 1:800 | 0.00125 | 3.02 3.07 1:400 0.0025 6.15 6.19 99.3 1:200 0.005 12.4 1:100 0.01 25.4 26.1 0.01 0.015 1:50 52.2 Cell culture supernatant of COLO201 Assay Expected y = 1351.7x + 0.113 Ratio Factor (x) ng/mL ng/mL 1:1600 0.000625 0.89 1:800 0.00125 1.82 1.72 1:400 0.0025 3.44 3.52 1:200 0.005 7.04 6.78 0.5 0.002 0.004 0.006 0.008 0.01 0.012 1:100 | 0.01 | 13.6

The linearity-of-dilution was good over the wide range of dilution of cell culture supernatant. Tim4-based ELISA Kit could detect sEVs included in the 0.1µL

Western blotting vs Tim4-based ELISA

Samples: Purified sEVs from COLO201 cell culture supernatant (sEVs were purified by the MagCapture™ Exosome Isolation Kit PS)





The sensitivity of the Tim4-based ELISA was 50 to 1,000 times higher than that of western blotting.