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Technical Report

Total protein amount does not reflect extracellular vesicle amount

- Comparison of EV yields evaluated by tetraspanin ELISA -

Introduction

There are many extracellular vesicle (EV) purification methods, ranging from conventional methods such as standard ultracentrifugation and density-gradient ultracentrifugation to use of kits based on various principles sold by many companies. It is known that the sample amount needed and the purity of EVs obtained differ depending on the purification method 1). It is therefore important to choose an appropriate purification method according to the purpose of the research. Several reports comparing purification methods have been published and can guide the choice of appropriate method.

A point to be noted, however, is that it is difficult to accurately compare methods under the circumstances that EVs recovered and impurity types vary with the purification method. For this reason, according to MISEV2018²⁾, it is important to use multiple, complementary methods for evaluation of the purification method itself and of the EVs obtained.

To compare EV yields, parameters including particle number, total protein amount, and total lipid amount are generally measured. For evaluating the particle number, methods using light scattering, such as Nano Tracking Analysis (NTA), and flow cytometry are known. For evaluating the protein amount, colorimetric methods such as BCA (bicinchoninic acid) assay and fluorescent detection are used. It should be noted that every method may be affected by contaminants in the original sample, resulting in overestimation of the EV amount²⁾. ELISA targeting tetraspanin EV markers has also been proposed as a method for EV quantification. As it is known that the expression level of tetraspanin markers on the surface of EVs varies depending on the culture conditions of the cells used for EV production,³⁾ caution should be exercised when comparing EV amounts prepared from different samples. However, the ELISA is useful as a method for comparing the purification efficiency of EVs purified by different methods from the same source.

This article describes the result of comparing the yield of EVs purified by ultracentrifugation with that of EVs purified using our EV purification kit, MagCapture™ Exosome Isolation Kit PS Ver.2 with sandwich ELISA kits targeting each of three EV markers, CD9, CD63, and CD81. In addition, particle number and total protein amount of the same sample were also measured, and the consistency of the results was verified.

Methods

Purification of EVs

EVs were purified from 1 mL of culture supernatant of bone marrow-derived MSC cultured in DMEM containing 10% EV depleted FBS by ultracentrifugation (UC method) or using MagCapture™ Exosome Isolation Kit PS Ver. 2 (Fujifilm Wako, 290-84103) (PS method). Solutions before and after purification, centrifuge supernatant, and bead flow through fraction were used as samples and their EV yields were compared.

UC method

- 1. Centrifuged at 110,000×g for 70 min.
- 2. PBS Wash × 1
- 3. Centrifuged at $110,000 \times g$ for 70 min.
- Suspended in 100 µL of PBS containing EV-Save™ Extracellular Vesicle Blocking Reagent (Fujifilm Wako, 058-09261).

PS method

- 1. Purified in accordance with the package insert of the kit.
- 2. Extracted twice with 50 μL of PBS/EDTA containing EV-Save $^{\text{TM}}$ to obtain 100 μL of purified product.

Tetraspanin ELISA

Three ELISA kits targeting EV markers were used to measure EV marker signal levels. Measurements were performed in accordance with the protocol attached to the kit. As a reference material, bone marrow MSC-derived EVs for which the number of particles was previously measured by NTA was simultaneously measured, and the EV concentration of each sample was calculated in arbitrary units (AU/mL) based on the standard curve obtained.

ELISA kits

- ① Capture antibody, anti-CD9 Ab; Detection antibody, anti-CD9 Ab CD9-Capture Human Exosome ELISA Kit (Fujifilm Wako, 296-83701)
- Capture antibody, anti-CD63 Ab; Detection antibody, anti-CD63 Ab CD63-Capture Human Exosome ELISA Kit (Fujifilm Wako, 290-83601)
- Capture antibody, anti-CD81 Ab; Detection antibody, anti-CD81 Ab CD81-Capture Human Exosome ELISA Kit (Fujifilm Wako, 292-83801)

NTA

NTA was performed using NanoSight NS300 (Malvern Panalytical).

Measurement of total protein amount

The total protein amount was measured using BCA assay, the highly sensitive method of Protein Assay BCA Reagent Kit (Fujifilm Wako, 297-73101).

Results

The sandwich ELISA targeting each of the EV markers, CD9, CD63, and CD81, detected stronger EV marker signal in the EV solution recovered by the PS method than that recovered by the UC method. Strong EV marker signal was detected in the supernatant obtained by ultracentrifugation, showing that recovery of EVs by the UC method was not as efficient as by the PS method (Figure 1).


Figure 1. Comparison of EV yields evaluated using tetraspanin ELISA

Next, the particle number in EV solution purified by each purification method was measured by NTA. Consistent with the result of ELISA, more particles were found in EV solution purified by the PS method (Table 1). However, the total protein amount measured was higher in the sample purified using the UC method (Table 2). As this result is contradictory to the above results, it is suggested that the protein amount in the sample purified by the UC method may have been higher because contaminant proteins derived from substances other than EVs were included in the total protein measurement.

Table 1. Comparison of particle numbers counted using NTA

	Number of particles
UC method	0.09 × 10 ¹⁰ particles/mL
PS method	0.21 × 10 ¹⁰ particles/mL

Table 2. Comparison of total protein amounts measured using BCA

	Protein amount
UC method	30.45 μg/mL
PS method	11.84 μg/mL

Conclusion

- EV yields were compared after purification using UC and PS, based on EV marker signal level evaluated by ELISA and particle number counted by NTA. The results showed that more EVs were recovered in the sample purified using the PS method.
- As the total protein amount may include contaminant proteins, comparison of the EV
 amount based on total protein amount may yield incorrect results.

References

- 1) Doyle, L. et al.: Cells, 8, 727 (2019).
- 2) Théry, C. et al.: J. EXTRACELL. VESICLES, 7, 1535750 (2018).
- 3) Ma, Y. et al.: Sci. Reports. 11, 13471 (2021)