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CYP450 Induction Assessment: 96 well cultures at a cell density of 0.7 million hepatocytes/mL (56,000 hepatocytes/well) were used in the CYP450 induction assessment. The hepatocytes were cultured as collagen-Matrigel® sandwich for 1 day followed by treatment duration of 48-72 hours for mRNA and 72 hours for activity using known enzyme inducers. Induction in CYP450 activity was assessed by quantifying respective metabolite formation by LC-MS/MS. Gene expression was quantified by RT-PCR. Values reflect mean and standard deviation of triplicate treatments (N=3).

Drug Metabolism Activity

Drug Metabolizing Enzyme	Substrate (μM)	Incubation Time (minutes)	Metabolite Quantified	Activity (pmol/minute/million cells)
CYP1A2	Phenacetin (100)	15	Acetaminophen	54.7
CYP2A6	Coumarin (50)	30	7-Hydroxycoumarin	37.0
CYP2B6	Bupropion (500)	15	Hydroxybupropion	10.8
CYP2C8	Paclitaxel (20)	15	6α-Hydroxypaclitaxel	2.4
CYP2C9	Diclofenac (25)	15	4-Hydroxydiclofenac	92.2
CYP2C19	S-Mephenytoin (250)	30	4-Hydroxymephenytoin	10.0
CYP2D6	Dextromethorphan (15)	15	Dextrorphan	25.7
CYP2E1	Chlorzoxazone (250)	15	6-Hydroxychlorzoxazone	34.7
CYP3A4	Midazolam (20)	10	1-Hydroxymidazolam	44.5
	Testosterone (200)	15	6β-Hydroxytestosterone	470.7
ECOD	7-Ethoxycoumarin (100)	30	7-Hydroxycoumarin	102.6

UGT
Sulfotransferase 7-Hydroxycoumarin (100) 30 7-Hydroxycoumarin glucuronide 395.0
Sulfotransferase 7-Hydroxycoumarin (100) 30 7-Hydroxycoumarin sulfate 21.0


Pravastatin Uptake Transporter Activity

Pravastatin (pmol/10 ⁶ Cells)	Pravastatin with Rifampin (pmol/10 ⁶ Cells)	% Inhibition by Rifampin
0.812	0.31	62 %


Pravastatin Uptake Assessment: 96 well cultures at a cell density of 0.7 million hepatocytes/mL (56,000 hepatocytes/well) were used in the Pravastatin Uptake Assessment. After approximately 6 hours in culture, the hepatocytes were treated with and without 20 μM Rifampin for a pre-incubation time of 30 minutes. Following pre-incubation, 25 μM pravastatin with and without rifampin was incubated for a duration time of 6 minutes. Values reflect the mean of triplicate treatments (N=3). The metabolites were identified and analyzed using LC-MS/MS.

Photomicrographs (100X, Phase Contrast)


Phase Contrast Day 2


Phase Contrast Day 5

Phase Contrast Day 9

Efflux Transporter Assessment

Monolayer Comments: HH1085/HH1085 has high attachment efficiency and a confluency of 90-100 % by 24 hours. This lot exhibits excellent morphology and remains intact for over 9 days in culture.

Efflux Transporter Assessment: The hepatocytes were cultured at a cell density of 0.7 million hepatocytes/mL in a 12-well plate as a collagen-Matrigel® sandwich. On day 5, the hepatocytes were treated with incubation medium containing 5 μM carboxy-2',7' dichlorofluorescein diacetate (CDFDA) and imaged on fluorescein isothiocyanate (FITC) fluorescent filter to assess bile canalicular formation.

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■ FUJIFILM Wako Chemicals U.S.A. Corporation 1600 Bellwood Road, Richmond, VA 23237, USA
■ FUJIFILM Wako Chemicals Europe GmbH Fuggerstraße 12, 41468 Neuss, Germany