# Wako Product Update

# **Protein Research Cell Culture**

Please visit the Wako Online Catalog http://search.wako-chem.com



WakoPURE system

No. 10

**Protein Research** 

2004

# 🗝 Index 🖂

Protein Research	
1.in vitro Protein-Synthesizing System	
WakoPURE system (#299-59501, 295-59503)	
Ni Agarose (#145-07981, 141-07983, 149-07984)	
WakoPURE MF-100K (#237-02231)	
WakoPURE Spin Empty Column (#234-02241)	
2.Simple SDS-PAGE Gel Staining	2
Quick-CBB PLUS (#174-00553, 178-00551)	
Molecular Weight Marker, High Range (#134-14501)	
Molecular Weight Marker, Middle Range (#131-14511)	
3. Proteome Research	3
a. Silver Stain MS Kit	
Silver Stain MS Kit (#299-58901)	
b. Negative Gel Stain MS Kit	4
Negative Gel Stain MS Kit (#293-57701)	
c. High purity Matrix at MALDI-TOFMS Analysis	
α-Cyano-4-hydroxycinnamic Acid [CHCA] (#037-19261)	
Sinapic Acid [SA] (#192-13361)	
2,5-Dihydroxybenzoic Acid [DHB] (#044-29101)	
Lysyl Endopeptidase®, MS Grade (#125-05061)	
Trypsin, from Porcine Pancreas, MS Grade (#202-15951)	

#### Cell Culture

1. Cell Freezing Medium	5
BAMBANKER <sup>™</sup> (#302-14681)	

ALPHA	BETICAL	INDEX
	page	Description
А	1	Ni Agarose
В	6	BAMBANKER™
С	5	CHCA
	2	Quick-CBB PLUS
	5	α-Cyano-4-hydroxycinnamic Acid
D	5	2,5-Dihydroxybenzoic Acid
	5	DHB
L	3, 5	Lysyl Endopeptidase®, MS Grade
м	2	Molecular Weight Marker, High Range
	2	Molecular Weight Marker, Middle Range
Ν	4	Negative Gel Stain MS Kit
	1	Ni Agarose
Q	2	Quick-CBB PLUS
\$	5	SA
	1	Silver Stain II Kit wako
	3	Silver Stain MS Kit
	5	Sinapic Acid
T	5	Trypsin, from Porcine Pancreas, MS Grade
W	1	WakoPURE MF-100K
	1	WakoPURE Spin Empty Column
	1	WakoPURE system

#### INDEX -

Online Catalog : http://search.wako-chem.com

# Protein Research

## 1. in vitro Protein-synthesizing System

Next-Generation in vitro Protein-synthesizing System

#### WakoPURE system

#### Individual Components bring out the Maximum Performance !

WakoPURE system is a proprietary protein synthesis & purification technology, which is an *in vitro* protein synthesizing system "reconstituted" from translation factors expressed in *Escherichia coli*.

It is a novel reconstituted system consisting of around 30 purified enzymes necessary for transcription, translation and energy recycling.

All the components involved in transcription and translation except for ribosomes are tagged with hexahistidine at the N or C terminus. Because all the factors are tagged and known components, synthesized protein can be easily purified by simply removing both the ribosome by ultrafiltration and the histidine-tagged factors and enzymes by metal affinity resin. protein of interest can be synthesized and purified about 3 hours.

#### [Features]

**No Tag Required** 

Because all the transcription/translation factors are tagged, synthesized protein can be easily purified. Native form, without any tags.

Time Saving

All you need is to only 3 steps; synthesis, affinity chromatography and ultrafiltration. It takes only 1 minute for handling and 2 hours for incubation and purification.

- Pure Protein Synthesis System WakoPURE system little contain contaminants such as protease and nuclease because of a novel reconstituted system.
- **Up to 100 Ug protein per 1 mL reaction** DHFR (dihydrofolate reductase) can be synthesized at a yield of 50 μg/mL in 1 hour.



#### Figure

- M, 6 : Molecular Weight Marker

  Negative Control (no DNA template)
  - I . Negative Control (no DINA template
  - 2 : Positive Control (DHFR control plasmid) at Step 1
    3 : Sample obtained at Step 1 was purified by ultrafiltration with WakoPURE MF-100K (Wako Cat. #237-02231)
  - 4 : Sample obtained at Step 1 was treated by Ni-Agarose
  - (Wako Cat.#145-07981) (removal of His-tagged factors) and filtered by WakoPURE Spin Empty Column (Wako Cat. #234-02241).
  - 5 : Sample obtained at Step 1 was treated by Ni-Agarose and purified by ultrafiltration with WakoPURE MF-100K

Result of DHFR after synthesis and various purification. (Apply samples on a 12.5% SDS-PAGE gel and electrophoresed, followed by staining the gel with Silver Stain II Kit Wako (Wako Cat. #291-50301)





#### [Note]

Plasmid DNA and PCR product can be used with WakoPURE system as a template DNA, and the template must contain the following components:

An initiation codon (ATG)/A stop codon (TAG, TGA or TAA)/ A T7 promoter sequence located upstream of the target gene/ Prokaryotic Shine-Dalgarno ribosome binding site (SD sequence) located 10 nt upstream of the initiation codon./ More than 6 nt must have at the upstream of the stop codon for the PCR product. No terminator sequence is necessary./ A T7 terminator located in the downstream of the stop codon for plasmid DNA.

Catalog No.	Description	Kit contents	Package Size
299-59501	Wake DLIDE exetons	Solution A/ Solution B/DHFR control vector	4 reactions
295-59503	WakoPURE system	(positive control)/ Universal Primer	16 reactions

This kit provides the necessary for procedures until protein synthesis at the step 1.

#### Related Products

Catalog No. (Package Size)	Description		
145-07981 (5mL), 141-07983 (10mL), 149-07984 (100mL)	Ni-Agarose		
237-02231 (20 EA)	WakoPURE MF-100K	ultrafiltration membrane	Available
234-02241 (20 EA)	WakoPURE Spin Empty Column	filtration column	3001

## 2. Simple SDS-PAGE gel staining

#### **Quick-CBB PLUS**

for electrophoresis

#### Cat.# 174-00553 250 mL 178-00551 1 L

Quick-CBB PLUS is a bottled protein staining kit, which allows simple and quick staining of proteins bands in polyacrylamide slab gels. In comparison with conventional Quick-CBB, fixing procedure is not required and organic solvents such as methanol and acetic acid are not used. The mixing procedure is now also unnecessary as all the required solutions for staining are contained in a single bottle. There are also other improvements such as coloration not occurring on the background. As in conventional Quick-CBB, destaining procedure is not required.

#### [Features]

FAQ

- · No fixing procedure and no organic solvents are necessary.
- · Simple and quick staining without the coloration on background







# COLUMN TWO IS NOT

#### Q1: When protein bands are detected a few minutes after the start of staining, can the staining procedure be terminated by washing the gel?

- A1: Quick-CBB PLUS usually stains protein bands in 10-20 minutes. The staining can be terminated at that time, but 30-60 minutes is recommended for proper staining.
- Q2: When protein bands can be clearly detected and background staining does not occur without washing after staining, is the washing procedure still required?
- A2 : The washing procedure is not imperative, but a clearer staining can be obtained by washing the gel. When protein bands are light after 60-minute of washing, a clearer staining can be obtained after overnight washing.
- Q3 : About disposal
- A3 : Quick-CBB PLUS does not contain any harmful substances, but is a very dark blue color. Therefore it should be collected in a container prepared for disposal.
- Q4 : As for the sensitivity of the staining
- A4 : Some protein bands visible to approx. 10 ng

Catalog No.	Description	Grade	Package Size	Storage
174-00553	Quick-CBB PLUS	for electrophoresis	250 mL	- RT
178-00551			1 L	

#### Related Products

Catalog No.	Description	Grade	Package Size	Storage
134-14501	Molecular Weight Marker, High Range	fan alaatuurkansia	1 mL	2 ~ 10
131-14511	Molecular Weight Marker, Middle Range	for electrophoresis	1 mL	

Wako Product Update No.10

2

3 4 5 6 7 8 9 10 11 12



Figure : Quick-CBB PLUS Stain

using 5~12% gel Sample: Lane 1, 2, 11 and 12 : Molecular Weight Marker, High Range (Wako Cat.#134-14501); Lane 3,4 and 10 : Molecular Weight Marker, Middle Range (Wako Cat.#131-14511) Lane 5, 6, 7, 8 and 9 : 10, 5, 2.5, 1.25 and 0.6µg of BSA, respectively

# Protein Research

# 3. Proteome Research



(100 ng each)



Figure 3. MALDI-TOF/MS of rabbit phosphorylase

The band was excised and treated with Lysyl Endopeptidase <sup>®</sup> (#125-02543). Following the in-gel digestion and preparation, the sample was analyzed on MALDI-TOF mass spectrometer. (These data were provided by Dr. Y. Wada at Osaka Medical Center, Japan.)

#### [Reference]

1. Shevchenko, A., et al., Anal.Chem., 68, 850 (1996)

Shake 20 min. Wash 1 min. × 2

Wash 1 min. × 3

De-staining mixture Leave 15 min.

The excised band

Mass Spectrometry

Developing Soln. Shake 3 ~ 10 min. Stopper Shake 1 min.

2. Farzin, G., et al., Electrophoresis, 20, 601 (1999)

Catalog No.	Description	Grade	Package Size	Storage
299-58901	Silver Stain MS Kit	for electrophoresis	20 tests	2 ~ 10

#### b. Negative Gel Staing MS Kit Negative Gel Staing MS Kit Cat. # 293-57701

It is known that protein bands separated by SDS/polyacrylamide gel electrophoresis (SDS-PAGE) can be visualized as transparent bands against the background of milky white gel stained by negative gel stain containing a Zn/imidazol reagent. We have improved the method using a new imidazol derivative reagent, which allows a clear and stable image of protein bands on the gel as sensitive as that by silver staining, in as little as 10 minutes. The staining technique is useful to obtain the clear and sensitive resolution pattern of the gel before immunoblotting as well as to excise and purify the band of interest from the gel without significant deterioration of amino acid residues for the subsequent studies of protein such as sequencing and mass analysis of peptide.



[Reference] Fernandez-patron, et al., Anal. Biochem., 224, 263 (1995)

Catalog No.	Description	Grade	Package Size	Storage
293-57701	Negative Gel Stain MS Kit	for electrophoresis	20 tests	Room Temperature

Negative Gel Stain MS Kit

Wako Product Update No.10 -

Online Catalog : http://search.wako-chem.com

# c. High-purity Matrix for MALDI-TOFMS analysis



3. Proteome Research

#### [Features]

- High purity matrix is mixed with a test sample and used for proteome analysis by mass spectrometry (MALDI-TOFMS).
- · High purity matrix gives good mass spec reading because of its high purity by recrystallization

#### Effects of recrystallization of $\alpha$ -Cyano-4-hydroxycinnamic Acid (CHCA)

When comparing with MALDI-TOFMS data for commercial CHCA, it was found that recrystallized CHCA (Wako Catalog No. 037-19261) gives clearer mass spec reading with less background noise. (These data were provided by Dr. Wada Y. at Osaka Medical Center, Japan)



#### [Product List]

Catalog No.	Product	Grade	Package Size	Storage
037-19261	α-Cyano-4-hydroxycinnamic Acid [CHCA]	for proteome research	5 × 50 mg	2~10°C
192-13361	Sinapic Acid [SA]		5 × 50 mg	
044-29101	2,5-Dihydroxybenzoic Acid [DHB]		5 × 50 mg	

#### Related Products

Catalog No.	Product	Grade	Package Size	Storage
125-05061	Lysyl Endopeptidase <sup>®</sup> , Mass Spectrometry Grade	for proteome	5 × 20 μg	20%
202-15951	Trypsin, from Porcine Pancreas, Mass Spectrometry Grade	research	5 × 20 μg	-20 C

## 1. Cell Freezing Medium





Cell Culture

Serum-Free Cell Freezing Medium



Wako Cat. No.

302-14681







Il number/vial 1.0 × 106





monkey B cell line Cell number/vial 1.0 × 106



10 months



-**80°**℃

100 80 60 40 20

BAMBANKER

Ce

Description

man B cell line I number/vial 1.0 × 10 <sup>6</sup>	PC Ce
9 months*	
%	100
%	80'
	601
	00
%	40'
	201
	20
Bambanker Serum Serum-free	0'





**BAMBANKER<sup>™</sup>** is manufactured by **HEE LYMPHOTEC Inc.** 

(Tokyo, Japan)

Listed products are intended for laboratory research use only, but not to be used for drug, food or human use.

Storage

Keep at 2 ~ 10℃

- Please visit our online catalog to search for other products from Wako ; http://search.wako-chem.com
  - This brochure may contain products that cannot be exported to your country due to regulations.

**Package Size** 

120 mL

Bulk quote requests for some products are welcomed. Please contact us.

#### Wako Pure Chemical Industries, Ltd.

http://www.wako-chem.co.jp 1-2, Doshomachi 3-Chome Chuo-Ku, Osaka 540-8605, Japan Tel: 81-6-6203-3741 Fax: 81-6-6201-5964 Online Cat.: http://search.wako-chem.com

#### Wako Chemicals USA, Inc.

http://www.wakousa. com Head Office: 1600 Bellwood Road, Richmond, VA 23237 Toll-Free (U.S. only): 1-877-714-1920 Tel: 1-804-714-1920/ Fax: 1-804-271-7791 Los Angeles Sales Office: 15625 Alton Parkway, Suite D, Irvine, CA 92618 Tel: 1-949-679-1700/ Fax: 1-949-679-1701

#### Wako Chemicals GmbH

04X15IBK

http://www.wako-chemicals.de Nissanstraße 2, D-41468 Neuss, Germany Tel: 49-2131-311-0 Fax: 49-2131-311100

6

BAMBANKER