

LIMULUS AMEBOCYTE LYSATE
ENDOSAFE® ENDOCHROME-K™
U.S. License No. 1197

MULTI-TEST VIAL FOR ENDOTOXIN (PYROGEN) DETECTION

INTENDED USE
Co-lyophilized Limulus amoebocyte lysate (LAL) and a synthetic color producing substrate, which is intended for quantitative detection of endotoxins by kinetic-chromogenic methods.

Co-lyophilized *Listeria amebocyt* lysate (LAL) and a synthetic color producing substrate, which is intended for quantitative detection of endotoxins by kinetic-chromogenic methods.

The LAL test is the most sensitive and specific means available to detect and measure bacterial endotoxin, a fever-producing byproduct of gram-negative bacteria commonly known as pyrogen. The basis of the test is that endotoxin produces an opacity and gelation in LAL that is easily recognized.^{1,5} The simplicity and economy of the LAL Test encourages the testing of in-process solutions and raw materials as well as end-product drugs, devices and biologics.⁶ The USP Bacterial Endotoxins Test <85> provides standard methods for validating the LAL Test as a replacement for the rabbit pyrogen test.⁹

Frederick Bang observed that bacteria caused intravascular coagulation in the American horseshoe crab, *Limulus polyphemus*.¹ In collaboration, Levin and Bang⁵ found that the agent responsible for the clotting phenomena resided in the crab's amoebocytes, or circulating blood cells, and that pyrogen (bacterial endotoxin) triggered the color and gel-forming reaction enzymatically.

Maximum Valid Dilution: USP <85> has listed endotoxin limits of 5 EU/kg for intravenous drugs and 0.2 EU/kg for intrathecal drugs.⁹ Specific limits for compendial items have been adopted.⁹ These limits may be used to determine the extent of dilution that may be used to overcome an interference problem without exceeding the limit endotoxin concentration.⁷ The Maximum Valid Dilution (MVD) is calculated by formulae presented in other pharmacopeia.⁹

Throughout the assay the tube or microplate reader has monitored the increase in absorbance. The reader records the time required for the absorbance to increase significantly over background, usually 0.050 to 0.200 OD units. This time is termed the ONSET TIME. The reader's software automatically produces a log/log correlation of the Onset Time of each standard with its corresponding endotoxin concentration. Standard curve features are then displayed and evaluated to determine if the analysis is valid.

An example of a standard series and recovery of an endotoxin spike of 0.5 EU/mL in sample products is presented below.

REPRESENTATIVE LINEAR ANALYSIS			
SPECIMEN	CSE (EU/mL)	MEAN ONSET TIME (sec)	RECOVERY (EU/mL)
STD 1	50.0	337	
STD 2	5.0	493	
STD 3	0.5	801	
STD 4	0.05	1400	
STD 5	0.005	3025	
AVG. SPL1	----	****	<0.005
AVG. PPC	0.5	822	0.706 (141%)

Fitted regression equation: Log(y) = -0.236*Log(x) + 2.879
Correlation Coefficient: r = -0.991

In this example, the positive product control (PPC) for each sample (SPL) yielded endotoxin recovery consistent with absence of inhibition. The negative control should be significantly lower than the lowest standard endotoxin concentration.

If the absolute value of the correlation coefficient is greater than 0.980, a polynomial regression model can be used to describe the standard curve. The polynomial regression is calculated for the user using Endoscan-V software¹⁰, or equivalent software. For additional information, please refer to the related software manuals and the section describing polynomial regression below.

POLYNOMIAL REGRESSION

A polynomial regression model may be fitted to the standard curve, providing the absolute value of the linear correlation coefficient is greater than, or equal to 0.980. Software packages such as Endoscan-V¹⁰ are capable of fitting the best polynomial function to the data.

NOTE: A polynomial standard curve can not be used for initial qualification assays. The FDA accepts linear regression models for such assays, per the regulatory guidelines.

The order of the polynomial regression is determined as follows:

Endoscan-V Software n-1

Where n = the number of standards utilized.

The example assay (above), analyzed using a fourth degree polynomial (five standards), would yield the following results:

REPRESENTATIVE POLYNOMIAL ANALYSIS			
SPECIMEN	CSE (EU/mL)	MEAN ONSET TIME (sec)	RECOVERY (EU/mL)
STD 1	50.0	337	
STD 2	5.0	493	
STD 3	0.5	801	
STD 4	0.05	1400	
STD 5	0.005	3025	
AVG. SPL1	----	****	<0.005
AVG. PPC	0.5	822	0.445 (89%)

Fitted regression equation: Log(y) = 0.0031*Log(x)⁴ + 0.0001*Log(x)³ + 0.011*Log(x)² -0.216*Log(x) + 2.837

Linear Correlation Coefficient: r = -0.991

LIMITATIONS OF PROCEDURE

Samples may be tested by LAL methods provided that no inhibition or enhancement conditions are present that can not be eliminated by an acceptable dilution (refer to MVD calculation) or sample-pretreatment, such as buffering. If the LAL method cannot be validated at a concentration within the maximum valid dilution, the LAL test cannot be substituted for the USP Bacterial Endotoxins Test.

PERFORMANCE CHARACTERISTICS

Linearity: The linearity of the standard curve within the concentration range used to determine endotoxin levels must be verified. No less than 3 endotoxin standards, spanning the desired concentration range, should be assayed at least in triplicate.¹⁰ The absolute value of the coefficient of correlation, r, shall be greater than or equal to 0.980.⁹

QUALITY CONTROL PROCEDURES FOR THE KINETIC-COLORIMETRIC METHOD

Follow the U.S. Pharmacopeia, 37th revision for end-product testing using Kinetic Colorimetric methods, including the attainment of a Positive Product Control within 50 – 200% of theoretical concentration. Standard curves must have a correlation coefficient ≤ -0.980.

Charles River Endosafe has developed a guide for initial qualification of kinetic incubating microplate readers.

EXPECTED VALUES

Endochrome-K™ LAL Reagent is standardized against the U.S. Reference Standard Endotoxin (RSE). Endotoxin can be quantified if the concentration is within the range of the standard curve. Water and materials derived from biological sources may contain measurable levels of endotoxin if purification efforts are incomplete. Determined endotoxin content should be compared to the endotoxin limit to assess its significance.

Using the appropriate conditions, Endosafe® Endochrome-K™ has an effective range from 100 - 0.001 EU/mL. Factors influencing the selection of the standard curve range include 1) the parameters of the analytical instrument, 2) the choice of regression models, and 3) the quality of supporting analytical reagents and labware.

BIBLIOGRAPHY

1. Bang, F.B. "A Bacterial Disease of Limulus Polyphemus." Bull. Johns Hopkins Hosp. 98, p.325 (1956).
2. Cooper, J.F. and Harbert, J.C. "Endotoxins as a Cause of Aseptic Meningitis after Radionuclide Cisternography." J. Nucl. Med., 16, p.809 (1976).
3. Cooper, J.F., Levin, J., and Wagner, H.N. "Quantitative Comparison of In Vitro and In Vivo Methods for the Detection of Endotoxin." J. Lab. Clin. Med., 78, p.138 (1971).
4. Hochstein, H.D. "The LAL Test versus the Rabbit Pyrogen Test for Endotoxin Detection: Update '87." Pharm. Technol., 11(6), p.124 (1987).
5. Levin, J. and Bang, F.B. "Clottable protein in Limulus: Its Localization and Kinetics of Its Coagulation by Endotoxin." Thromb. Diath. Haemorrh., 19, p.186 (1968).
6. McCullough, K.Z. "Process Control: In-process and Raw Material Testing Using LAL." Pharm. Technol., 12(5) p.40 (1988).
7. Weary, M.E. "Understanding and setting endotoxin limits." J. Parent. Sci. & Tech., 44:1, p. 16 (1990).
8. Cooper, J.F. "Resolving LAL Test Interferences." J. Parent. Sci. & Tech., 44:1, p.13 (1990).
9. Bacterial endotoxins test <85>. In The U.S. Pharmacopeia, 37th rev. United Book Press, Inc., Baltimore, MD
10. Reference Guide for Endoscan-V, Charles River Laboratories, Inc., 1023 Wappoo Road, Suite 43B, Charleston, SC, 29407 USA

Manufactured By: CHARLES RIVER ENDOSAFE
Div. of Charles River Laboratories, Inc.
1023 Wappoo Road, Suite 43B
Charleston, SC 29407, USA
PHONE NUMBER: 843-402-4900
FAX NUMBER: 843-766-7576

PIR17003