Code No. 294-36731 (50 tests)

PYROSTAR™ Neo Package Insert

[Summary and History]

Lipopolysaccharide (LPS) also known as endotoxin, is a component of the outer membrane of Gram-negative bacteria. Since endotoxin causes serious symptoms such as fever and shock in the human body, the testing of bacterial endotoxins in parenteral drugs and medical devices is important.

The bacterial endotoxins test is performed using a lysate reagent that applies the blood coagulation of horseshoe crab. PYROSTAR™ Neo is a chromogenic endotoxin detection reagent in which a recombinant clotting factors, chromogenic substrate, and a buffer component are lyophilized. Horseshoe crab factor C, factor B and clotting enzyme are prepared by recombinant technologies.

Since natural horseshoe crabs are not used as a raw material, you can contribute to the protection of horseshoe crabs by switching from the conventional horseshoe crab-derived lysate reagent to PYROSTAR™ Neo.

[Features]

1. Since the main component is recombinant protein, the difference in reactivity between reagent lots is small.
2. Since (1→3)- β -glucan sensitive factor G is not contained, there is no false positive reaction due to (1→3)- β -glucan.
3. PYROSTAR™ Neo is a reagent for Kinetic Chromogenic Assay by onset time method using a microplate and can measure many samples at the same time in a wide quantitative range.

[Principle]

PYROSTAR™ Neo is a lyophilized reagent that contains three horseshoe crab recombinant proteins, factor C, factor B, and proclotting enzyme, with the chromogenic substrate Boc-Thr-Gly-Arg-pNA. The reaction mechanism of PYROSTAR™ Neo in the presence of endotoxin is summarized in Figure 1. Endotoxin in the test sample activates factor C, and the subsequent cascade reaction releases *p*-nitroaniline from the chromogenic substrate. Endotoxin concentration in the test sample is quantified based on the relationship between the time required for the absorbance to reach a preset threshold value and the endotoxin concentration.


Figure 1 : Reaction Mechanism of PYROSTAR™ Neo

[Kit Contents]

PYROSTAR™ Neo

- A lyophilized reagent in which recombinant *Limulus polyphemus* factor C, factor B and proclotting enzyme, and additives such as chromogenic substrate and buffer are co-lyophilized.

[Test Procedure]

I. Materials & Equipment Required but NOT Supplied

1. Microplate reader and software
Capable of maintaining a temperature at 37°C during kinetic readings, and processing absorbance data by onset time method.
2. 96-well microplates* (BioCleanPlate Wako™ Code : 293-35221)
3. Pipettes
4. Pipette tips* (e.g. BioCleanTip Wako® Extend S II Code : 294-35011, 1000 II Code : 298-35031, 200 II Code : 291-35021)
5. Repeating pipette (e.g. Eppendorf Multipette M4)
6. Repeating pipette tips* (e.g. Eppendorf Combitips advanced Biopur)
7. Test tubes for dilution* (with aluminum caps are recommended)
8. Vortex mixer (e.g. Vial Tube Mixer VTM-252 II Code : 291-36241)
9. Water for Bacterial Endotoxins Test*
10. Standard Endotoxin (e.g. US Pharmacopeia Reference Standard Endotoxin)

* : Glassware such as test tubes should be sterilized by dry heat at 250°C for 30 minutes or more, and water for bacterial endotoxins test should be used. When using plastic wares such as 96-well microplates and pipette tips, make sure there is no endotoxin contamination or interference with the measurement.

II. Preparation of the reagents and test samples

1. Reconstitution of PYROSTAR™ Neo
 - 1) Collect the lyophilized product into the bottom of the vial by tapping on a firm surface. Small amounts of the product on the rubber stopper will not affect testing.
 - 2) Remove the aluminum cap from the vial.
 - 3) Slowly pull up the rubber stopper. The inside of the vial forms a vacuum, therefore avoid dispersion of the product. Do not touch the inside of the vial or the inside part of the rubber stopper when pulling up the rubber stopper.
 - 4) Dispense 2.7 mL of water for bacterial endotoxins test into the vial. Do not allow the rim of the vial to become wet.
 - 5) Gently swirl the vial to dissolve the lyophilized product completely. Be careful not to make bubbles in the solution or agitate the vial vigorously.

Notes : Dissolve the reagent just before use and use it all at once. If the dissolved reagent needs to be stored, freeze it at -30°C or lower and use it within 2 weeks. If refrigerate, keep it at 4°C or lower within 4 hours. Do not freeze and thaw repeatedly.

2. Preparation of Standard Endotoxin Solution

Dissolve the endotoxin standard solution by the method specified for the standard and dilute it with water for bacterial endotoxins test.

3. Dilution of the test samples

Dilute the test samples as needed.

4. pH adjustment of samples

If the pH of a mixture of the sample solution and PYROSTAR™ Neo solution is outside the range of 7.2 to 8.0, adjust the pH of the sample solution with a diluted sodium hydroxide or diluted hydrochloric acid solution.

III. Assay procedure

1. Recommended microplate reader system settings

Temperature : 37°C

Assay mode : Kinetic Onset time

Onset OD : 0.015

Auto mix : Once

Wavelength*1 : 405 nm for main and 650 nm for reference

Reading time*2 : 60 minutes

*1 : Two-wavelength measurement is recommended.
The main wavelength can be in the range of 370 nm to 410 nm, and the reference wavelength can be used in the range of 490 nm to 700 nm. Sensitivity varies slightly depending on the wavelength.

*2 : The reading time can be shortened depending on the required sensitivity.

2. Dispense 50 µL each of water for bacterial endotoxins test (negative control), standard solutions and samples into the corresponding wells of a microplate.

3. Gently swirl the PYROSTAR™ Neo solution again and quickly add 50 µL to each well using a repeating pipette.

4. Immediately after the addition of PYROSTAR™ Neo solution is completed, transfer the microplate into the microplate reader and start measurement. Be sure to agitate the plate at the start of measurement.

IV. Calculation of endotoxin concentration

A calibration curve is created using the onset time obtained from the standard endotoxin, and the endotoxin concentration of the sample is calculated from the onset time of the sample and the calibration curve.

The logarithm of endotoxin concentration can be used as the x-axis of the calibration curve and the logarithm of onset time as the y-axis.


Figure 2 : Example of calibration curve of PYROSTAR™ Neo

【Expected Performance】

When used under the recommended conditions,

1. 0.001 EU/mL endotoxin can be detected within 60 minutes.
2. A linear calibration curve ($|r| \geq 0.980$) can be created in the range of 0.001 to 50 EU/mL.

【Warning and General Precautions】

1. PYROSTAR™ Neo does not correspond to the amebocyte lysate specified in bacterial endotoxins test of the Pharmacopoeia. When used as an alternative method for bacterial endotoxins test of pharmacopoeia, confirm that the accuracy, precision, sensitivity, specificity, etc. are equal or better compared to the bacterial endotoxins test using the amebocyte lysate.
2. PYROSTAR™ Neo is extremely sensitive to endotoxins. Be careful not to cause endotoxin contamination and always use aseptic techniques.
3. Since PYROSTAR™ Neo detects endotoxin on the same principle as the horseshoe crab-derived lysate reagent, the test sample may interfere with the reaction as in the case of the lysate reagent. Please evaluate the interfering factors in the same way as the bacterial endotoxins test.
4. Discard product if the color has changed to yellow or any insoluble matter is found after dissolution.

【Storage】

2 to 10°C

【Expiration Date】

Printed on the label

【Package】

1 multi-test vial (for 2.7 mL, 50 tests)

【References】

1. Bacterial Endotoxins Test and alternative methods using recombinant protein reagents for endotoxin assay, 〈G4-4-180〉 General Information, The Japanese Pharmacopoeia 18th edition, 2021
2. 4.01 Bacterial Endotoxins Test, The Japanese Pharmacopoeia 18th edition, 2021
3. Guidance for Industry, Pyrogen and Endotoxins testing : Questions and Answers, U.S. Department of Health and Human Services, Food and Drug Administration, 2012

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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