[LBIS Rat Insulin ELISA Kit (U-E type)]

Cat #636-05581

Please, read this instruction carefully before use.

This kit is manufactured by FUJIFILM Wako Pure Chemical Corporation. Use only the current version of Instruction Manual enclosed with the kit!

1. Intended use

LBIS Rat Insulin ELISA Kit (U-E type) is a highly sensitive sandwich ELISA system for quantitative measurement of rat insulin. This kit is intended for research use only.

2. Storage and expiration

When the complete kit is stored at 2 °C - 8 °C (Do not freeze), the kit is stable until the expiration date shown on the label on the container. Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

3. Introduction

Insulin is a peptide hormone secreted from B cells of islet of Langerhans in the pancreas with a molecular weight of about 5800 and pI 5.4. It is consisted of 2 chains, A and B. It has 3 disulfide bonds formed between A6 and A11, A7 and B7, and A20 and B19. Insulin exists as a dimer molecule in acidic to neutral solution without Zn ion, and as a hexamer including two Zn ions in neutral solution if Zn ions are present. Main targets of insulin are liver, muscle, and adipose tissue. Insulin actions in these targets are as follows. In the liver, it promotes glycogenesis, protein synthesis, fatty acid synthesis, carbohydrate utilization, and inhibition of gluconeogenesis. In the muscle, it promotes membrane permeability for carbohydrates, amino acids and K ion, glycogenesis, protein synthesis, while inhibits protein degradation. In the adipose tissue, it promotes membrane permeability for glucose and fatty acid synthesis. A precursor of insulin, called proinsulin with a single polypeptide chain, is first synthesized in the cell, then sulfide bonds are formed, and finally by enzymatic cutting at two sites, active insulin and c-peptide (connecting peptide) are formed. Potency of an insulin preparation was originally determined by bioassay. However, whole body bioassay inevitably shows poor precision owing to individual variation.

WHO issued 1st International Standard for human insulin in 1986 which has the potency of 26 IU/mg (0.038 mg/IU). In the same year, 1st International Standard of bovine insulin, the potency of which is 25.7 IU/mg, and Porcine insulin 1st International Standard, 26 IU/mg, were provided. Before these standards, in 1974, 1st International Reference Preparation of human insulin for immunoassay was provided as 3 IU/ampoule. Based on the above data, if the biological activity of insulin per molecule is the same among various animal species, potencies of animal insulin might be calculated from their molecular weights. But, so far, we do not have experimental proof about this. As the molecular weights of insulin of various animals are nearly the same, and the differences are within 1 %, there may be no critical fault if we think that the general potency of insulin is 26 IU/mg.

Rat and mouse have two molecular species of insulin, type 1 and type 2. Amino acid sequences of these molecular species are same between rat and mouse. But as their ratios are different between these two animal species, it is recommended to use standard preparation derived from each animals.

4. Assay principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS Rat Insulin ELISA Kit (U-E type), biotinylated anti-insulin antibody, and standard or sample are incubated in monoclonal anti-insulin coated wells to capture insulin bound with biotinylated anti-insulin antibody. After 2 hours' incubation and washing, HRP (horse radish peroxidase)-conjugated streptavidin is added, and incubated for 30 minutes. After washing, HRP-conjugated streptavidin remaining in wells are reacted with a chromogen (TMB) for 30 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is proportional to insulin concentration. The standard curve is prepared by plotting absorbance against standard insulin concentrations. Insulin concentrations in unknown samples are determined using this standard curve.

5. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.

LBIS Rat Insulin ELISA Kit (U-E-type)

- Use clean laboratory glassware.
- Avoid contact with the acidic stop solution and chromogen (TMB) containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the stop solution because it is sulfuric acid. The stop solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Unused samples and used tips should be rinsed in 1 % formalin, 2 % glutaraldehyde, or more than 0.1 % sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- The materials must not be pipetted by mouth.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20 °C - 25 °C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30 %.

6. Reagents supplied

| Components | State | Amount |
|---|-----------------------------------|------------------|
| (A) Anti-Insulin-coated plate | Use after washing | 96 wells/1 plate |
| (B) Standard Insulin solution (5 ng/mL) | Concentrated. Use after dilution | 500 µL/1 vial |
| (C) Buffer solution | Ready for use. | 60 mL/1 bottle |
| (D) Biotinylated anti-insulin antibody | Concentrated. Use after dilution. | 200 µL/1 vial |
| (E) HRP-conjugated streptavidin | Concentrated. Use after dilution. | 200 µL/1 vial |
| (F) Chromogen (TMB) | Ready for use. | 12 mL/1 bottle |
| (H) Stop solution Be careful! | Ready for use. | 12 mL/1 bottle |
| (I) Wash stock solution (10×) | Concentrated. Use after dilution. | 100 mL/1 bottle |
| Plate seal | — | 3 sheets |
| Instruction Manual | — | 1 copy |

7. Equipments or supplies required but not supplied Use as a check box

Deionized water (or Distilled water). Test tubes for preparation of standard solution series. Glassware for dilution of Wash stock solution (10×) (a graduated cylinder, a bottle) Pipettes (disposable tip type). One should be able to deliver 10 µL precisely, and another for 100 µL. Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 100 µL. Paper towel to remove washing buffer remaining in wells. A vortex-type mixer. A shaker for 96 well-plate (600 rpm - 1200 rpm) An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. A 96 well-plate reader (450 nm ± 10 nm, 620 nm: 600 nm - 650 nm) Software for data analysis.

8. Preparation of reagents

- ◆ Bring all reagents of the kit to room temperature (20 °C - 25 °C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Standard Insulin solution (5 ng/mL)]

Make a serial dilution of master standard (5 ng/mL) solution to prepare each standard solution (39 pg/mL - 2500 pg/mL). *Unit reduction for µIU/mL is 26 IU/mg. (Refer to 3. Introduction.)

LBIS Rat Insulin ELISA Kit (U-E-type)

| Volume of standard solution | (C) Buffer solution | Concentration (pg/mL) | Concentration (μ IU/mL) |
|----------------------------------|---------------------|-----------------------|------------------------------|
| Original solution: 100 μ L | 100 μ L | 2500 | 65 |
| 2500 pg/mL solution: 100 μ L | 100 μ L | 1250 | 32.5 |
| 1250 pg/mL solution: 100 μ L | 100 μ L | 625 | 16.3 |
| 625 pg/mL solution: 100 μ L | 100 μ L | 313 | 8.13 |
| 313 pg/mL solution: 100 μ L | 100 μ L | 156 | 4.06 |
| 156 pg/mL solution: 100 μ L | 100 μ L | 78 | 2.03 |
| 78 pg/mL solution: 100 μ L | 100 μ L | 39 | 1.02 |
| 0 (Blank) | 100 μ L | 0 | 0 |

[(D) Biotinylated anti-insulin antibody]

Prepare working solution by dilution of (D) with the buffer solution (C) to **1:100**.

10 mL of the diluted solution is enough for 96 wells.

[(E) HRP-conjugated streptavidin]

Prepare working solution by dilution of (E) with the buffer solution (C) to **1:100**.

10 mL of the diluted solution is enough for 96 wells.

[(I) Wash stock solution (10 \times)]

Dilute 1 volume of the concentrated Wash stock solution (10 \times) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example: 100 mL of concentrated washing buffer (10 \times) and 900mL of deionized water (or distilled water).

【Storage and stability】

[(A) Anti-Insulin coated plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2 °C - 8 °C. The strip will be stable until expiration date.

[(B) Standard Insulin solution (5 ng/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C) Buffer solution] & [(F) Chromogen (TMB)]

If not opened, store at 2 °C -8 °C. It maintains stability until expiration date. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Biotinylated anti-insulin] & [(E) HRP-conjugated streptavidin]

Unused working solution (already diluted) should be disposed.

[(H) Stop solution]

Close the stopper tightly and store at 2 °C - 8 °C. It maintains stability until expiration date.

[(I) Wash stock solution (10 \times)]

The rest of undiluted buffer: if stored tightly closed at 2 °C - 8 °C, it is stable until expiration date.

Dispose any unused diluted buffer.

9. Technical tips

- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2 °C - 8 °C.
- Time the reaction from the pipetting of the reagent to the first well.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The chromogen (TMB) should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of stop solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the anti-insulin-coated plate is module type of 8 wells \times 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30 %, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step

of incubation.

10. Preparation of samples

This kit is intended to measure insulin in rat serum, plasma (preferably obtained with heparin), culture medium and tissue/cell extracts. The necessary sample volume for the standard procedure is 10 μ L. Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results.

Hemolytic and hyperlipemic serum samples are not suitable.

*** To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 40 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Dilution of a sample should be made in a test tube using buffer solution prior to adding them to wells. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

Storage and stability

Insulin in samples will be inactivated if stored at $2^{\circ}\text{C} - 8^{\circ}\text{C}$. If it is necessary to store sample in refrigerator ($2^{\circ}\text{C} - 8^{\circ}\text{C}$), add aprotinin at final concentration of 100 KIU/mL - 500 KIU/mL. (KIU: kallikrein inhibitor unit). If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C . Avoid repeated freeze-thaw cycles.

11. Assay procedure

Remove the cover sheet of the anti-Insulin-coated plate after bringing up to room temperature.

- (1) Wash the anti-Insulin-coated plate (A) by filling the wells with 300 μ L of washing buffer and discard 4 times(*②), then strike the plate upside-down onto several layers of paper towels to remove residual buffer in the wells.
- (2) Pipette 100 μ L of biotinylated anti-insulin antibody to all wells. Shake the plate gently on a plate shaker (*③).
- (3) Pipette 10 μ L of sample to the designated sample wells.
- (4) Pipette 10 μ L of standard solution to the wells designated for standards.
- (5) Shake the plate gently on a plate shaker(*③).
- (6) Stick a plate seal (*④) on the plate and incubate for 2 hours at $20^{\circ}\text{C} - 25^{\circ}\text{C}$.
- (7) Discard the reaction mixture and rinse wells as step (1).
- (8) Pipette 100 μ L of HRP-conjugated streptavidin to all wells, and shake as step (5).
- (9) Stick a plate seal (*④) on the plate and incubate the plate for 30 minutes at $20^{\circ}\text{C} - 25^{\circ}\text{C}$.
- (10) Discard the reaction mixture, and then wash the plate as step (1).
- (11) Pipette 100 μ L of chromogen (TMB) to wells, and shake as step (5).
- (12) Stick a plate seal (*④) on the plate and incubate the plate for 30 minutes at $20^{\circ}\text{C} - 25^{\circ}\text{C}$.
- (13) Add 100 μ L of the stop solution to all wells and shake as step (5).
- (14) Measure the absorbance of each well at 450 nm (reference wavelength, 620*nm) immediately using a plate reader.


*Refer to the page 7 for notes of *②, *③ and *④.

12. Calculations

- (1) Prepare a standard curve using semi-logarithmic or two-way logarithmic section paper by plotting absorbance* (Y-axis) against insulin concentration (pg/mL) on X-axis. *Absorbance at 450 nm minus absorbance at 620 nm.
- (2) Using the standard curve, read the insulin concentration of a sample at its absorbance*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.

* We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

Physiological or pathological situation of animals should be judged comprehensively taking other examination results into consideration.

Rat Insulin assay standard curve (an example)
Absorbance may change due to assay environment.

13. Performance characteristics

- Assay range

The assay range of the kit is 39 pg/mL - 2500 pg/mL.

If some samples show absorbance more than that of 2500 pg/mL standard, please repeat the assay after proper dilution of samples.

- Specificity

The antibodies used in this kit are specific to insulin. Cross-reactivity of the kit is shown below.

| Substances | Cross-reactivity | Substances | Cross-reactivity |
|----------------|------------------|-----------------|------------------|
| Rat C-peptide | — | Mouse insulin | + |
| Rat proinsulin | + | Mouse C-peptide | — |
| Human insulin | + | | |

*Cross-reactivity at Conc. 2.5 ng/mL.

- Precision of assay

Within assay variation (3 samples, 8 replicates assay) Mean CV was within 10 %.

- Reproducibility

Between assay variation (3 samples, 4 days, assayed in triplicates) Mean CV was within 10 %.

- Recovery test

Standard insulin was added in 3 concentrations to 2 serum samples and were assayed.

The recoveries were 94.1 % - 101 %.

- Dilution test

Two serum samples were serially diluted by 3 steps.

The dilution curves showed excellent linearity with $R^2 > 0.992$.

14. Reference assay data

Insulin assay data in fasted rats, assayed in duplicates.

Mean assay value: 245 pg/mL Standard deviation: 72 pg/mL

Rat strain: SD, Number of animals: 8 males, 6 weeks of age.

After fasting for 24 hours. Serum was used.

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for insulin levels independently.

15. Trouble shooting

- Low absorbance in all wells

Possible explanations:

1) The standard or samples might not be added.

2) Reagents necessary for coloration such as biotinylated anti-insulin antibody, HRP-conjugated

LBIS Rat Insulin ELISA Kit (U-E-type)

streptavidin, or chromogen (TMB) might not be added.

3) Wrong reagents related to coloration might have been added. Wrong dilution of biotinylated anti-insulin antibody or HRP-conjugated streptavidin.

4) Contamination of peroxidase enzyme inhibitor(s).

5) Influence of the temperature under which the kits had been stored.

6) Excessive hard washing of the well plate.

7) Addition of chromogen (TMB) soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

● Intense coloration in all wells including blank

Possible explanations:

1) Improper or inadequate washing. Change washing frequency from 4 times to 5-8 times at the constant stroke after the reaction with HRP-conjugated streptavidin.)

2) Overdeveloping. Incubation time with chromogen (TMB) should be decreased before addition of stop solution.

3) Too high incubation temperature. Adjust the temperature to 20 °C - 25 °C.

● High coefficient of variation (CV)

Possible explanation:

1) Improper or inadequate washing.

2) Improper mixing of standard or samples.

3) Pipetting at irregular intervals.

● Q-1: Can I divide the plate to use it for the other testing?

A-1: Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon

● Q-2: I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2: When we manufacture 96 well-plate, we insert preservation stabilizer in wells.


Summary of assay procedure : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

Bring the well-plate and **all reagents to 20 °C - 25 °C for 2 hours.**

Wash stock solution (10×) concentrate must be diluted to 10 times by deionized water (or distilled water) that returned to 20 °C - 25 °C.

Standard Insulin solution dilution example:

Biotinylated anti-insulin antibody(D) : Dilute to **100 times** by using (C) Buffer solution and use.

LBIS Rat Insulin ELISA Kit (U-E-type)

| | | | |
|--------------------------|--|--------|----|
| <input type="checkbox"/> | Anti-Insulin coated plate | | |
| <input type="checkbox"/> | ↓Washing 4 times (*②) | | *⑥ |
| <input type="checkbox"/> | Biotinylated anti-insulin antibody | 100 μL | |
| <input type="checkbox"/> | ↓Shaking (*③) | | |
| <input type="checkbox"/> | Samples/Standards | 10 μL | |
| <input type="checkbox"/> | ↓Shaking (*③), Incubation for 2 hours at room temp. (Standing(*④)) | | |
| <input type="checkbox"/> | HRP-conjugated streptavidin(E) | | |
| <input type="checkbox"/> | Dilute to 100 times by using (C) Buffer solution and use. (Dilute reagents during the first reaction.) | | |
| <input type="checkbox"/> | ↓Washing 4 times (*②) | | *⑥ |
| <input type="checkbox"/> | HRP-conjugated streptavidin | 100 μL | |
| <input type="checkbox"/> | ↓Shaking (*③), Incubation for 30 minutes at room temp. (Standing(*④)) | | |
| <input type="checkbox"/> | ↓Washing 4 times(*②) | | *⑥ |
| <input type="checkbox"/> | Chromogen (TMB) (After dispense, the color turns to blue depending on the concentration.) | 100 μL | |
| <input type="checkbox"/> | ↓Shaking (*③), Incubation for 30 minutes at room temp. (Standing(*④)) | | |
| <input type="checkbox"/> | Stop solution (After dispense, the color turns to yellow depending on the concentration.) | 100 μL | |
| <input type="checkbox"/> | ↓Shaking (*③)(Immediately shake.) | | |
| <input type="checkbox"/> | Measurement of absorbance (450 nm, Ref 620 nm) immediately (Ref. wave cancels the dirt in the back of plate.) | | *⑤ |

*② After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume: 300μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with HRP-conjugated streptavidin. Standard of plate-washing pressure: 5 mL/min -25 mL/min. (Adjust it depending on the nozzle's diameter.)

*③ Guideline of shaking: 600 rpm - 1200 rpm for 10 seconds × 3 times.

*④ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

*⑤ 600 nm - 650 nm can be used as reference wavelength.


*⑥ After removal of wash buffer, immediately dispense the next reagent.

Worksheet example

| | Strip 1&2 | Strip 3&4 | Strip 5&6 | Strip 7&8 | Strip 9&10 | Strip 11&12 |
|---|------------|-----------|-----------|-----------|------------|-------------|
| A | 2500 pg/mL | Sample 1 | Sample 9 | Sample 17 | Sample 25 | Sample 33 |
| B | 1250 pg/mL | Sample 2 | Sample 10 | Sample 18 | Sample 26 | Sample 34 |
| C | 625 pg/mL | Sample 3 | Sample 11 | Sample 19 | Sample 27 | Sample 35 |
| D | 313 pg/mL | Sample 4 | Sample 12 | Sample 20 | Sample 28 | Sample 36 |
| E | 156 pg/mL | Sample 5 | Sample 13 | Sample 21 | Sample 29 | Sample 37 |
| F | 78 pg/mL | Sample 6 | Sample 14 | Sample 22 | Sample 30 | Sample 38 |
| G | 39 pg/mL | Sample 7 | Sample 15 | Sample 23 | Sample 31 | Sample 39 |
| H | 0 (Blank) | Sample 8 | Sample 16 | Sample 24 | Sample 32 | Sample 40 |

LBIS Rat Insulin ELISA Kit (U-E-type)

Assay worksheet

LBIS Rat Insulin ELISA Kit (U-E-type)

[Storage condition] Store the kit at 2 °C - 8 °C (Do not freeze).

[Term of validity] Expiration date is indicated on the container.

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