**FUJIFILM**

**Wako**

Code No. 178-00551 (1 L)  
174-00553 (250 mL)

## Quick-CBB PLUS

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Quick-CBB PLUS stains protein bands on polyacrylamide gel with CBB G-250. This staining method is simple and fast, and does not require a fixing process. Therefore waste of reagents such as methanol and acetic acid is reduced. Additionally, the protein bands are clearly visible in contrast with low background.

### 【Features】

1. Quick-CBB PLUS does not use organic solvents such as methanol or acetic acid.
2. The protocol is simple and fast.
3. Stains protein bands clearly with low background.

### 【Procedure】

1. Rinse the SDS-PAGE gel with deionized water for 5 minutes in a staining tray. Repeat 2 more times. Drain off the deionized water from the tray.
2. Add 25 mL of Quick-CBB PLUS to the tray and shake it for 30-60 minutes until the protein bands are visible.
3. Drain off the solution and rinse the gel with distilled water for 60 minutes.

### 【Usage note】

1. Use a tray with a smooth inner surface.
2. Immerse the gel completely in staining solution while staining.

### 【Cautions for handling】

1. If the Quick-CBB PLUS appear to cause precipitate due to low temperature, warm it to resolubilize before use.
2. Quick-CBB PLUS does not contain any toxic substances. However, the solution is blue, so disposal in a specific vessel may be necessary.

### 【Reference】

Blakesley, R., W. and Boezi, A., J.: *Anal. Biochem.*, **82**, 580 (1977)

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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