For Electrophoresis Quick-CBB

CBB(Coomassie Brilliant Blue) R-250 has been widely employed as a stain for detection of proteins after resolution with PAGE. It is typically necessary, however, to destain the gel, which is an inconvenient and time-consuming process.

Quick-CBB resolves this problem, allowing for the rapid staining of protein bands within an hour. Background staining does not occur, eliminating the need for destaining.

Features

1. Involves simple preparation and simple procedures.
2. After electrophoresis, staining is completed within 50 minutes.

Reagents

One kit includes :

Staining Solution A (contains CBB R-250)	1L×1
Staining Solution B	1L×1

Preparation

Volumes of reagents below are calculated based on a 75mm×75mm×1.0mm sized polyacrylamide slab gel.

1. **Fixing solution** — Prepare 200 mL of solution by mixing 100 mL of methanol and 20 mL of acetic acid with 80 mL of deionized water.
2. **Staining solution**—Prepare 60 mL of solution by mixing 30 mL of Staining Solution A with 30 mL of Staining Solution B.

[Notes]

1. The volume of reagents used is based on the size of polyacrylamide slab gel and on the size of the tray used for staining. In general, reagents prepared are 5 to 10 times the volume of the gel. Ensure that the gel is soaked completely in solution throughout the procedure.
2. Prepare solution just prior to staining the gel.
3. Use superior quality solvents (methanol and acetic acid) and use only deionized or distilled water.
4. Vessels used should be cleaned thoroughly prior to use.

Staining Procedure

1. **Fixing:** Soak gel in a tray containing 100 mL of fixing solution and shake well for 10 minutes. Repeat this step two times.
2. **Staining:** Remove fixing solution and pour 60 mL of staining solution into tray. Shake for 30 minutes.
3. **Rinsing:** Transfer gel into another tray and wash the stained gel with deionized water to remove the excess stain.

[Notes]

1. The volume of reagents used is based on the size of polyacrylamide slab gel and on the size of the tray used for staining. In general, reagents prepared are 5 to 10 times the volume of the gel. Ensure that the gel is soaked completely in solution throughout the procedure.
2. Prepare solution just prior to staining the gel.
3. Use superior quality solvents (methanol and acetic acid) and use only deionized or distilled water.
4. Vessels used should be cleaned thoroughly prior to use.

Technical notes

1. **Next day staining** — Store electrophoresed gel in fixing solution until staining the following day.
2. **Confirmation of staining** — It is difficult to assess the degree of staining during the procedure due to the dark color of the stain. For easier viewing, the gel may be inclined and pressed every 10 minutes using disposable gloves.
3. **Destaining background** — Destaining is not necessary, but if one chooses to do so, the gel may be shaken for 10 minutes in 100 mL of 7% acetic acid.
4. **Reusing staining solution**— Stain may be reused after adding additional staining solution (approximately 10 mL), though it is not recommended due to decreased staining sensitivity.
5. **Disposal** —
 - a) Fixing solution contains methanol and should be disposed of in a proper container after collecting.
 - b) Staining solution does not contain any hazardous material but should be collected in a disposable container to avoid unwanted staining or spillage.
6. **Storage of gel** — After rinsing the stained gel with deionized water, shake for 5 minutes in 5% glycerin solution. Dry with a gel dryer and seal when storing. If a gel dryer is not used, seal in a plastic bag.
7. **Precipitation** — CBB R-250 will precipitate if stored at a low temperature. If this occurs, dissolve in warm water (30 ~ 40°C).

Storage

Store at room temperature.

Reference

Stephano, J.L.,Gould,M.,Rojas-Galicia,L.: Advantages of Picrate Fixation for Staining Polypeptides in Polyacrylamide Gels., *Anal. Biochem.*, **152** : 308–313(1986).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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