Code No. 295-50201 (For 50 tests)

**DNA Extractor<sup>®</sup> Kit**

This product is in compliance with sodium Iodide method and designed for use in extracting residual host cell DNA in biological samples such as vaccines and biopharmaceuticals, and DNA in serum. It offers high recovery of DNA even from samples containing only a very small amount of DNA. The entire extraction step is completed with a single tube for 60-90 minutes. The extracted DNA can be quantified by qPCR or threshold systems from Molecular Devices, LLC.

There are two protocols for DNA extraction using this kit. Protocol #1 is used for DNA extracting from clean samples such as purified proteins or serum. Protocol #2 is used for crude samples such as cell lysate, solution containing high-concentration protein or lipid, and when the extracted DNA is quantified by threshold systems from Molecular Devices, LLC.

Difference between the two protocols are indicated in the table below.

	Sample of interest	Pretreatment step of sample	Number of washing steps	"Washing Solution (A)" in the kit
Protocol #1	Clean or Serum	Not required	1 time	Not required
Protocol #2	Crude or Threshold systems	Optional	2 times	Required

**【Features】**

- In compliance with sodium Iodide method
- High quality and recovery
- Completed in a single tube for 60-90 minutes.

**【Principle of sodium Iodide method for DNA extraction】**

An chaotropic reagent, Sodium Iodide, and an anionic detergent participate in solubilization of proteins and lipids contained in biological samples. After addition of 2-propanol to the mixture, nucleic acids are co-precipitated with glycogen as a carrier, while other components remain soluble in the solution phase.

**【Kit contents】**

Code No.	Components	Size
292-50211	Sodium Iodide Solution	26 mL × 1
299-50221	Sodium N-Lauryl Sarcosinate Solution	1.2 mL × 1
296-50231	Washing Solution (A)	42 mL × 1
293-50241	Washing Solution (B)	40 mL × 2
290-50251	Glycogen Solution	0.1 mL × 1

**【Storage】**

Store at 2-10°C.

**【Additional materials required for 50 tests】****〈Reagents〉**

- 1) 2-Propanol 20 mL
- 2) Sterile distilled water 6 mL

**〈Instruments〉**

- 1) Microcentrifuge
- 2) Vortex mixer
- 3) 50 × 2 mL or 1.5 mL of centrifuge tube
- 4) Micropipette
- 5) Micropipette tips
- 6) Heating block

**【Precautions for Use】**

- The color of "Washing Solution (A)" may be slightly different between product lots. e.g. faint yellow or colorless. It does not affect quality of DNA extraction.
- If precipitate is observed in "Sodium Iodide Solution" or "Sodium N-Lauryl Sarcosinate Solution", warm it or them at 50°C until the precipitate is no longer visible.
- Use DNA-free, DNase-free and sterilized reagents and Instruments.
- DNA concentration cannot be accurately measured by absorbance if sodium Iodide remains in the extracted DNA. Measure the amount of DNA by qPCR or threshold systems.

**【Protocol #1】****〈Preparation of Reagents〉****Solution I-A**

Add 6 mL of sterile distilled water, 1 mL of "Sodium N-Lauryl Sarcosinate Solution" and 65  $\mu$ L of "Glycogen Solution" into the bottle of 26 mL of "Sodium Iodide Solution". Mix it with a vortex mixer.

Note : Solution I-A can be stored at 2-10°C in the dark for approximately 5 months.

**Solution II**

Add 2  $\mu$ L of "Glycogen Solution" into the bottle of 40 mL of "Washing Solution (B)". Mix it with a vortex mixer immediately.

Notes :

- It does not affect quality of DNA extraction if white precipitates appear in Solution II.
- Solution II can be stored at 2-10°C in the dark for approximately a week at least. If preparing a small volume of Solution II, use required amount of "Glycogen Solution" and "Washing Solution (B)". For example, add 1  $\mu$ L of "Glycogen Solution" into 20 mL of "Washing Solution (B)" in a sterilized tube.

**〈DNA Extraction Procedure〉**

1. Dispense 100  $\mu$ L of sample solution to a 2 mL or 1.5 mL blank centrifuge tube.
2. Add 300  $\mu$ L of Solution I-A into the tube and mix it with a vortex mixer.
3. Warm the tube at 60°C for 15 minutes in a heating block.
4. Remove the tube from the heating block. Add 400  $\mu$ L of 2-Propanol

into the tube and mix it with a vortex mixer.

5. Leave the tube at room temperature for 15 minutes.
6. Centrifuge the tube with  $10,000\times g$  at room temperature for 15 minutes.  
A faint white pellet appears in the tube.
7. Remove supernatant from the tube by using a pipette and/or decanting carefully until the liquid volume in the tube is less than approximately  $100\ \mu\text{L}$ .  
Notes :
  - Be careful not to remove the pellet from the tube.
  - Centrifuge the tube with  $10,000\times g$  at room temperature for 15 minutes again if the pellet floats.
  - In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.
8. Add 1 mL of Solution II into the tube and mix it by inverting gently. Ensure that the pellet is detached from the tube wall.
9. Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes.
10. Remove as much supernatant as possible from the tube by using a pipette and/or decanting carefully.  
Notes :
  - Be careful not to remove the pellet from the tube.
  - Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes again if the pellet floats.
  - In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.
11. Dry the pellet by warming it at  $60^\circ\text{C}$  in a heating block, vacuum dryer or natural drying.  
Notes :
  - Too dry pellet is difficult to be dissolved. Stop drying the pellet when it still is wet.
  - Much residual liquid in the tube may cause low qPCR efficiency because Solution II contains ethanol, which can inhibit qPCR.
  - The pellet contains DNA and glycogen.
12. Dissolve the pellet with sterile distilled water or buffer, and proceed with qPCR.

#### **[Protocol #2]**

##### **⟨Preparation of Reagents⟩**

###### 1. Solution I-B

$65\ \mu\text{L}$  of "Glycogen Solution" into the bottle of 26 mL of "Sodium Iodide Solution". Mix it with a vortex mixer.

Note : Solution I-B can be stored at  $2-10^\circ\text{C}$  in the dark for approximately 5 months.

###### 2. Solution II

Add  $2\ \mu\text{L}$  of "Glycogen Solution" into the bottle of 40 mL of "Washing Solution (B)". Mix it with a vortex mixer immediately.

Notes :

- It does not affect quality of DNA extraction if white precipitates appear in Solution II.
- Solution II can be stored at  $2-10^\circ\text{C}$  in the dark for approximately a week at least. If preparing a small volume of Solution II, use required amount of "Glycogen Solution" and "Washing Solution (B)". For example, add  $1\ \mu\text{L}$  of "Glycogen Solution" into 20 mL of "Washing Solution (B)" in a sterilized tube.

##### **⟨Pretreatment Step of Sample⟩**

This step is optional and can be skipped. In case of extracting DNA from samples containing much protein and/or lipid, do the pretreatment according to ⟨Optional Protocol for Preparation of Samples⟩ mentioned later.

##### **⟨DNA Extraction Procedure⟩**

1. Dispense  $400-500\ \mu\text{L}$  of sample solution to a 2 mL blank centrifuge tube.
2. Add  $20\ \mu\text{L}$  of "Sodium N-Lauryl Sarcosinate Solution" into the tube and mix it with a vortex mixer.
3. Add  $500\ \mu\text{L}$  of Solution I-B into the tube and mix it with a vortex mixer.
4. Warm the tube at  $40^\circ\text{C}$  for 15 minutes in a heating block.
5. Remove the tube from the heating block. Add  $900\ \mu\text{L}$  of 2-Propanol into the tube and mix it with a vortex mixer.
6. Leave the tube at room temperature for 15 minutes.
7. Centrifuge the tube with  $10,000\times g$  at room temperature for 15 minutes.  
A faint white pellet appears in the tube.
8. Remove supernatant from the tube by using a pipette and/or decanting carefully until the liquid volume in the tube is less than approximately  $100\ \mu\text{L}$ .  
Notes :
  - Be careful not to remove the pellet from the tube.
  - Centrifuge the tube with  $10,000\times g$  at room temperature for 15 minutes again if the pellet floats.
  - In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.
9. Add  $800\ \mu\text{L}$  of "Washing Solution (A)" into the tube and mix it by inverting gently.  
Ensure that the pellet is detached from the tube wall.
10. Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes.
11. Remove supernatant from the tube by using a pipette and/or decanting carefully until the liquid volume in the tube is less than approximately  $100\ \mu\text{L}$ .  
Notes :
  - Be careful not to remove the pellet from the tube.
  - Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes again if the pellet floats.
  - In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.
12. Add 1.5 mL of Solution II into the tube and mix it by inverting gently.
13. Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes.
14. Remove as much supernatant as possible from the tube by using a pipette and/or decanting carefully.  
Notes :
  - Be careful not to remove the pellet from the tube.
  - Centrifuge the tube with  $10,000\times g$  at room temperature for 5 minutes again if the pellet floats.
  - In case that residual liquid remains on the tube wall, place the tube upside down on a paper filter to remove it.
15. Dry the pellet by warming it at  $60^\circ\text{C}$  in a heating block, vacuum

dryer or natural drying.

Notes :

- Too dry pellet is difficult to be dissolved. Stop drying the pellet when it still is wet.
- Much residual liquid in the tube may cause low qPCR efficiency because Solution II contains ethanol, which can inhibit qPCR.
- The pellet contains DNA and glycogen.

16. Dissolve the pellet with sterile distilled water or buffer and proceed with qPCR.

Or dissolve the pellet with Zero Calibrator in Total DNA Assay Kit from Molecular Devices, LLC. and quantify the extracted DNA by Threshold systems.

#### 〈Optional Protocol for Preparation of Samples〉

■ In case that protein concentration in the sample solution is less than 2 mg/mL.

1. Add the reagents below into the sample solution to the final concentration.

Reagents	Final concentration
Sodium Dodecyl Sulfate	0.10%
Dithiothreitol	50 mmol/L

2. Warm it at 55°C for 30 minutes in a heating block.

3. Proceed with 〈DNA Extraction Procedure〉 in Protocol #2.

■ In case that protein concentration in the sample solution is more than 2 mg/mL.

1. Add the reagents below into the sample solution to the final concentration.

Reagents	Final concentration
Sodium Dodecyl Sulfate	0.10%
Dithiothreitol	50 mmol/L
Sodium Chloride	150 mmol/L
EDTA	1 mmol/L

2. Add 20  $\mu$ g of Proteinase K per 1 mg of protein in the sample solution in to the sample solution.

3. Adjust pH to be approximately 7.5 with Tris.

4. Warm it at 55°C in a heating block for an hour at least.

5. Proceed with 〈DNA Extraction Procedure〉 in Protocol #2.

#### 【FAQ】

Q. Peak wavelength appeared at 230 nm when extracted DNA concentration was measured by absorbance.

A. Sodium Iodide remains in the extracted DNA. Measure the amount of DNA by qPCR or Threshold systems.

Q. DNA recovery rate was too low.

A. A part of DNA pellet may have been removed at steps for removing supernatant. Remove the supernatant with a pipette carefully.

#### 【Reference】

Ishizawa, M., Kobayashi, Y., Miyamura, T. and Matuura, S. (1991) Simple Procedure of DNA isolation from human serum, *Nucleic Acids Res.*, **19**, 5792.

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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