(90 × 190mm Size)

FUJIFILM

Wako

Code No. 295-55201 (for 1,000 cm²)
Code No. 291-55203 (for 5,000 cm²)

ImmunoStar[®] Reagents

(for chemiluminescence detection)

ImmunoStar reagents from Wako are designed for a simple and highly sensitive immunoblotting utilizing detection by enhanced chemiluminescence. Detection levels comparable to those reached with radioactive labels are achieved by use of a unique enhancer. Sensitivity is further enhanced by use of Wako's ABC Solution (Streptavidin and Biotin-conjugated peroxidase complex).

[Features]

1. The immunoblotting assay using these chemiluminescent reagents can detect sub-pg amounts of the target protein within 1-3 minute(s) of exposure, which is 10 to 100 times more sensitive than those utilizing chromogenic detection.
2. The chemiluminescence emission continues for several hours, enabling the detection of additional small amounts of the target protein by prolonging the exposure.
3. The membrane is reusable after washing.

[Contents]

Reagent	1,000 cm ²	5,000 cm ²
Luminescence Solution A	70 mL	330 mL
Luminescence Solution B	70 mL	330 mL
Luminescence Solution C	30 mL	120 mL

[Materials and apparatuses to be prepared for immunoblotting]

- 1) Primary antibody
The titer or dilution rate of the primary antibody should be determined before use.
- 2) Blocking solution
Usually, 5% skim milk solution is used for blocking of the blotting membrane. For detection of phosphotyrosin, BSA or gelatin is recommended as the blocker.
- 3) Secondary antibody
Approximately 0.5 μg/mL of the following listed antibodies (available from FUJIFILM Wako Pure Chemical Corporation.), can be used as the secondary antibody. The dilution rate for each of the secondary antibodies should be determined by an appropriate experiment such as dot blot immunoanalysis of IgG before use.

- 1/8 -

Code No.	Description	Package Size
010-14031	Anti Mouse IgG(H + L), Goat, Affinity Purified, Biotin-conjugated	1 mg
013-14021	Anti Rabbit IgG(H + L), Goat, Affinity Purified, Biotin-conjugated	1 mg
014-16131	Anti Rat IgG(H + L), Goat, Affinity Purified, Biotin-conjugated	1 mg

- 4) ABC solution
ABC solution, composed of streptavidin and biotin-conjugated peroxidase, is useful for intensifying the chemiluminescent signal of the target protein in immunoblotting assays. Dilute ABC solution (Wako 017-15881), 1 : 100, for use as the working solution.
- 5) Wash solution
0.05% (w/v) Tween20 and 50 mmol/L TBS (pH 7.2)
- 6) Diluent
50 mmol/L TBS (pH 7.2)
- 7) Membrane
Nitrocellulose membrane or PVDF. The signal on nitrocellulose membrane is higher than on PVDF.
- 8) Detection apparatuses
X-ray film or luminescence image analyzer such as LAS 1000, (Fujifilm Corp. Tokyo, Japan)

[Immunoblot assay]

1. Preparation of working solution

Just before use, prepare the chemiluminescent solution by mixing Chemiluminescence Reagents A, B and C at a ratio of 6 mL, 6 mL, and 2 mL, respectively, per 100 cm² of membrane.

2. Chemiluminescent staining

- 1) Prepare a protein blotted membrane, 100 cm², by either electrophoretic transfer or dot blotting.
- 2) Immerse the membrane in 50 mL of wash solution for 5 minutes.
- 3) Introduce the membrane into 12 mL of the blocking solution, and incubate for 1 hour at 37 °C.
- 4) Wash the membrane with 50 mL of wash solution for 5 minutes. Repeat.
- 5) Incubate the membrane in 10 mL of the primary antibody solution for 1 hour at 37 °C, or under pre-determined conditions for the antibody.
- 6) Wash the membrane in 50 mL of wash solution for 5 minutes. Repeat twice.
- 7) Incubate the membrane in 10 mL of the biotin-conjugated

- 2/8 -

secondary antibody in 10mL of the antibody solution for 20 minutes at 37 °C.

- 8) For signal amplification using ABC solution, repeat step 4), followed by incubation of the membrane in 10mL of ABC solution for 20 minutes at room temperature or 37 °C.
- 9) Repeat step 6).
- 10) Introduce the membrane into 14 mL of the mixed chemiluminescent solution prepared just before use, allow it to react with the chemiluminescent solution for 1 minute. Then wrap the membrane with colorless plastic wrap. Avoid forming bubbles on the membrane.
- 11) Expose the wrapped membrane to x-ray film for approximately 1 minute and develop the film, or apply the wrapped membrane to an image analyzer such as a LAS 1000, (Fujifilm Corp. Tokyo, Japan). When using x-ray film, the exposure time may vary from 0.5 to 60 minutes, depending on the intensities of the signal and backgrounds. As for using image analyser LAS 1000, the wrapped membrane is exposed for initial 5 minutes to estimate optimized exposure time for the sample.

3. Restaining of used membrane

The membrane once immunologically stained can be re-used after removal of the antibodies used for the detection. The procedure is described below.

- 1) Place the membrane in a solution containing 7 mol/L guanidine HCl, 50 mmol/L glycine, 0.05 mol/L EDTA-3Na, 0.1 mol/L KCl, and 20 mmol/L 2-mercaptoethanol. Place the solution and membrane on a rocker for 10 minutes.
- 2) Wash the membrane with 50 mL of wash solution for 10 minutes. Repeat.
- 3) The membrane is ready to be used in the immunoblotting assay protocol described above.

[Notes]

- 1) Before beginning the procedure, calculate the volumes of the solutions needed in each step according to the membrane size (cm²), based on the description in the protocol.
- 2) To conserve the primary antibody, spread several drops of the antibody solution on parafilm, cover the membrane with it and allow the membrane to react with the antibody.
- 3) Use of peroxidase inhibitors such as sodium azide in the solutions and reaction media should be avoided.

[Trouble Shooting]

Possible causes	Treatments
1. No signal detected <ul style="list-style-type: none"> • Failure blotting • Proteolytic digestion of sample 	Check the transfer process Add protease inhibitor in sample preparation

- | | |
|------------------------------------|--|
| • Property of the primary antibody | Determination the titer or optimal conditions for the immunoreaction |
|------------------------------------|--|

2. Low signal

- | | |
|------------------------------|-----------------------------------|
| • Insufficient sample | Increase the sample quantity used |
| • Insufficient exposure time | Prolong the exposure time |
| • The sample causes as in 1 | See the treatment in 1 |

3. Signal Too high

- | | |
|---|-----------------------------|
| • Excess sample amount | Decrease the sample amount |
| • Conditions for immunoreaction of the primary antibody | Optimize the immunoreaction |

4. High background

- | | |
|---|---|
| • Conditions for immunoreaction of the primary antibody | Optimize the immunoreaction |
| • Failure of blocking | Prolong the blocking time or try other blocker |
| • Contamination of apparatuses used | Use clean apparatuses |
| • Longer exposure | Diminish the exposure time or expose the membrane after 5 min on the reaction |
| • Insufficient washing in step 11) | Repeat the wash several times |

5. Non specific signal

- | | |
|--|--|
| • The property of the primary antibody | Optimize the immunoreaction |
| • Insufficient blocking | Prolong the blocking or try another blocking reagent |

Storage 2-10 °C, in the dark

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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