FUJIFILM

Wako

Code No. 295-58501 (50 assays)

for the research of forensic medicine

DNA Extractor FM Kit

[Preface]

Identification of individual by DNA analyses is powerful method in forensic medicine. However, samples for forensic DNA analyses obtained in the scene are often insoluble materials such as body hairs, nails, salivary stains and blood stains, which usually contain quite limited amounts of DNA. Therefore, it is sometimes difficult to examine DNA due to complicated hard isolation procedures or its insufficient recovery. This kit provides a rapid procedure for isolation of DNA from such the rigid materials, overcoming difficulties described above by combining rapid solubilization of DNA from the rigid materials with our original NaI DNA extraction procedure without use of harmful organic solvents such as phenol and chloroform.

[Assay Materials]

Body hair (2 pieces of 1 cm long or the less), nail (0.5 mg or the less), saliva stain, blood stain

[Features]

1. Solubilizes a variety of tough hairs as body hair completely within several ten minutes.
2. Useless of harmful organic solvents such as phenol and chloroform.
3. Avoids loss of DNA form limited amounts of sample DNA unlike solid phase DNA extraction by non-specific absorption of DNA to solid support.
4. Attains DNA purification with nearly 100% recovery by using NaI-DNA Extraction.
5. Ready to use DNA preparation obtained by the kit in DNA amplification as the polymerase chain reaction.
6. Provides a set of essential reagents for the procedure with the kit.

[Kit contents] (For 50 assays)

Lysis Solution	9.5 mL × 1
Enzyme-activated Reagent (EAR)	80 mg × 1
Reconstitution Solution for EAR	500 µL × 1
Protease	10 mg × 1
Sodium Iodide Solution	12.5 mL × 1
Washing Solution (A)	50 mL × 1
Washing Solution (B)	50 mL × 1

[Storage]

2-10 °C

[Materials and apparatuses prepared]

2-Propanol
Sterilized distilled water
Microcentrifuge
1.5mL-microfugetubes
Microdispensors

[Procedure]

1. Preparation of reagent solutions and their storage

Following working solutions should be prepared after check of the solutions provided, in which whether precipitin arose during storage. The working solutions prepared with the solutions containing precipitin may significantly reduce potency of solubilization to sample materials. Use such the solutions after precipitin is dissolved by warming at 37-60 °C.

- 1) Reconstitution of EAR (EAR Solution): Dissolve 80 mg of EAR in 0.5 mL of the reconstitution solution provided. On the storage, keep it at -20 °C. The solution is stable for 3 months even though with repeated freezing-sawing.
- 2) Preparation of Enzyme Solutions: Dissolve 10 mg of Protease provided in 0.55 mL of sterilized distilled water (SDW) (**Enzyme Solution**) except the case for nail.
For nail, use three-fold Diluted Enzyme Solution with SDW. The solution stored at -20 °C is stable for 3 months even though repeated freezing-sawing.

2. Procedures

- 1) A variety of body hairs¹
 - (1) Rinse a hair with 1 mL Ethanol by shaking centrifuge tube.
 - (2) Take out the hair from the tube, and dry it by removing ethanol with soft tissue papers.
 - (3) Cut a hair into 2 pieces of 1 cm in length. More than 2 pieces of hair may be inappropriate owing to inhibitory effect of melanin in the body hairs to DNA amplification reaction.
 - (4) Put the body hairs into a 1.5 mL microfuge tube, and add 190 µL of Lysis Solution, 10 µL of EAR Solution and 10 µL of Protease Solution^{*1} in this order, followed by vortexing.^{*2}
 - (5) Incubate the mixture at 56-60 °C for 20-40 min.^{*3}
 - (6) Centrifuge the tube at 10,000-15,000 rpm (9,600-21,600 × g) for a few minutes.^{*4}
 - (7) Carefully, transfer 200 µL of the supernate by pipette to a new tube, letting insoluble materials left in the tube.^{*5}

*1 The reagent solutions pre-mixed before addition to the tube can be used, but need vortexing to be homologous because of formation of precipitate on mixing.

*2 Make sure that sample materials sink under the solution.

*3 In the middle of and at the end of incubation, make sure resolution of sample material by vortexing.

*4 Do not use flash centrifugation.

*5 Insoluble materials contain melanin that may inhibit DNA amplification reaction.

- (8) Add 250 μ L of Sodium Iodide Solution, followed by addition of 450 μ L of 2-Propanol to the tube, and mixing well.
 - (9) Centrifuge the tube at 12,000-15,000 rpm (13,800-21,600 \times g) for 10 min at 4 °C.
 - (10) Remove the supernate, leaving the pellet^{*6} in the tube.
 - (11) Add 1 mL of Washing Solution (A) to the tube, and vortex the tube so vigorously that the pellet may detach from wall of the tube.^{*7}
 - (12) Centrifuge at 12,000-15,000 rpm (13,800-21,600 \times g) for 5-10 min at 4 °C.
 - (13) Carefully remove the supermate from the tube, leaving the pellet in the tube.
 - (14) Add 1 mL of Washing Solution (B) to the tube, and vortex well so as to detach the pellet from the wall.
 - (15) Centrifuge at 12,000-15,000 rpm (13,800-21,600 \times g) for 5-10 min at 4 °C.
 - (16) Carefully remove the supermate from the tube, leaving the pellet in the tube.^{*8}
 - (17) Dry the pellet under a reduced pressure or by warming.
 - (18) Dissolve the pellet in an appropriate volume of TE Buffer from 20 to 50 μ L.^{*9}
- 2) Blood stain and saliva stainⁱⁱ
- (1) Excise blood or saliva stain to approximately 5 \times 5 mm, and process maximally 2 sheets of approximately 5 \times 5 mm for an extraction.
 - (2) Put the excised stain in a 1.5 mL microfuge tube, and add 190 μ L of Lysis Solution, 10 μ L of EAR Solution and 10 μ L of Enzyme Solution, in this order^{*1}, followed by vortexing.^{*2}
 - (3) Incubate for 3-4 hours at 56-60 °C.^{*3}
 - (4) Remove the solution from the tube as much as possible, and transfer it to a new tube.^{*5}
 - (5) Add the 2.5 volume of NaI and 4.5 volume of 2-Propanol to the 2 volume of supernate.
 - (6) Centrifuge at 12,000-15,000 rpm (13,800-21,600 \times g) for 10 min at 4 °C.
 - (7) Carefully, remove the supernate from the tube, leaving the pellet^{*6} in the tube.
 - (8) Wash the pellet with 1 mL of Washing Solution (A) by adding it to the tube and vortexing.^{*10}

^{*6} Pellet can be seen due to glycogen as a carrier.

^{*7} Although DNA-amplification may be less reproducible, washing with the solution (A) from (11)-(13) can be omitted.

^{*8} At this stage, the pellet may be transparent, and hard to see.


^{*9} DNA can be amplified with 1/2-1/10 volume of sample, depending upon the samples. Some staining reagents for body may be inhibitory for DNA amplification.

^{*10} Although DNA-amplification may be less reproducible, washing with the solution (A) from (8)-(10) can be omitted.

^{*11} Nail of 1 mm squares weigh about 0.5 mg. Smaller pieces of nail, faster resolution by the solution. Thus, cut the nail into small pieces more than 2. Nail polish and the like should be clean up before extraction.

^{*12} Although DNA-amplification may be less reproducible, washing with the solution (A) from (11)-(13) can be omitted.

- (9) Centrifuge at 12,000-15,000 rpm (13,800-21,600 × g) for 5-10 min at 4 °C.
- (10) Carefully, remove the supernate from the tube, leaving the pellet in the tube.*⁷
- (11) Add 1mL of Washing Solution (B) and vortexing vigorously.
- (12) Centrifuge the tube at 12,000-15,000 rpm (13,800-21,600 × g) for 5-10 min at 4 °C.
- (13) Carefully, remove the supernate from the tube, leaving the pellet in the tube.*⁸
- (14) Dry the pellet under a reduced pressure or by warming.
- (15) Dissolve the pellet in an appropriate volume of TE Buffer from 20 to 50 μ L.*⁹
- 3) Nailⁱⁱⁱ
- (1) Cut about 0.5 mg of nail into small pieces of 0.5 × 1 mm.*¹¹
- (2) Put the pieces in a 1.5 mL microfuge tube, and add 190 μ L of Lysis Solution, 10 μ L of EAR Solution and 10 μ L of Diluted Enzyme Solution, in this order, followed by vortexing.*²
- (3) Incubate the tube at 56-60 °C.*³
- (4) Add 10 μ L of Diluted Enzyme Solution, to the tube 2 times every 4 hour interval, and then continue incubation over night.
- (5) Certify resolution of nail by vortexing the tube.
- (6) Centrifuge the tube at 10,000-15,000 (9,600-21,600 × g) for a few minutes.*⁴
- (7) Carefully, remove and transfer 200 μ L of the supernate to a new tube, leaving insoluble materials in the tube.
- (8) Add 250 μ L of Sodium Iodide Solution and 450 μ L of 2-Propanol to the new tube in this order, followed by vortexing.
- (9) Centrifuge the tube at 12,000-15,000 rpm (13,800-21,600 × g) for 10 min at 4 °C.
- (10) Carefully, remove the supernate from the tube, leaving the pellet*⁵ in the tube.*⁶
- (11) Wash the pellet with 1 mL of Washing Solution (A) by adding it to the tube and vortexing.*¹²
- (12) Centrifuge at 12,000-15,000 rpm (13,800-21,600 × g) for 5-10 min at 4 °C.
- (13) Carefully, remove the supernate from the tube, leaving the pellet in the tube.*⁷
- (14) Add 1 mL of Washing Solution (B) and vortexing vigorously.
- (15) Centrifuge the tube at 12,000-15,000 rpm (13,800-21,600 × g) for 5-10 min at 4 °C.*⁸
- (16) Carefully, remove the supernate from the tube, leaving the pellet in the tube.
- (17) Dry the pellet under a reduced pressure or by warming.
- (18) Dissolve the pellet in an appropriate volume of TE Buffer from 20 to 50 μ L.*⁹

ii

Blood or saliva stain (5mm×5mm)

←190 μL Lysis solution
←10 μL EAR solution
←10 μL Enzyme solution

Incubate for 3-4h at 56-60°C

Cent. 9,600-21,600 × g for a few minutes

Supernate (For example 200 μL)

←250 μL NaI solution }
←450 μL 2-propanol } * Sup 2 : NaI 25 : 2-propanol 4.5

Cent. 13,600-21,600 × g for 10min, 4°C

Pellet

←1mL Washing solution (A)

Cent. 9,600-21,600 × g for 5-10min, 4°C

Pellet

←1mL Washing solution (B)

Cent. 9,600-21,600 × g for 5-10min, 4°C

Pellet

Dry by warming or reduced pressure

←20-50 μL TE Buffer

DNA solution

iii

Nail (less than 0.5mg with 0.5×1 mm squares)

←190 μL Lysis solution

←10 μL EAR solution

←10 μL Diluted enzyme solution, twice every 4 hour interval

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