(105 × 230mm Size)

FUJIFILM

Wako

Code No. 299-58901 (For 20 sheets)

Silver Stain MS Kit

This silver staining kit is produced for mass spectrometric sequencing of protein separated by polyacrylamide gel electrophoresis. By the silver staining, protein is rarely modified chemically due to omitting treatment of glutaraldehyde and is detected at several nanograms level on the electrophoretic gel.

[Contents]

1. Enhancing Stock Solution	200 mL × 1 vial
2. Staining Stock Solution	200 mL × 1 vial
3. Developing Stock Solution	100 mL × 1 vial
4. Developing Powder	20 g × 1 vial
5. Stopper	200 mL × 1 vial
6. De-staining Solution A	50 mL × 1 vial
7. De-staining Solution B	50 mL × 1 vial

[Notes]

1. Use deionized or distilled water.
2. Before staining, ensure that all labware and equipment to be used are clean and have a smooth surfaces.
3. Be careful not to damage the gel, and use disposable gloves during the procedure to avoid direct contact with the gel.
4. Ensure that the gel is completely immersed in each solution throughout the procedure.
5. At Procedure 5, ensure that the gel is shaken gently but enough to avoid sticking to the tray during subsequent step.
6. The used staining solution may form explosive silver amide if it is left as it is. As soon as staining is completed, add sodium chloride or 1N hydrochloric acid to the used staining solution until it becomes white and discard it.

[Staining Procedure]

<Preparation>

The following solutions are prepared for staining a 100 × 100 mm polyacrylamide gel in 1 mm thickness.

1. Fixing solution-1 : Mix 50 mL methanol, 5 mL acetic acid and 45 mL deionized water.
2. Fixing solution-2 : Mix 50 mL methanol with 50 mL deionized water.
3. Enhancing solution : Mix 10 mL Enhancing Stock Solution with 90 mL deionized water.
4. Staining solution : Mix 10 mL Staining Stock Solution with 90 mL deionized water.
5. Developing solution : Mix 5 mL Developing Stock Solution and 1 g Developing Powder with 95 mL deionized water.

<Staining Procedure>

The following solutions are prepared for staining a 100 × 100 mm polyacrylamide gel in 1 mm thickness.

1. Soak the gel in 100 mL Fixing solution-1 in the tray and shake it gently for 20 minutes.
2. After discarding the solution, shake the gel gently for 10 minutes in 100 mL Fixing solution-2.
3. After discarding the solution, shake the gel gently for 10 minutes in 100 mL deionized water.
4. After discarding the solution, shake the gel gently for 1 minutes in 100 mL Enhancing solution.
5. After discarding the solution, shake the gel gently for 1 minutes in 100 mL deionized water. Repeat this step two times.
6. After discarding the solution, shake the gel gently for 20 minutes in 100 mL Staining solution.
7. After discarding the solution^{*1}, shake the gel gently for 1 minute^{*2} in 100 mL deionized water. Repeat this step two times.

*1···The used staining solution may form explosive silver amide if it is left as it is. As soon as staining is completed, add sodium chloride or 1N hydrochloric acid to the used staining solution until it becomes white and discard it.

*2···Do not shake the gel in deionized water for more than 1 minute, or detection sensitivity may decrease in a detection step.

8. After discarding the solution, shake the gel gently for 3 to 10 minutes until protein bands become moderate intensity in 100 mL Developing solution. During the development, change the Developing solution when the solution become turbid.
9. After protein bands become moderate intensity, immediately add 10 mL Stopper to the tray and shake gently for 1 minutes.
10. After discarding the Developing solution, shake the gel gently for 1 minutes in 100 mL deionized water. Repeat this step 3 times.
11. For mass spectrometric sequencing, excited band from the gel can be de-stained by treating with 0.1 mL mixture of De-staining Solution A and De-staining Solution B in equal quantity for 15 minutes at room temperature.

[Technical Notes]

1. When intensity's of the protein bands are weaker, shorten the washing time shorter than those describe in 5 and 7 of staining Procedure.
2. For double staining following CBB staining – Decolorize the background completely by soaking the CBB stained gel for at least 1 hour in Fixing solution-1. Then proceed on with staining procedure 2.
3. Background staining may be caused by one of the following reasons :
 - a. Insufficient rinsing in staining procedure 5 and 7.
 - b. Insufficient shaking in staining procedure 6 and 8.

4. In the case of overstaining – Destain using following method. Avoid destaining chromatic picture at the same time.
 - a. Preparation of destaining solution
Dissolve 3.2 g of $\text{Na}_2\text{S}_2\text{O}_3$ and 40 mg of Na_2CO_3 in 200 mL deionized water.
 - b. Soak the gel in a tray containing 200 mL destaining solution and shake for several minutes. Discard the solution and wash the gel several times within 300 to 400 mL deionized water.
5. Storage of gel – After rinsing the staining gel with deionized water, shake it for 5 minutes in 1% glycerin solution.
Dry with a gel dryer and seal it when storing. If a gel dryer is not used, the gel with the deionized water should be sealed and stored in a plastic bag without contact with air or hands directly.

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