Dietary Fiber Assay Kit

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Dietary Fiber is defined in Japan as products determined by an enzymatic-gravimetric method of Prosky *et al.* (1992) (AOAC 985.29), which is universally applied, or products undecomposed by a series of enzymatic treatment with α -amylase, protease and amyloglucosidase such as polysaccharide and lignin are defined as Dietary Fiber in "Analytical Methods for Nutrition Labeling in Japan".

This is a kit to determine Dietary Fiber with a simplified Prosky method on the Enzyme treatment.

【Kit Contents】

| | |
|---|----------------|
| ① Thermostable α -Amylase Solution | 20 mL × 1 vial |
| ② Protease Solution | 20 mL × 1 vial |
| ③ Amyloglucosidase Solution | 20 mL × 1 vial |
| ④ Diatomaceous Earth | 100 g × 1 vial |

【Principle】

(1) Determination of total Dietary Fiber Contents

Digestible components in samples, such as starch and protein, are hydrolyzed with α -amylase, protease, and amyloglucosidase, precipitate undecomposed residual macromolecule by adding quadruple volume 95 % Ethanol, and then trapped with a glass filter crucible containing diatomaceous earth.

Total Dietary Fiber Contents are determined by deducting weight of protein and ash content from the dry weight of trapped residue on a glass filter crucible.

(2) Fractional Analysis of Soluble Dietary Fiber Contents and Insoluble Dietary Fiber Contents

Digestible components in samples, such as starch and protein, are hydrolyzed with α -amylase, protease, and amyloglucosidase, then filtrated with a glass filter crucible containing diatomaceous earth.

Insoluble Dietary Fiber Contents are determined by deducting weight of protein and ash content from the dry weight of trapped insoluble residue on a glass filter crucible.

Meanwhile, add a quadruple volume 95% Ethanol into the filtrate, precipitate undecomposed residual soluble macromolecule, and trap with a glass filter crucible containing diatomaceous earth. Soluble Dietary Fiber Contents are determined by deducting weight of protein and ash content from the dry weight of the trapped precipitated residue on a glass filter crucible.

Total Dietary Fiber Contents are defined as the sum of Insoluble Dietary Fiber Contents and Soluble Dietary Fiber Contents.

(3) Determination of Low-molecular Soluble Dietary Fiber

Specific Low-molecular Soluble Dietary Fibers, which are not precipitated even by addition of a quadruple volume of 95% Ethanol, are determined by high-performance liquid chromatography.

After digestible components such as starch and protein in samples are hydrolyzed with α -amylase, protease, and amyloglucosidase, the undecomposed residual macromolecule are precipitated by adding a quadruple volume of 95% Ethanol. Then filter it with a glass filter crucible containing diatomaceous earth.

Proteins, organic acids and inorganic salts in the filtrate are removed with ion-exchange resin, and then the solution is applied to high-performance liquid chromatography. Low-molecular Soluble Dietary Fiber Content is determined by the ratio of the peak area of the sample with that of glucose (internal standard) on the chromatogram.

Meanwhile, Insoluble and Soluble Macromolecular Dietary Fiber Content is determined by deducting weight of protein and ash content from the dry weight of the trapped precipitated residue on a glass filter crucible.

Total Dietary Fiber Contents is defined as sum of Low-molecular Soluble Dietary Fiber Contents and Insoluble and Soluble Macromolecular Dietary Fiber Contents.

【Features】

- Easy Operation and good working efficiency since this is an enzyme reagent kit with a simplified Prosky method on the enzyme treatment.

Assay Method

(1) Determination of Total Dietary Fiber Contents

【Materials required but not supplied】

〈Reagents〉

1. MES : Special Grade
2. TRIS : Special Grade
3. Sodium Chloride : Special Grade
4. Calcium Chloride Dihydrate : Special Grade
5. 95% (v/v) Ethanol : Special Grade
6. 78% (v/v) Ethanol : Mix 95% (v/v) Ethanol and water [4 : 1 (volume ratio)]
7. Acetone : Special Grade
8. Analytical reagents for Kjeldahl method : Concentrated sulfuric acid, hydrogen peroxide, Proteolysis-accelerator, sodium hydroxide, boric acid, Indicator, Sulfuric Acid Standard Solution, etc.

〈Appliance〉

1. 500 mL tall beaker, graduated
2. Analytical balance : measurable to the tenth of a milligram
3. Water bath
4. Thermostatic bath : Bath in which a 500 mL tall beaker can continuously be shaken while maintaining a temperature of 60°C

5. Glass filter crucible : Gooch type 2G (pore size 40-50 mm), about 4 cm filter diameter (2G2) is preferable. Wash thoroughly after purchase and heat at $525 \pm 5^\circ\text{C}$ for 1 hour.
6. Filtration equipment : such as a suction pump with a filtering flask (or multiple aspirator)
7. Adapter : for attachment of a glass filter crucible to a filtering flask
8. Electric furnace : adjustable to $525 \pm 5^\circ\text{C}$
9. Dry oven : adjustable to $105 \pm 5^\circ\text{C}$ and $130 \pm 5^\circ\text{C}$
10. Desiccator : with silica gel as a desiccant
11. pH meter
12. Kjeldahl apparatus : flask, heating apparatus, ammonia distillation apparatus, titrator, etc.

[Assay Method]

<Measurement setup>

[Glass filter crucible]

Add 1.0 g of Diatomaceous Earth into a glass filter crucible before use, and wash the Diatomaceous Earth carefully with water, 78 % Ethanol and 95 % Ethanol, in that order until the leakage of microparticles of diatomaceous earth is no longer observed, and obtain a homogeneous Diatomaceous Earth layer. Heat at 130°C for 1 hour and stand to cool in a desiccator. Determine the constant weight to the tenth of a milligram and store in a desiccator until use.

[Sample Preparation]

Low moisture food items such as grains, beans, nuts and seeds are ground evenly. Ground foods are then sifted through a 30-mesh (0.50 mm) sieve. The sifted part is used as the sample. However, if the food contains lipids at a concentration of 5 % or more in the foods such as soybean, nuts and seeds, defatting is conducted in advance.

High moisture food stuffs such as potatoes, vegetables, fungi are homogenized and freeze-dried (or dried overnight at 70°C). Then they are ground evenly and sifted through a 30-mesh (0.50 mm) sieve.

Food items with high sugar content such as fruits are homogenized and sifted through a 4-mesh (4.75 mm) sieve. In case the food has high acid, it affects the pH. You should do pH adjustment before the enzymatic treatment. In addition, loss in weight due to drying and defatting is determined.

[Buffer]

Preparation of 50 mmol/L MES-Tris Buffer, containing 10 mmol/L Na and 3 mmol/L Ca (pH 6.3, 24°C) :

- Preparation of 50 mmol/L MES solution : Weigh 9.76 g of MES (the final concentration : 50 mmol/L), 0.584 g of NaCl (the final concentration : 10 mmol/L) and 0.441 g Calcium Chloride Dihydrate (the final concentration : 3 mmol/L) and add them into a 1 L beaker, graduated. Then add deionized water to make the 1 L solution.
- Preparation of 50 mmol/L TRIS solution : Weigh 6.055 g of TRIS (the final concentration: 50 mmol/L), 0.584 g of

Sodium Chloride (the final concentration: 10 mmol/L) and 0.441 g of Calcium Chloride Dihydrate (the final concentration : 3 mmol/L) and add them into a 1 L beaker, graduated. Then add deionized water to make the 1 L solution.

- Preparation of 50 mmol/L MEW-TRIS Buffer : Mix the 50 mmol/L MES solution and 50 mmol/L TRIS solution and adjust the pH to 6.3 at 24°C . [MES solution : TRIS solution = 64 : 36 (v/v)]

<Assay>

1) Digestion with thermostable α -amylase

Accurately weigh 1 g of a sample to the tenth of a milligram. Take two samples (W1, W2) for one analysis. The weight difference between the two samples shall be within 20 mg. Use one sample to determine undigested protein and the other to determine ash content. Add a sample to 40 mL of 50 mmol/L MES-TRIS buffer (pH 6.3, 24°C) in a 500 mL tall beaker and disperse it. Check the pH to make sure it is 6.3 ± 0.1 (24°C). Add 0.2 mL of thermostable α -amylase and cover the tall beaker with aluminum foil. Let the mixture react in a boiling water bath for 30 min. while shaking it every 5 min.

2) Digestion with protease and amyloglucosidase

Add 10 mL of water while washing the wall of the tall beaker, cool to 60°C , and add simultaneously 0.2 mL of protease and 0.2 mL of amyloglucosidase. Let the mixture react in a thermostatic bath for 30 min. while continuous shaking.

3) Ethanol Precipitation

Add a quadruple volume of 95% Ethanol (preheated at 60°C) to the enzyme treated solution obtained, and allow to stand for 1 hour at room temperature to precipitate enzyme undigested macromolecule.

4) Filtration and Wash

Conduct a suction filtration of Ethanol Precipitation solution using a glass filter crucible, which has a homogeneous Diatomaceous Earth layer by washing with 78% Ethanol just before use, to obtain a residue. Under suction, rinse the residue in the tall beaker with 78% Ethanol (20 mL \times 3 times), 95% Ethanol (10 mL \times 2 times), and acetone (10 mL \times 2 times), and add the washings to the residue in the glass filter crucible.

5) Drying and Weighing

Dry the residue with the glass filter crucible at $105 \pm 3^\circ\text{C}$ overnight. Weigh the residue to a tenth of a milligram after allowing to cool in a desiccator. Residues for protein determination and that for ash determination are set at R1 and R2, respectively.

6) Protein Determination

Scrape the total amount of residues for protein determination together with Diatomaceous Earth and determine the nitrogen content in the residue by Kjeldahl protein determination. Multiply the obtained nitrogen content by 6.25 to obtain the protein content (P1).

7) Ash determination

Conduct incineration treatment of the residues for ash determination with the glass filter crucible at $525 \pm 5^\circ\text{C}$ for 5 hours, allow to cool in a desiccator, weigh to the tenth of a milligram and obtain the ash content (A1).

8) Blank Test

Perform the above-mentioned procedure without samples and obtain the blank values r1, r2, p1 and a1.

〈Calculation〉

Calculate Total Dietary Fiber Content using the following formula :

$$\text{Total Dietary Fiber Content (g/100 g)} \\ = \frac{(R1 + R2)/2 \times \{1 - (P1/R1 + A1/R2)\} - B}{(W1 + W2)/2} \times 100 \\ B(g) = (r1 + r2)/2 \times \{1 - (p1/r1 + a1/r2)\}$$

In the case of dried samples, it is converted to Dietary Fiber Content before drying using the following formula :

$$\text{Dietary Fiber Content before drying (g/100 g)} \\ = \text{Total Dietary Fiber Content} \times (1 - D/100) \\ D : \text{Loss on drying (\%)}$$

In the case of defatted sample, it is converted to Dietary Fiber Content before defatting using the following formula :

$$\text{Dietary Fiber Content before defatting (g/100 g)} \\ = \text{Total Dietary Fiber Content} \times (1 - F/100) \\ F : \text{Loss on defatting (\%)}$$

(2) Fractional Analysis of Soluble Dietary Fiber Contents and Insoluble Dietary Fiber Contents

[Materials required but not supplied]

〈Reagents〉 Same as the above-mentioned (1)
〈Appliance〉 Same as the above-mentioned (1)

[Assay Method]

〈Measurement setup〉

[Glass filter crucible] : Same as the above-mentioned (1)
[Sample Preparation] : Same as the above-mentioned (1)
[Buffer] : Same as the above-mentioned (1)

〈Assay〉

- 1) Digestion by thermostable α -amylase :
Same as the above-mentioned (1)
- 2) Digestion by protease and amyloglucosidase :
Same as the above-mentioned (1)
- 3) Filtration of soluble fraction and insoluble fraction
Conduct a suction filtration of the enzyme treated solution obtained using a glass filter crucible, which has a homogeneous Diatomaceous Earth layer by washing with 78 % Ethanol just before use, and divide into residue and filtrate. Rinse the residue in the tall beaker with about 10 mL of water, and add the washings to the residue in the glass filter crucible. Combine the washings and

the filtrate. The filtrate is the fraction containing Soluble Dietary Fiber, and the residue is the fraction containing Insoluble Dietary Fiber.

4) Ethanol Precipitation, Filtration and Wash of fraction contains soluble dietary fiber

Add the quadruple volume of 95% Ethanol (preheated at 60°C) in the filtrate obtained at the procedure of 3) and leave it for 1 hour at room temperature to precipitate enzyme undigested macromolecule. Conduct a suction filtration of Ethanol Precipitation solution using a glass filter crucible, which has a homogeneous Diatomaceous Earth layer by washing with 78% Ethanol just before use, to obtain a residue. Under suction, add 78% Ethanol (20 mL \times 3 times), 95% Ethanol (10 mL \times 2 times), and acetone (10 mL \times 2 times) successively to the residue in the glass filter crucible.

5) Filtration and Wash of fraction contains insoluble dietary fiber

Under suction, add 95% Ethanol (10 mL \times 2 times), and acetone (10 mL \times 2 times) successively to the residue in the glass filter crucible obtained in procedure 3).

6) Drying and Weighing

Dry the residue with the glass filter crucible at $105 \pm 3^\circ\text{C}$ overnight. Weigh the residue to a tenth of a milligrams with an analytical balance after allowing to cool in a desiccator.

R1 : Residue for protein determination containing Soluble Dietary Fiber fraction
R2 : Residue for ash determination
R3 : Residue for protein determination containing Insoluble Dietary Fiber fraction
R4 : Residue for ash determination

7) Protein Determination

R1 and R3, which are residues for protein determination, are scraped with Diatomaceous Earth, and then nitrogen content in the residue is determined with Kjeldahl method. Multiplying the obtained nitrogen content by 6.25 gives protein contents P1 and P2, respectively.

8) Ash Determination

Conduct incineration treatment of the residues for ash determination (R2 and R4) with the glass filter crucible at $525 \pm 5^\circ\text{C}$ for 5 hours, allow to cool in a desiccator, weigh to the tenth of a milligram and obtain ash contents A1 and A2, respectively.

9) Blank Test

Perform the above-mentioned procedure without samples and obtain the blank values r1, r2, r3, r4, p1, p2, a1 and a2.

〈Calculation〉

Calculate Soluble Dietary Fiber, Insoluble Dietary Fiber, and Total Dietary Fiber Content are using the following formulas :

Soluble Dietary Fiber Content (g/100 g)

$$= \frac{(R1 + R2)/2 \times \{1 - (P1/R1 + A1/R2)\} - Bs}{(W1 + W2)/2} \times 100$$

$$Bs(g) = (r1 + r2)/2 \times \{1 - (p1/r1 + a1/r2)\}$$

Insoluble Dietary Fiber Content (g/100 g)

$$= \frac{(R3 + R4)/2 \times \{1 - (P2/R3 + A2/R4)\} - Bi}{(W1 + W2)/2} \times 100$$

$$Bi(g) = (r3 + r4)/2 \times \{1 - (p2/r3 + a2/r4)\}$$

In the case of dried samples, it is converted to Dietary Fiber Content before drying using the following formula :

$$\begin{aligned} &\text{Dietary Fiber Content before dried (g/100 g)} \\ &= \text{Soluble or Insoluble Dietary Fiber Content} \times \\ &\quad (1 - D/100) \\ &\quad D : \text{Loss on drying (\%)} \end{aligned}$$

In the case of defatted sample, it is converted to Dietary Fiber Content before defatting using the following formula :

$$\begin{aligned} &\text{Dietary Fiber Content before defatted (g/100 g)} \\ &= \text{Soluble or Insoluble Dietary Fiber Content} \times \\ &\quad (1 - F/100) \\ &\quad F : \text{Loss on defatting (\%)} \end{aligned}$$

Total Dietary Fiber Content (g/100 g)

$$= \text{Soluble Dietary Fiber Content} + \text{Insoluble Dietary Fiber Content}$$

(3) Determination of Low-molecular Soluble Dietary Fiber such as indigestible dextrin, polydextrose**[Materials required but not supplied]**

〈Reagents〉 Required reagents besides (1)

9. LabAssay™ Glucose (Wako Catalog No. 638-50971 (500 tests))

〈Appliance〉 Required appliances besides (1)

13. rotary evaporator

14. membrane-filter : 0.45 μm

15. high-performance liquid chromatography (HPLC) : degasifier, pump, differential refractive index detector, etc.

16. analytical column: gel filtration or ligand exchange column

17. ion-exchange column: one-to-one mixture of OH and H form resin, or equivalent product

[Assay Method]

〈Measurement setup〉

[Glass filter crucible] : Same as the above-mentioned (1)

[Sample preparation] : Same as the above-mentioned (1)

[Buffer] : Same as the above-mentioned (1)

〈Assay〉

1) Digestion by thermostable α-amylase :

Same as the above-mentioned (1)

2) Digestion by protease and amyloglucosidase :

Same as the above-mentioned (1)

3) Ethanol Precipitation :

Same as the above-mentioned (1)

4) Filtration and Wash

Conduct a suction filtration of Ethanol Precipitation solution using a glass filter crucible, which has a homogeneous Diatomaceous Earth layer by washing with 78% Ethanol just before use, to obtain residue and filtrate. Under suction, rinse the residue in the tall beaker with 78% Ethanol (20 mL×3 times), 95% Ethanol (10 mL×2 times), and acetone (10 mL×2 times), and add the washings successively to the residue in the glass filter crucible. Combine the washings until 95% Ethanol and the filtrate. The filtrate is concentrated and ethanol in the filtrate is removed with a rotary evaporator to adjust the volume of the solution to 100 mL. The solution is the fraction of Low-molecular Soluble Dietary Fiber.

5) Drying and Weighing : Same as the above-mentioned (1)

6) Protein Determination :

Same as the above-mentioned (1)

7) Ash determination : Same as the above-mentioned (1)

8) Blank Test : Same as the above-mentioned (1)

9) Removal of protein, organic acids and inorganic salts from the fraction of Low-molecular Soluble Dietary Fiber Elute 50 mL of the fraction of Low-molecular Soluble Dietary Fiber obtained in procedure 4) with 50 mL of ion exchange resin packed column (glass tube, I.D. 20 mm×300 mm) at SV1.0 (flow rate : 50 mL/hr). Elute with distilled water and obtain 200 mL of eluate.

10) High-Performance Liquid Chromatography

The eluate obtained in procedure 9) is concentrated with a rotary evaporator. The obtained solution is adjusted to an appropriate concentration with distilled water (ex. to Brix 5), filtered with a membrane-filter, which pore diameter is 0.45 μm and applied to HPLC.

11) Glucose Determination^{Note)}

Determine the glucose content in the eluate obtained in procedure 9) using LabAssay™ Glucose (Wako Catalog No. 638-50971 (500 tests)).

〈Calculation〉

Calculate the Low-molecular Soluble Dietary Fiber, Insoluble and soluble Macromolecular Dietary Fiber Content, and Total Dietary Fiber Content using the following formulas :

Low-molecular Soluble Dietary Fiber Content (g/100 g)

$$= \frac{\text{Peak area of Dietary Fiber}}{\text{Peak area of Glucose}} \times \text{Glucose content in eluent (g)} \times \frac{100}{(W1 + W2)/2} \dots (I)$$

Insoluble and Soluble Macromolecular Dietary Fiber Content (g/100 g)

$$= \frac{(R1 + R2)/2 \times \{1 - (P1/R1 + A1/R2)\} - B}{(W1 + W2)/2} \times 100 \dots (II)$$

$$B(g) = (r1 + r2)/2 \times \{1 - (p1/r1 + a1/r2)\}$$

Total Dietary Fiber Content (g/100 g)

= Low-molecular Soluble Dietary Fiber Content
+ Insoluble and Soluble Macromolecular Dietary Fiber Content (III)

In case of dried samples, in formulas (I)-(III), they are converted to Dietary Fiber Content before drying using the following formula :

$$\text{Dietary Fiber Content in a sample before dried (g/100 g)} \\ = \text{Dietary Fiber Content} \times (1 - D/100)$$

D : Loss on drying (%)

In case of defatted samples, in formulas (I)-(III), they are converted to Dietary Fiber Content before defatting using the following formula :

$$\text{Dietary Fiber Content in a sample before defatted} \\ (\text{g/100 g}) = \text{Dietary Fiber Content} \times (1 - F/100)$$

F : Loss on Defatting (%)

Note : In another approach, a known weight of internal standard substance can be added to the eluate obtained in procedure 9), and the Dietary Fiber Content can be calculated from the correlation between the peak area of the eluate and that of the internal standard substance on the chromatogram. Use a substance other than glycerin as an internal standard substance, because glycerin is contained in the enzyme reagent.

[Acknowledge]

This kit shows measured values of the same degree as the conventional Prosky's method in various foods. Please note the measured value may be low in foods which resistant starch is added.

[Storage Condition]

Keep at 2 ~ 10°C in a dark place.

[Package Size]

100 tests

[Measurement Outline]

1. Sample dissolved in MES/Tris buffer
↓
2. Digestion with thermostable α -amylase
↓
3. Digestion protease and amyloglucosidase
↓ for 30 min. at pH 6.3 and 60°C
4. Ethanol Precipitation
↓ Leave 1 hour.
5. Filtration
↓

6. Wash with ethanol and acetone
↓

7. Dry overnight
↓

8. Kjeldahl Protein Determination and Ash Determination
↓

9. Total Dietary Fiber Determination

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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