

90 × 210mm サイズ

**FUJIFILM**

**Wako**

Code No. 018-22381 (200  $\mu$ g)  
014-22383 ( 1 mg)  
012-22384 ( 5 mg)

**Anti DYKDDDDK tag, Monoclonal Antibody**  
**抗DYKDDDDKタグ, モノクローナル抗体**

The DYKDDDDK was consisted of 8 amino acids and has been used as the tag for detection and/or purification of recombinant protein.

Monoclonal anti DYKDDDDK tag antibody was established by fusing mouse myeloma cells and spleen cells of mice immunized with synthetic DYKDDDDK peptide which was coupled with keyhole limpet hemocyanin (KLH). This antibody was purified from mouse ascites.

This product is for laboratory use only ; use in any such application is the responsibility of the user.

**[Concentration]**

Describing at label.

**[Formulation]**

1 × PBS aqueous solution with 50% glycerol.

**[Clone No.]**

1E6

**[Subclass]**

IgG<sub>2b</sub>

**[Specificity]**

Specific for the DYKDDDDK peptide.

**[Working Dilution]**

Western blot 1 : 1,000 - 1 : 10,000

Immunoprecipitation 5  $\mu$ g/20  $\mu$ L 10% Protein A or Protein G beads slurry

**[Storage]**

Keep at -20 °C in the dark.

Do not freeze and thaw cycle.

**[Package]**

Code No.	Packaging
018-22381	200 $\mu$ g
014-22383	1 mg
012-22384	5 mg

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : + 49-2131-3111-0  
Facsimile : + 49-2131-311100  
<http://www.wako-chemicals.de>

1801KA1