(90 × 210mm Size)

FUJIFILM

Wako

Code No. 291-73501

For research use only

Mouse/Rat PYY ELISA Kit *Wako*

< Competition assay >

Package Size : 96 tests

[Introduction]

PYY is a gastrointestinal peptide hormone consisting of 36 amino acid residues. It is mainly released from endocrine cells localizing in the intestinal tract after food ingestion, and inhibits appetite through the vagus nerve and blood. The ELISA kit is used for the quantitative determination of mouse/rat PYY in plasma and serum samples.

[Kit Contents]

1. Antibody Coated Plate	1 plate
2. Mouse/Rat PYY Standard	12.5 ng
3. Biotinylated Mouse/Rat PYY	1 vial
4. Anti Mouse/Rat PYY, Rabbit	8.5 mL
5. HRP Labeled Streptavidin Solution	12 mL
6. Enzyme Substrate Solution	12 mL
7. Stopping Solution	12 mL
8. Buffer Solution	25 mL
9. Washing Solution (concentrated)	50 mL
10. Adhesive Foil	3 sheets

[Principle of the assay]

This ELISA kit for determination of PYY in mouse/rat plasma or serum samples is based on a competitive enzyme immunoassay using the combination of an antibody to mouse/rat PYY and a biotin-streptavidin affinity system. The 96 wells plate is coated with goat anti rabbit IgG, to which biotin labeled antigen, standard antigen or samples and rabbit anti mouse/rat PYY antibody are added for the competitive immunoreaction. After incubation and plate washing, horse radish peroxidase (HRP) labeled streptavidin (SA) is added to form a HRP labeled SA-biotinylated antigen-antibody complex on the surface of the wells. Finally, HRP enzyme activity is determined by 3,3',5,5'-Tetramethylbenzidine (TMB) and the concentration of mouse/rat PYY is calculated.

[Procedures]

(1) Preparations

1. Preparation of the standards :
Reconstitute the Mouse/Rat PYY Standard (lyophilized 12.5 ng/vial) with 1 mL of Buffer Solution, which affords a 12.5 ng/mL standard solution. The reconstituted standard solution (0.1 mL) is diluted with 0.2 mL of Buffer Solution which yields a 4.17 ng/mL standard solution. Continue the dilution procedure to prepare a series of one each of 1.39, 0.46, and 0.15 ng/mL standard solutions. The Buffer Solution itself is used as 0 ng/mL.
2. Preparation of the labeled antigen solution :
Reconstitute the Biotinylated Mouse/Rat PYY with 7 mL of Buffer solution.

- 1/8 -


3. Preparation of the washing solution :
Dilute 50 mL of the Washing Solution (concentrated) to 950 mL with distilled or deionized water.
4. The other reagents are ready-to-use.

(2) Assay Procedure

1. Bring all the reagents to room temperature (20-30 °C) before starting assay.
2. Add 350 μ L of washing solution to each well and aspirate the washing solution in the wells. Repeat this washing procedure two more times (total 3 times). Finally, invert the plate and firmly tap it onto an absorbent surface, such as paper toweling, to ensure blotting free most of the residual washing solution.
3. Pipette 50 μ L of the labeled antigen solution into each well first, then introduce 25 μ L of each of standard solution (0, 0.15, 0.46, 1.39, 4.17, 12.5 ng/mL) or samples and finally add 75 μ L of Anti Mouse/Rat PYY, Rabbit into the wells.
4. Cover the plate with adhesive foil and incubate it at 4 °C for 18 hours (\pm 1 hour) without shaking and then 30 minutes more at room temperature with shaking (100-150 rpm).
5. After incubation, take off the adhesive foil, aspirate and wash the wells 5 times with 350 μ L of washing solution. Finally, invert the plate and firmly tap it onto an absorbent surface, such as paper toweling, to ensure blotting free most of the residual washing solution.
6. Add 100 μ L of HRP Labeled Streptavidin Solution into each well.
7. Cover the plate with adhesive foil and incubate it at room temperature for 1 hour with shaking (100-150 rpm).
8. Take off the adhesive foil, aspirate and wash the wells 5 times with 350 μ L of washing solution. Finally, invert the plate and firmly tap it onto an absorbent surface, such as paper toweling, to ensure blotting free most of the residual washing solution.
9. Add 100 μ L of Enzyme Substrate Solution into each well ; cover the plate with adhesive foil and keep it for 30 minutes at room temperature in a dark place without shaking for the color reaction.
10. Add 100 μ L of Stopping Solution into each well to stop the color reaction.
11. Shake the plate 5-10 seconds to mix the contents and immediately read the optical absorbance of the wells at 450 nm.

(3) Standard Curve

The standard curve of PYY is obtained by measurement of the series dilution of the standard PYY solution with the kit, as is shown in the figure.

- 2/8 -

[Capability]

- (1) Dynamic range
0.15 ~ 12.5 ng/mL
- (2) Reproducibility
Intra-assay CV (%) 3.1 ~ 9.8
Inter-assay CV (%) 4.2 ~ 14.2
- (3) Cross-reactivity
Specific to Mouse, Rat
No cross reactivity to Human PYY

[Notes on use]

- 1) EDTA-2Na additive blood collection tube is recommended for the plasma collection. It is recommended that serum or plasma samples should be tested as soon as possible after separation. If the sample is tested later, they should be aliquoted and frozen below -30 °C (for long term storage, it is recommended that the sample should be stored in a -70 °C deep freezer). Avoid repeated freezing and thawing of samples. Samples should be kept in an ice bath after thawing before the assay and used within 60 minutes.
- 2) Standard and labeled antigen solutions should be prepared immediately before use. The plate can be used for twice if separated. In that case, the reconstituted reagents (standard and labeled antigen solution) should be stored at or below -30 °C.
- 3) During storage of washing solution (concentrated) at 2-8 °C, precipitates may be observed occasionally, however they will be dissolved when diluted.
- 4) As pipetting operations may affect accuracy and precision of the assay, pipette solutions especially for standard solutions or samples precisely into each well of the plate. In addition, use clean test tubes or vessels in assay and use a new tip for each sample or standard solution and for each step of preparation of the standard diluting solution to avoid cross contamination.
- 5) Perform all the determinations in duplicate.
- 6) If a high range sample (over than 12.5 ng/mL) will be tested or predicted, dilute this sample with assay buffer, before starting the assay procedure.
- 7) To quantitate accurately, always run a standard curve when measuring samples.
- 8) The color reaction should be carried out under light proof conditions.
- 9) Read optical absorbance of the reaction solution in wells as soon as possible after stopping the color reaction.
- 10) Protect reagents from strong light (e.g. direct sunlight) during storage and assay.

[Storage condition]

Keep at 2~10 °C

[Shelf life]

18 months after production

[Package size]

96 tests

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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