

90 × 190mm サイズ

**FUJIFILM**

**Wako**

<for Research Use Only>

Code No. 162-26031 (1 mL用)  
168-26033 (1 mL用 × 5)

**Protease Inhibitor Cocktail Set V  
(EDTA free) (× 100)**

**プロテアーゼ阻害剤カクテル セット V,  
(EDTA 不含) (× 100)**

This protease inhibitor cocktail is the mixture of four protease inhibitors for the inhibition of serine and cysteine proteases, but not metalloproteases.

**Appearance :** Lyophilized powder

**Reconstitution :** Reconstitute each vial with 1 mL of distilled water to obtain a 100 × stock solution.

※Precaution  
Centrifuge the vial prior to opening.

**Usage :** 100 × stock solution dilute to 1 × solution with appropriate solution or sample solution before use.

※Precaution  
Since not all organisms contain the same level of endogenous proteases, it may sometimes be necessary to increase the concentration of inhibitors.

**Contents (as 100 × stock solution) :**

Inhibitor	M. W.	Concentration	Inhibition Target
AEBSF Hydrochloride	239.69	50 mmol/L	Serine Proteases
Aprotinin, from Bovine Lung	Approx. 6.5kDa	15 μmol/L	Serine Proteases and Esterases
E-64	357.41	0.1 mmol/L	Cysteine Proteases
Leupeptin Hemisulfate	475.59	0.1 mmol/L	Cysteine Proteases and Trypsin-like Proteases

**[Expiration date]** Indicated on the label.

**[Storage]** Store at -20 °C.  
After reconstitution, please use the reconstituted solution as soon as possible.

For extended storage, store at -20 °C and avoid freeze-thaw cycles. For long time storage is not recommended.

- 1/2 -

**[Package]**

Code No.	Packaging
162-26031	1 mL用
168-26033	1 mL用 × 5

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : +81-6-6203-3741  
Facsimile : +81-6-6201-5964  
<http://www.wako-chem.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : +1-804-271-7677  
Facsimile : +1-804-271-7791  
<http://www.wakousa.com>

**FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
D-41468 Neuss  
Germany  
Telephone : +49-2131-311-0  
Facsimile : +49-2131-311100  
<http://www.wako-chemicals.de>

1802KA1

- 2/2 -