FUJIFILM

Wako

Code No. 197-17571 (100 mL)
193-17573 (100 mL × 4)

StemSure[®] hPSC Medium Δ

StemSure[®] hPSC Medium Δ is appropriate for maintaining human pluripotent stem cells (hPSCs), human ES cells and human iPSC cells, in a feeder-free, serum-free and animal-free environment. This medium is low in protein, and does not contain albumin.

This product is for laboratory research use only.

Procedures :

<Preparation of complete StemSure[®] hPSC Medium >

1. Thaw StemSure[®] hPSC Medium Δ at 2-10 °C for several hours or overnight.

NOTE : Do not thaw the frozen medium at 37 °C. After thawing, store this medium at 2-10 °C in the dark and use within a week. Avoid repeating freeze-thaw.

2. Add 35-100 ng/mL basic fibroblast growth factor (bFGF) (Code No. 064-05381, 068-05384) (final concentration) to the thawed StemSure[®] hPSC Medium Δ, and gently mix well to prepare the complete StemSure[®] hPSC Medium (sshPSC medium).

NOTE : The concentration of 35 ng/mL is recommended as a final concentration of bFGF. When you transfer hPSCs cultured in other feeder-free medium to the sshPSC medium, the same concentration of bFGF used might be suitable for the sshPSC medium also.

3. Before use, warm the sshPSC medium to room temperature.
NOTE : Do not warm the sshPSC medium in a water bath.

<Transferring of hPSCs from the culture on feeder layers to the feeder-free culture >

NOTE : When hPSCs cannot be directly transferred from feeder-dependent culture with mouse embryonic fibroblasts (MEFs) to feeder-free culture which uses the sshPSC medium, a preculture in MEF-conditioned medium [Xu, C., et al.: Nat. Biotechnol., 19, 971, (2001)], prior to using the sshPSC medium might be effective to transfer smoothly.

Preparation of Matrigel[®]-coated plates

1. Thaw Matrigel[®] at 4 °C. To avoid gelation, do not store at room temperature.
2. Dilute Matrigel[®] in cold DMEM/Ham's F-12 and gently mix well by pipetting (300 μL of Matrigel[®] in 25 mL of D-MEM/Ham's F-12).
3. Coat 6 well-plate with 1.4 mL of the diluted Matrigel[®] per well (9.5 cm²).
4. Incubate at room temperature for at least 1 hour before use.

Passaging

1. Aspirate the culture medium of a feeder-dependent cultured well.
2. Add 0.4 mL of CTK solution [Suemori, H., et al.: Biochem. Biophys. Res. Commun., 345, 926, (2006).] per well of 6-well plate.
3. Incubate at 37 °C with 5% CO₂ for 5 minutes.
4. Aspirate the CTK solution, and add 2 mL per well of the sshPSC medium and suspend the colonies in small clumps of about 100 cells by pipetting with a P1000 micropipette.
5. Transfer the cell suspension to a 15mL tube, and add 6 mL of the sshPSC medium and gently mix by inverting the tube several times.
6. Spin at 1,000 rpm (about 170 × g) at room temperature for 3 minutes.
7. Remove the supernatant, and add 2 mL per well of ROCKi + medium[※] and resuspend the pellet.
NOTE : Do not break up the clumps to a smaller number of cells, especially single cells.
※ROCKi + medium : sshPSC medium with 10 μmol/L Y-27632, which is an inhibitor of p160-Rho-associated coiled kinase (ROCK).
8. Aspirate the coating solution from the Matrigel[®]-coated wells.
9. Transfer 2 mL of the cell suspension to the well.
NOTE : Do not plate the cells at low-density. This will lead to the decrease in the viability of the cells and the increase in the differentiation of the plated cells.
NOTE : After the transfer to the feeder-free condition, MEFs will remain present for the first 1-3 passages. To eliminate MEFs from the culture, the cell suspension is plated to a non-coated or a gelatin-coated well, and incubate at 37 °C with 5% CO₂ for about 1 hour. Residual MEFs, not hPSC clumps, adhere to the well. The suspension containing hPSC clumps is harvested from the well, and then plated to a Matrigel[®] coated-well.
10. Incubate at 37 °C with 5% CO₂.
11. Replace the medium daily with 2 mL of sshPSC medium per well.
12. Passage the cells at 5 to 7 days.

<Passaging in sshPSC medium >

1. Aspirate the culture medium, and rinse the cells once with D-PBS(-).
2. Remove D-PBS(-), and add 1 mL of StemPro[®] Accutase, TrypLE[™] Select or TrypLE[™] Express per well.
NOTE : Do not use trypsin, Dispase[®] and CTK solution.
3. Incubate at 37 °C with 5% CO₂ for 5 minutes.
4. Add 2 mL of the ROCKi + medium, and suspend the colonies to single cells by using a P1000 micropipette.
5. Transfer the cell suspension to a 15 mL tube.
6. Spin at 1,000 rpm (about 170 × g) at room temperature for 3 minutes.
7. Aspirate the supernatant and suspend the pellet with 2 mL of the ROCKi + medium.
8. Count the viable cell numbers.
9. Aspirate the coating solution from the Matrigel[®]-coated wells (see Preparation of Matrigel[®]-coated plates).

10. Add 2 mL of the ROCKi + medium per well, and plate 1×10^5 cells per well.

NOTE : Do not plate the cells at low-density. This will lead to the decrease in the viability of the cells and the increase in the differentiation of the plated cells.

NOTE : In early passages, when the cells do not adapt to the feeder-free culture, $2-3 \times 10^5$ cells per well should be plated. The optimal split ratio is different in hPSC clones. The optimal plating cell number should be adjusted carefully for reproducible culture. In general, the cells should reach 80-90% confluent of the well at 4 or 5 days after plating.

11. Incubate at 37 °C with 5% CO₂.
12. The next day (at day 1), replace the medium with 2 mL of the sshPSC medium per well, not the ROCKi + Medium.
13. Replace the medium daily with 2 mL of the sshPSC medium per well. The replacement of medium is unnecessary on day 2.
14. Passage the cells when cells grow to cover approximately 80-90% of the well, usually every 4-6 days.

Storage : Store at -20 °C in the dark.
Avoid repeating freeze-thaw. After thawing, store at 2-10 °C in the dark and use within a week.

Package :

Code No.	Packaging
197-17571	100 mL
193-17573	100 mL × 4

Related Reagents

- D-MEM/Ham's F-12 with L-Glutamine, Phenol Red, HEPES and Sodium Pyruvate (Code No.042-30555)
- Fibroblast Growth Factor (basic), Human, recombinant, Animal-derived-free (Code No. 064-05381, 068-05384)
- Y-27632 (Code No. 257-00511, 253-00513, 251-00514)
- D-PBS(-) (Code No. 045-29795)
- Matrigel® hESC-Qualified Matrix (Code No. 643-55461, Corning No. 354277)
- 6 Well Clear TC-Treated Multiple Well Plates (Code No. 649-01111, Corning No. 3516)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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