(90 × 210mm Size)

FUJIFILM

Wako

Code No. 017-25151 (2.5 mL)
013-25153 (2.5 mL × 5 vials)
011-25154 (200 μL)

For immunochemistry

Anti DYKDDDDK tag Antibody magnetic beads

Anti DYKDDDDK (Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys) tag Antibody magnetic beads is the slurry of anti DYKDDDDK monoclonal antibody-immobilized magnetic beads. This product can be used for the purification and immunoprecipitation of DYKDDDDK-tagged recombinant proteins. DYKDDDDK-tagged recombinant proteins can be eluted by competitive elution with DYKDDDDK tag Peptide under neutral pH condition after immunoprecipitation.

This product is intended for laboratory research use only.

[Formulation]

10 mg/mL Anti DYKDDDDK tag Antibody magnetic beads, 1 × TBS (pH 7.4), 50w/v% glycerol, 0.05w/v% sodium azide.

[Immobilized antibody]

Anti DYKDDDDK mouse monoclonal antibody (1E6)

[Subclass of Immobilized antibody]

IgG2b

[Assay]

50 μL/immunoprecipitation

[Storage]

- 20 °C

[Additional required materials]

Magnetic stand, Centrifuge, Rotator, Vortex mixer, Tube mixer, Centrifuge tube, 1.5 mL microcentrifuge tube, Cell lysis buffer (see below), PBS (-), TBS, Protease inhibitors, Phosphatase inhibitors.

Recommended cell lysis buffer

- RIPA Buffer (Code No. 182-02451)
50 mmol/L Tris-HCl, pH 8.0, 150 mmol/L Sodium chloride, 0.5w/v% Sodium Deoxycholate, 0.1w/v% Sodium Dodecyl Sulfate, 1.0w/v% NP-40 substitute
- Cell Lysis Buffer M (Code No. 038-21141)
20 mmol/L Tris-HCl, pH 7.4, 200 mmol/L Sodium Chloride, 2.5 mmol/L Magnesium Chloride, 0.05w/v% NP-40 substitute

When inhibitors are added to a cell lysis buffer, please use it in accordance with manufacturer's instruction.

[Protocol]

The immunoprecipitation method for DYKDDDDK-tagged recombinant proteins (from mammalian cells)

1. The preparation of cell lysates

A. Adherent cells

- 1) Culture the cell lines expressing DYKDDDDK-tagged recombinant proteins ($1 \times 10^6 \sim 10^7$).

- 1/8 -

- 2) Remove the culture medium and wash the cells twice with PBS (-).
- 3) Dissociate the cells using a cell scraper or trypsin-EDTA solution from dishes, and transfer the cell suspension into a new centrifuge tube. In the case of using trypsin-EDTA solution, add the medium including serum to the cell suspension immediately for suppressing damage by trypsin.
- 4) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 5) Remove the supernatant.
- 6) After adding PBS (-) and resuspending the cell pellet with PBS (-), centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C. Repeat step 5) to 6).
- 7) Remove the supernatant.
- 8) Add 1 mL of cell lysis buffer to the cell pellet and suspend cells by pipetting or vortex mixer.
- 9) Incubate on ice for 10 minutes.
- 10) Centrifuge the cell lysate at $20,000 \times g$ for 10~20 minutes at 4 °C.
- 11) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. (Cell lysate)*1

*1 Please keep on ice until you use it for the research.

B. Suspension cells

- 1) Culture the cell lines expressing DYKDDDDK tagged recombinant proteins ($1 \times 10^6 \sim 10^7$).
- 2) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 3) Remove the supernatant.
- 4) After adding PBS (-) and resuspending the cell pellet with PBS (-), centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C. Repeat step 3) to 4).
- 5) Remove the supernatant.
- 6) After suspending the cell pellet with 1mL of PBS (-), transfer the cell suspension into a new 1.5 mL micro centrifugation tube.
- 7) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C and remove the supernatant.
- 8) Add 1 mL of cell lysis buffer to the cell pellet and suspend cells by pipetting or vortex mixer
- 9) Incubate on ice for 10 minutes.
- 10) Centrifuge the cell suspension at $20,000 \times g$ for 10~20 minutes at 4 °C.
- 11) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. (Cell lysate)*1

*1 Please keep on ice until you use it for the research.

2. Immunoprecipitation

- 1) Resuspend anti DYKDDDDK tag antibody magnetic beads by using vortex mixer.
- 2) Transfer 50 μL of beads slurry into a new 1.5 mL micro centrifuge tube.
- 3) Add 1mL of cell lysis buffer into 1.5 mL microcentrifuge tube including beads slurry.
- 4) Place this tube on magnetic stand for 30 seconds to 1 minute.
- 5) After confirming the separation of beads and supernatant, discard the supernatant carefully not to aspirate any beads. (Prewashed beads)

- 2/8 -

- 6) Add appropriate volume of cell lysate (200~1000 μ L) into 1.5 mL microcentrifuge tube including prewashed beads and resuspend the beads by vortex mixer.
- 7) Incubate with rotation for 4 hours at 4 °C.
- 8) Place this tube on magnetic separation stand after incubation.
- 9) After confirming the separation of beads and supernatant, discard the supernatant carefully not to aspirate any beads.
- 10) Add 1 mL of cell lysis buffer into 1.5 mL microcentrifuge tube including magnetic beads and resuspend by vortex mixer.
- 11) Place this tube on magnetic separation stand for 30 seconds to 1 minute.
- 12) After confirming the separation of beads and supernatant, discard the supernatant carefully not to aspirate any beads.
- 13) Repeat 10)-12) steps twice. (Antigen binding beads)

3. The elution of DYKDDDDK-tagged recombinant proteins

Please select an appropriate elution method in accordance with goal and property of protein.

A. The elution by using DYKDDDDK peptide (non-denaturing condition)

- 1) Prepare the DYKDDDDK peptide (Code No. 044-30951) stock solution (5 mg/mL final concentration). Dissolve 5 mg of DYKDDDDK peptide in 1 mL of 1 \times TBS (10 mmol/L Tris-HCl, pH 7.4, 150 mmol/L Sodium Chloride).
- 2) Prepare the diluted solution. Dilute DYKDDDDK peptide stock solution in 1 \times TBS at a concentration of 150-500 μ g/mL.
- 3) Add 50 μ L of above diluted solution of DYKDDDDK peptide into 1.5 mL microcentrifuge tube including antigen binding beads.
- 4) Incubate above sample at room temperature or 4 °C with shaking by using tube mixer etc.
- 5) After incubation, spin down the 1.5 mL microcentrifuge tube and place it on magnetic stand for 30 seconds to 1 minute.
- 6) After confirming the separation of beads and supernatant, transfer the supernatants to a new 1.5 mL microcentrifuge tube.
- 7) Store the collected solution on ice until further experiments. In the case of long storage, store at -20 °C.

B. The elution by using 0.1 mol/L Glycine-HCl, pH 3.0 (acidic condition)

- 1) Add 50 μ L of 0.1 mol/L Glycine-HCl (pH 3.0) into microcentrifuge tube including antigen binding beads.
- 2) Incubate the samples at room temperature for 10 minutes with shaking by using tube mixer etc.
- 3) After incubation, spin down the 1.5 mL microcentrifuge tube and place it on magnetic separation stand for 30 seconds to 1 minute.
- 4) After confirming the separation of beads and supernatant, transfer the supernatant to a new 1.5 mL microcentrifuge tube.
- 5) Neutralize with neutralization buffer. (Example : 0.5 μ L of 1M CAPS, pH 11)
- 6) Store the solution on ice until further experiments. In the case of long storage, store at -20 °C.

C. The elution by using SDS sample buffer (denaturing condition)

When 2 \times SDS sample buffer includes reducing agent such as 2-mercaptoethanol or DTT, heavy chain (50 kDa) and light chain (25 kDa) may dissociate from anti DYKDDDDK tag antibody and are eluted into a recovered sample.

- 1) Add 50 μ L of 2 \times SDS sample buffer (0.125 mol/L Tris-HCl, pH 6.8, 4w/v% SDS, 20w/v% Glycerol, 10 w/v% 2-mercaptoethanol or 100 mM DTT, 0.01w/v% Bromophenol Blue) into microcentrifuge tube including antigen binding bead and resuspend the samples by using vortex mixer.
- 2) Incubate the samples at 95 °C for 5 minutes.
- 3) Place it at room temperature for 5 minutes.
- 4) After incubation, spin down the 1.5 mL microcentrifuge tube and place it on magnetic separation stand for 30 seconds to 1 minute.
- 5) After confirming the separation of beads and supernatant, transfer the supernatants to a new 1.5 mL microcentrifuge tube.
- 6) Store the collected solution on ice until further experiments. In the case of long storage, store at -20 °C.

D. The elution by using 4w/v%SDS solution (denaturing condition)

This elution method inhibits co-elution of immunoglobulin heavy chain (50kDa) and light chain (25kDa) into eluate.

- 1) Prepare 4w/v%SDS solution and add 50 μ L 4w/v% SDS solution into microcentrifuge tube including antigen binding beads.
- 2) Incubate the samples at room temperature for 10 minutes with shaking by using tube mixer etc.
- 3) After incubation, spin down the 1.5 mL microcentrifuge tube and place it on magnetic separation stand for 30 seconds to 1 minute.
- 4) After confirming the separation of beads and supernatant, transfer the supernatants to a new 1.5 mL microcentrifuge tube.
- 5) Store the collected solution on ice until further experiments. In the case of long storage, store at -20 °C.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : +1-804-271-7677
 Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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