(90 × 210mm Size)

FUJIFILM

Wako

Code No. 136-18341 (2 mL)

MagCapture™ Tamavidin® 2-REV
MagCapture™ タマビジン® 2-REV

MagCapture™ Tamavidin® 2-REV is the slurry of Tamavidin® 2-REV immobilization magnetic beads.

Tamavidin® 2-REV is a genetic mutant of Tamavidin® 2 which is identified as a novel avidin-like protein from mushroom. Tamavidin® 2-REV is a useful tool for the purification of biotinylated molecules such as nucleic acids, chemicals, antibodies and enzymes by using competitive biotin elution method.

[Formulation]

10 mg/mL Tamavidin® 2-REV magnetic beads, 1 × TBS (pH 7.4), 50 w/v% glycerol, 0.05 w/v% sodium azide.

[Source/ Tamavidin® 2-REV]

E. coli expressed Tamavidin® 2-REV.

[Assay]

Please optimize the condition for your samples and experiments.

[Additional materials]

Magnetic separation stand, Centrifuge, Rotator, Vortex mixer, Tube mixer, 1.5 mL centrifugation tube, Washing buffer (see below).

Recommended washing buffer :

- PBS(-)
- PBS-T
- TBS-T
- Cell Lysis Buffer M (Code No. 038-21141)
20 mmol/L Tris-HCl, pH 7.4, 200 mmol/L Sodium Chloride,
2.5 mmol/L Magnesium Chloride, 0.05 w/v% NP-40 substitute.

【Protocol】

Immobilization protocol

1. Resuspend the product.
2. Add 50 μ L of beads suspension into a new 1.5 mL centrifugation tube.
3. Wash beads with an appropriate washing buffer for your experiment^{※1}.
4. After washing the magnetic beads, add the sample solution including biotinylated molecules.
5. Resuspend well with rotator for 1 hour at 4 °C.
6. Place the tube on a magnetic stand for a while.
7. Remove the supernatant with pipette.
8. Wash beads with an appropriate washing buffer for three times.
9. Remove the supernatant with pipette [A].

- 1/4 -

Elution protocol

1. Prepare the excess biotin solution ; 2 mmol/L biotin/PBS.
 2. Add 15 μ L excess biotin solution to [A].
 3. Mix well and incubate for 15 minutes at 4 °C .
 4. Place the tube on a magnetic stand.
 5. Recover the supernatant.
- * When the recovery amount is low, repeat steps 2 to 5 once again.

※ 1 Washing protocol

1. Set the tube containing 50 μ L of beads suspension on a magnetic stand for several minutes.
2. Remove the supernatant by pipetting on a magnetic stand.
3. Add 1 mL of an appropriate washing buffer and resuspend by vortex mixer.
4. Spin down and place the tube on a magnetic stand for several minutes.
5. Repeat steps 2 to 4 proper times.
6. Remove the supernatant by pipetting on a magnetic stand.

[Storage]

Store at -20 °C.

[References]

- 1) Takakura, Y. *et al.*, "Tamavidin 2-REV : An engineered tamavidin with reversible biotin-binding capability.", *J. Biotechnol.*, **164**(1), 19 (2013).
- 2) Ukekawa, R. *et al.*, "A useful new tool for isolation of long non-coding RNA-binding proteins.", *Medical Science Digest*, **40**(7), 355 (2014).

[License]

Tamavidin® is the registered trademark of Japan Tobacco Inc.. FUJIFILM Wako Pure Chemical Corporation. has sold Tamavidin® under the sales license from Japan Tobacco Inc..

[Package]

2 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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