(90×210mm Size)

FUJIFILM

Wako

Code No. 145-09681 (10 mL, NET 5 mL)
141-09683 (100 mL, NET 50 mL)
149-09684 (2 mL, NET 1 mL)

for genetic research

Ni-NTA Agarose HP

Ni-NTA Agarose is designed for affinity purification of 6× His fusion protein. Nickel is used for a wide range of protein purification, because it has high affinity for 6× His fusion proteins. Ni-NTA Agarose consists of a metal ion and a chelate that is coated with 6% cross-linked agarose coupled to Ni^{2+} ions via a nitrilotriacetic acid (NTA) ligand. Therefore this product can be used for purification of denatured proteins. In addition, Fast Protein Liquid Chromatography (FPLC), is a method of liquid chromatography often used to purify mixtures of proteins. This product is available to these purifications also.

This product is for laboratory use only; use in any such application is the responsibility of users.

[Matrix] Highly cross-linked 6% agarose

[Ligand] Nitrilotriacetic acid (NTA)

[Metal ion capacity] $\geq 15 \mu\text{mol Ni}^{2+}/\text{mL gel}$

[Bead size] 50~150 μm

[Storage buffer] 50% slurry in 20% ethanol

[Maximum total pressure] 0.6MPa

[Maximum liner flow rate] 1,800cm/h or 60 mL/min

[Recommended buffers]

For native conditions :

- Equilibration buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 10 mmol/L Imidazole, pH 8.0
- Wash buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 20 mmol/L Imidazole, pH 8.0
- Elution buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 250 mmol/L Imidazole, pH 8.0

For denaturing conditions :

- Equilibration buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 10 mmol/L Imidazole, 8 mol/L Urea, pH 8.0
- Wash buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 20 mmol/L Imidazole, 8 mol/L Urea, pH 8.0
- Elution buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 250 mmol/L Imidazole, 8 mol/L Urea, pH 8.0

Note : It is important to avoid agents such as EDTA or citrate in the above buffers.

- 1/6 -

[Protocol]

1. Procedure for batch purification of His-tagged proteins

1-A Resin equilibration

- 1) Thoroughly suspend Ni-NTA resin and immediately transfer an appropriate amount of Ni-NTA resin to a tube.
- 2) Centrifuge at $500 \times g$ for 5 minutes.
- 3) Remove the supernatant. Be careful not to aspirate the resin.
- 4) Add 10 resin-bed volumes of Equilibration buffer and suspend with a vortex mixer.
- 5) Centrifuge at $500 \times g$ for 5 minutes.
- 6) Remove the supernatant. Be careful not to aspirate the resin.

1-B Sample application

- 1) Add the prepared protein extract.
- 2) Incubate with rotation for 30~60 minutes.
- 3) Centrifuge at $500 \times g$ for 5 minutes.
- 4) Remove the supernatant. Be careful not to aspirate the resin.

1-C Washing

- 1) Add 10 resin-bed volumes of Wash buffer and suspend with a vortex mixer.
- 2) Centrifuge at $500 \times g$ for 5 minutes.
- 3) Remove the supernatant. Be careful not to aspirate the resin.
- 4) Repeat 1)~3) step 2 times.

1-D Elution

- 1) Elute the His-tagged proteins using a resin-bed volume of Elution Buffer.
- 2) Incubate with rotation for 10 minutes.
- 3) Centrifuge at $500 \times g$ for 5 minutes.
- 4) Transfer the supernatant into a new tube. Be careful not to aspirate the resin.
- 5) Repeat 1)~4) step 2 times.
- 6) Store the collected supernatant at 2~8 °C.

2. Procedure for gravity purification of His-tagged proteins

2-A Resin equilibration

- 1) Thoroughly suspend Ni-NTA resin. Immediately add an appropriate amount of the Ni-NTA resin to an appropriate column.
- 2) Remove the storage buffer up to the vicinity of the resin surface by gravity flow.
- 3) Add 5 resin-bed volumes of Equilibration buffer.
- 4) Remove Equilibration buffer up to the vicinity of the resin surface by gravity flow.
- 5) Repeat 3)~4) step 2 times.

2-B Sample application

- 1) Add the prepared protein extract.
- 2) Remove the sample buffer up to the vicinity of the resin surface by gravity flow as much as possible.

- 2/6 -

2-C Washing

- 1) Add 10 resin-bed volumes of Wash buffer to the resin.
- 2) Remove the Wash buffer up to the vicinity of the resin surface by gravity flow.
- 3) Repeat 1)~2) step 2 times.

2-D Elution

- 1) Add a resin-bed volume of Elution buffer and elute the His-tagged proteins by gravity flow.
- 2) Store the fraction on ice.
- 3) Repeat 1)~2) step 2 times.
- 4) Store the collected fraction at 2~8 °C.

[Compatibility of reagents]

Ni-NTA Agarose HP is compatible with the following compounds up to the concentrations given.

Reducing agents	30 mmol/L Reduced glutathione
	20 mmol/L β -Mercaptoethanol
	10 mmol/L DTT
	10 mmol/L DTE
	0.3% SDS
Denaturing agents	8 mol/L Urea
	6 mol/L Guanidine Hydrochloride
Detergents	2% Triton® X-100
	2% Tween® 20
Other additives	50% Glycerol
	1 mmol/L EDTA
	20% Ethanol

Triton is a registered trademark of Dow Chemical Company.
Tween is a registered trademark of ICI Americas, Inc.

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