FUJIFILM

Wako

Code No. 143-09763 (10 mL, Net 5 mL)
141-09764 (100 mL, Net 50 mL)

for genetic research
Ni-NTA Agarose

Ni-NTA Agarose is designed for affinity purification of 6×His fusion protein. Ni-NTA Agarose consists of a metal ion and a chelate that is coated with 6% cross-linked agarose coupled to Ni²⁺ ions via a nitrilotriacetic acid (NTA) ligand. 6×Histidine tag binds to Ni²⁺ ions. Protein purification with Ni-NTA agarose is simple and this affinity purification is used for a wide range of proteins : soluble proteins, insoluble proteins, and membrane proteins.

[Matrix]	Cross-linked 6% agarose
[Ligand]	Nitrilotriacetic acid (NTA)
[Protein binding capacity]	> 50mg/mL gel*
[Bead size]	50 ~ 150 μm
[Storage buffer]	50% slurry in 30% ethanol

*Binding capacity is different from each target protein

[Recommended buffers]

For native conditions :

- Equilibration buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 10 mmol/L Imidazole, pH 8.0
- Wash buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 20 mmol/L Imidazole, pH 8.0
- Elution buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 250 mmol/L Imidazole, pH 8.0

For denaturing conditions :

- Equilibration buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 10 mmol/L Imidazole, 8 mol/L Urea, pH 8.0
- Wash buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 20 mmol/L Imidazole, 8 mol/L Urea, pH 8.0
- Elution buffer : 50 mmol/L Sodium Dihydrogenphosphate, 300 mmol/L Sodium Chloride, 250 mmol/L Imidazole, 8 mol/L Urea, pH 8.0

Note : It is important to avoid agents such as EDTA or citrate in the above buffers.

[Protocol]

1. Procedure for batch purification of His-tagged proteins

1-A Resin equilibration

- 1) Thoroughly suspend Ni-NTA resin and immediately transfer an appropriate amount of Ni-NTA resin to a tube.
- 2) Centrifuge at 500×g for 5 minutes.
- 3) Remove the supernatant. Be careful not to aspirate the resin.

- 4) Add 10 resin-bed volumes of Equilibration buffer and suspend with a vortex mixer.
- 5) Centrifuge at 500×g for 5 minutes.
- 6) Remove the supernatant. Be careful not to aspirate the resin.

1-B Sample application

- 1) Add the prepared protein extract.
- 2) Incubate with rotation for 30 ~ 60 minutes
- 3) Centrifuge at 500×g for 5 minutes.
- 4) Remove the supernatant. Be careful not to aspirate the resin.

1-C Washing

- 1) Add 10 resin-bed volumes of Wash buffer and suspend with a vortex mixer.
- 2) Centrifuge at 500×g for 5 minutes
- 3) Remove the supernatant. Be careful not to aspirate the resin.
- 4) Repeat 1) ~ 3) step 2 times.

1-D Elution

- 1) Elute the His-tagged proteins using a resin-bed volume of Elution Buffer.
- 2) Incubate with rotation for 10 minutes.
- 3) Centrifuge at 500×g for 5 minutes.
- 4) Transfer the supernatant into a new tube. Be careful not to aspirate the resin.
- 5) Repeat 1) ~ 4) step 2 times.
- 6) Store the collected supernatant at 2 ~ 8°C.

2. Procedure for gravity purification of His-tagged proteins

2-A Resin equilibration

- 1) Thoroughly suspend Ni-NTA resin. Immediately add an appropriate amount of the Ni-NTA resin to an appropriate column.
- 2) Remove the storage buffer up to the vicinity of the resin surface by gravity flow.
- 3) Add 5 resin-bed volumes of Equilibration buffer.
- 4) Remove Equilibration buffer up to the vicinity of the resin surface by gravity flow.
- 5) Repeat 3) ~ 4) step 2 times.

2-B Sample application

- 1) Add the prepared protein extract.
- 2) Remove the sample buffer up to the vicinity of the resin surface by gravity flow as much as possible.

2-C Washing

- 1) Add 10 resin-bed volumes of Wash buffer to the resin.
- 2) Remove the Wash buffer up to the vicinity of the resin surface by gravity flow.
- 3) Repeat 1) ~ 2) step 2 times.

2-D Elution

- 1) Add a resin-bed volume of Elution buffer and elute the His-tagged proteins by gravity flow.
- 2) Store the fraction on ice.

- 3) Repeat 1) ~ 2) step 2 times.
4) Store the collected fraction at 2 ~ 8°C.

Note : If it is necessary to prepare a fresh sample, please run the protocol at 2-8°C.

【Compatibility of reagents】

Ni-NTA Agarose HP is compatible with the following compounds up to the concentrations given.

Reducing agents	30 mmol/L Reduced glutathione
	20 mmol/L β -Mercaptoethanol
	10 mmol/L DTT
	10 mmol/L DTE
	0.3% SDS
Denaturing agents	8 mol/L Urea
	6 mol/L Guanidine Hydrochloride
Detergents	2% Triton [®] X-100
	2% Tween [®] 20
Other additives	50% Glycerol
	1 mmol/L EDTA
	20% Ethanol

Triton is a registered trademark of Dow Chemical Company.
Tween is a registered trademark of ICI Americas, Inc.

【Storage】

2~10°C

【Packages】

Code No.	Packages
143-09763	10 mL (Net 5 mL)
141-09764	100 mL (Net 50 mL)

【Related products】

Code No.	Product Name	Packages
146-09731	Ni-NTA Cartridge	1 Cartridge (5 mL)
142-09733		1 Cartridge (5 mL) × 5
149-09684	Ni-NTA Agarose HP	2 mL (NET 1 mL)
145-09681		10 mL (NET 5 mL)
141-09683		100 mL (NET 50 mL)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
http://www.wako-chem.co.jp

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
http://www.wakousa.com

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