For Tissue Optical Clearing Reagent SCALEVIEW-S Trial Kit

[Introduction]

Dr. Atsushi Miyawaki et al. developed a sorbitol-based optical clearing method, *Scale S*.

SCALEVIEW-S is a high-performance clearing method that blessed with effective clearing and tissue/signal preservation. It is also noteworthy that this method is able to preserve highly fluorescence signal of fluorescent proteins and ultrastructure of cell membranes.

SCALEVIEW-S is a simple and reproducible method of accurate visualization of biological tissue.

[Kit contents]

This kit includes 6 components.

- | | |
|-------------------|-------------------|
| (1) SCALEVIEW-S0 | 100 ml × 1 bottle |
| (2) SCALEVIEW-S1 | 100 ml × 1 bottle |
| (3) SCALEVIEW-S2 | 100 ml × 1 bottle |
| (4) SCALEVIEW-S3 | 100 ml × 1 bottle |
| (5) SCALEVIEW-S4 | 100 ml × 1 bottle |
| (6) SCALEVIEW-SMt | 100 ml × 1 bottle |

Note : This kit needs an optional deScale solution (Code No. 041-34425).

[Storage]

Store in the dark at 2-10°C

[Additional required materials]

Reagents

- 1) Paraformaldehyde (Code No. 160-16061)
- 2) 1 × PBS (-) (Code No. 164-25511)
- 3) Agarose
- 4) deScale Solution (Code No. 041-34425)

Equipments :

- 1) 5 ml/15 ml/30 ml/50 ml conical centrifuge tube
- 2) Culture dish (100 mm, 60 mm, and 35 mm diameter) (for each sample)
- 3) Seesaw shaker
- 4) Fluorescence microscope and Recommend Objective lenses for SCALEVIEW-S treated samples.
 1. Multiphoton microscope OLYMPUS, model : FVMPE-RS
 - Multiphoton dedicated objectives :
 - OLYMPUS, model : XLPLN10XSVM : 10x, NA=0.6, WD=8 mm
 - OLYMPUS, model : XLSLPLN25XGMP : 25x, NA=1.0, WD=8 mm
 - 2. Confocal microscope : OLYMPUS, model : FV3000
 - Objectives :
 - OLYMPUS, model : UPLSAP010X2 : 10x, NA=0.4, WD=3.1 mm
 - OLYMPUS, model : UCFLFLN20X : 20x, NA=0.7, WD=0.8-1.8 mm

[Procedure]

(Fixation)

A case for making mouse brain transparent according to reference (1).

- 1) Transcardially perfuse an anesthetized mouse with 4% paraformaldehyde (PFA)/PBS (pH 7.6 ~ 7.8).
- 2) Remove the whole brain and subject it to post-fixation in 4% PFA/PBS at 4°C for 72 hrs.
- 3) Wash the samples in PBS.
- 4) Prepare vibratome sections (optional, 0.2-3 mm thick)

(Clearing) (1-2 mm brain samples)

- 5) Transfer the sample into 5 ml of SCALEVIEW-S0 solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 30 min.
- 6) Transfer the sample into 5 ml of SCALEVIEW-S1 solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 30 min.
- 7) Transfer the sample into 5 ml of SCALEVIEW-S2 solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 30 min.
- 8) Transfer the sample into 5 ml of SCALEVIEW-S3 solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 30 min.
- 9) Transfer the sample into 5 ml of deScale Solution and incubate it with gentle shaking (10 rpm) at 4°C for 3 hours × 2 times.
- 10) Transfer the sample into 5 ml of SCALEVIEW-S4 solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 12-24 hours.
- 11) The transparency can be evaluated by eye at this stage. If the sample is successfully cleared, the brain sample should look amber under a light source.
- 12) Transfer the sample into 5 ml of SCALEVIEW-SMt solution and incubate it with gentle shaking (10 rpm) at 35-40°C for 1 hour.


Figure 1. Transmittance images of 1mm brain slice before and after clearing with SCALEVIEW-S.

(Imaging)

- 13) Mount onto an imaging chamber. Use appropriate amount of SCALEVIEW-SMt solution to mount the sample. Observe the SCALEVIEW-S cleared samples using fluorescence microscopy.


Figure 2. Confocal imaging of adult hemisphere brain (Thy1-YFP-H mouse).

Note :

- For thicker brain slices (>2 mm) or hemisphere mouse brain whole mouse brain, incubation time in each step can become longer and SCALEVIEW-S solutions in each step can increase amount (e.g., hemisphere brain 25 ml, whole brain 40 ml each).
- SCALEVIEW-SMt mounted samples will not be solidified. If necessary for each sample, 1.5-3.0% melting agarose can be added to minimize movement artifacts during imaging.

[References]

- 1) Hama, H. et al. : *Protocol Exchange* (2016).
- 2) Hama, H. et al. : *Nat Neurosci*, **18**, 1518 (2015).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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