Code No. 294-80201(for 100mL)
290-80203(for 1L)

ScreenFectTM UP-293

ScreenFectTM UP-293 is a transfection reagent kit which is highly efficient for production of antibodies and proteins in Expi293FTM or FreeStyleTM 293-F cells. ScreenFectTM UP-293 Booster, which is added to transfected cells, enhances remarkably protein expression level in cells.

【Kit contents】

	For 100mL	For 1L
ScreenFect TM UP-293 Transfection Reagent	0.2mL	2x1mL
ScreenFect TM UP-293 Dilution Buffer	6.7mL	67mL
ScreenFect TM UP-293 Booster	0.5mL	5mL

【Storage】

Store at 2-10°C. Do not freeze.

【Additional required reagents】

	For Expi293 TM Expression System	For FreeStyle TM 293 Expression System
Cells	Expi293F TM Cells	FreeStyle TM 293-F Cells
Medium	Expi293 TM Expression Medium	FreeStyle TM 293 Expression Medium
Dilution Buffer	-	Opti-MEM TM I Reduced Serum Medium

【Attention】

- Please be careful of bacteria contamination under experiment.
- After storage of several months without using the reagent, slight precipitation in ScreenFectTM UP-293 Transfection Reagent might occurs. Vortex it well before use.
- Optimization study of amounts of ScreenFectTM UP-293 Transfection Reagent, plasmid DNA and cell culture time is required for successful protein expression.

Transfection protocol for 30 ml culture volume

①

One day prior to transfection, prepare culture cells to be appropriate cell numbers for transfection using a 125mL or 250mL of Erlenmeyer shake flask.
Note: Recommended seeding cell density is >0.8x10⁶ viable cells/mL.

②

Prepare cell suspension

Just before transfection, prepare 28mL of cell suspension with the cell densities below in a 125 mL of Erlenmeyer shake flask.

【For Expi293F™ Cells】

2.7x10⁶ viable cells/mL medium (7.5x10⁷ viable cells/28mL of medium)

【For FreeStyle™ 293-F Cells】

1.1x10⁶ viable cells/mL medium (3.0x10⁷ viable cells/28mL of medium)


Note: cell viability must not be < 90 %.

③

Mix ScreenFect™UP-293 Transfection Reagent well with a vortex mixer.

④

Prepare a diluted transfection reagent


【For Expi293F™ Cells】

Add 60µL of ScreenFect™UP-293 Transfection Reagent into 940µL of ScreenFect™UP-293 Dilution Buffer in a sterile microfuge tube and mix it using pipette strokes.

【For FreeStyle™ 293-F Cells】

Add 30µL of ScreenFect™UP-293 Transfection Reagent into 970µL of Opti-MEM™ I Reduced Serum Medium in a sterile microfuge tube and mix it using pipette strokes.

⑤

Prepare a diluted DNA solution


【For Expi293F™ Cells】

Add 30µg of plasmid DNA with a density of 1µg/µL into 970µL of ScreenFect™UP-293 Dilution Buffer in a sterile microfuge tube and mix it using pipette strokes.

【For FreeStyle™ 293-F Cells】

Add 30µg of plasmid DNA with a density of 1µg/µL into 970µL of Opti-MEM™ I Reduced Serum Medium in a sterile microfuge tube and mix it using pipette strokes.

⑥


Prepare a lipoplex solution

Mix the diluted transfection reagent and the diluted DNA solution with rapid pipette strokes and spin down it. Incubate it for 5-20minutes at room temperature.

Note: Do not vortex the mixed solution.

⑦

Add the 2mL of the lipoplex solution into the cell suspension in the Erlenmeyer Shaker Flask with slowly swirling.

Note: The final volume should be approximately 30mL in the flask.

⑧

Incubate the cells at 37°C(8% CO₂), rotating at 125rpm on an orbital shaker for 16hours.

⑨

Add 150µL of ScreenFect™UP-293 Booster to the cells.

⑩

Continue to incubate the cells at 37°C(8% CO₂), rotating at 125rpm on an orbital shaker .

Collect cells or culture supernatants containing target proteins in 48-120 hours after transfection.

Note: The incubation period required for successful transfection depends on the target protein and cell viability.

Expi293F™, Freestyle™ and Opti-MEM™ are registered trademarks of Thermo Fisher Scientific.

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	For Expi293 TM Expression System	For FreeStyle TM 293 Expression System
細胞	Expi293F TM Cells	FreeStyle TM 293-F Cells
培地	Expi293 TM Expression Medium	FreeStyle TM 293 Expression Medium
希釈バッファー	-	Opti-MEM TM I Reduced Serum Medium

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