## For Tissue Optical Clearing Reagent CUBIC Trial Kit

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Dr. Etsuo A. Susaki *et al.* developed a simple, efficient, and scalable brain / body clearing method and imaging / computational analysis pipeline, CUBIC (clear, unobstructed brain/body imaging cocktails and computational analysis). CUBIC offers a fluorescent-preserving, high-performance and device-free tissue clearing method based on hydrophilic Solutions (ScaleCUBIC Solutions). It enables reproducible whole-organ and whole-body clearing and rapid 3D imaging with LSFM. CUBIC also provides processing and analysis of 3D images for extracting biological information.

Therefore, CUBIC presents a platform of whole-organ/body imaging with single-cell resolution and image informatics, which enables a wide range of users to perform experiments targeting cellular and organ layers with multiple samples.

### 【Kit contents】

This kit includes 4 components.

(1) ScaleCUBIC-1 Solution	125 mL×2 bottles
(2) ScaleCUBIC-2 Solution	125 mL×2 bottles
(3) Mounting Solution 1	40 mL×1 bottle
(4) Mounting Solution 2	40 mL×1 bottle

### 【Storage】

Stored in the dark

### 【Additional required materials】

Reagents :

- 1) Paraformaldehyde (Code No.160-16061)
- 2) 1×PBS (-) (Code No. 164-25511)

Equipment :

- 1) 5 mL / 15 mL / 30 mL / 50 mL conical centrifuge tubes
- 2) Seesaw shaker

### 【Procedure】

(Solution preparation)

1. 50% ScaleCUBIC-1 Solution (1:1 mixture)  
Dissolve ScaleCUBIC-1 Solution in distilled water.
2. 50% ScaleCUBIC-2 Solution (1:1 mixture)  
Dissolve ScaleCUBIC-2 Solution in distilled water.

(Fixation)

A case for making mouse organs transparent according to reference 1), 2) and 3).

- 1) Transcardially perfuse an anesthetized mouse with 10 mL ~ of cold PBS with Heparin (10 U/mL) followed by 20 mL ~ of cold 4% paraformaldehyde (PFA)/PBS (pH 7.4).

- 2) Dissect the organs and subject them to post-fixation in 4% PFA/PBS at 4°C for 18-24 hrs.
- 3) Wash the samples in PBS.
- 4) Prepare thick tissue sections (optional, 2 mm thick for example).

### 【Clearing】 (2 mm section samples)

Delipidation and decolorization

- 1) Transfer the sample into 2.5 mL of 50% ScaleCUBIC-1 Solution in 5 mL tube and incubate it with gentle shaking (30 rpm) at RT for 6-24 hrs.
- 2) Transfer the sample into 2.5 mL of ScaleCUBIC-1 Solution and incubate it with gentle shaking (30 rpm) at 37°C for 2 days.  
If needed, transfer the sample into 2.5 mL of fresh ScaleCUBIC-1 Solution and continue to incubate it with gentle shaking (30 rpm) at 37°C for an additional few days until the samples sufficiently cleared (dependent on tissue and sample size). Replace the ScaleCUBIC-1 Solution every 2 days.
- 3) Transfer the sample into 25 mL of PBS in 30 or 50 mL tube and incubate it with gentle shaking (30 rpm) at RT to 37°C at least for 2 hrs ×3 times (dependent on sample size).
- 4) Transfer the sample into 2.5 mL of 50% ScaleCUBIC-2 Solution in 5 mL tube and incubate it with gentle shaking (30 rpm) at RT for 1 day.  
\*The sample will become swollen.

RI match

- 5) Transfer the sample into 2.5 mL of ScaleCUBIC-2 Solution in 5 mL tube and incubate it with gentle shaking (30 rpm) at RT to 37°C for 1-2 days. \*Because ScaleCUBIC-2 Solution is highly viscous, manual mixing with inversion of the tube will be occasionally needed several times.

\*\*Sample size will be recovered to the original size.


Figure 1. Transmittance images of mouse 2 mm slices before and after clearing with ScaleCUBIC Solutions.

### 【Clearing】 (whole organ samples)

Delipidation and decolorization

- 1) Transfer the sample into 12-15 mL of 50% ScaleCUBIC-1 Solution in 30 mL tube and incubate it with gentle shaking (30 rpm) at RT for 1 day.
- 2) Transfer the sample into 12-15 mL of ScaleCUBIC-1 Solution in 30 mL tube and incubate it with gentle shaking (30 rpm) at 37°C for 2 days.
- 3) Transfer the sample into 12-15 mL of fresh ScaleCUBIC-1 Solution in 30 mL tube and continue to incubate it with gentle shaking (30 rpm)

at 37°C for additional 6-8 days (dependent on tissue and sample size).  
Replace the ScaleCUBIC-1 Solution every 2 days.

- 4) Transfer the sample into ~ 50 mL of PBS in 50 mL tube and incubate it with gentle shaking (30 rpm) at RT to 37°C at least for 2 hrs × 3 times. It is recommended to wash the sample once for 2 hrs, once overnight and again once for 2 hrs.
- 5) Transfer the sample into 12-15 mL of 50% ScaleCUBIC-2 Solution in 30 mL tube at RT for 1 day. \*The sample will become swollen.

RI match

- 6) Transfer the sample into 12-15 mL of ScaleCUBIC-2 Solution in 30 mL tube at RT to 37°C for 2 days. \*Because ScaleCUBIC-2 Solution is highly viscous, manual mixing with inversion of the tube will be occasionally needed several times.

\*\*Sample size will be recovered to the original size.

Note :

- For thicker slices (>2 mm) or whole organs, incubation time in each step can become longer and ScaleCUBIC Solutions in each step can increase the amount.


Figure 2. Transmittance images of mouse whole organs before and after clearing with ScaleCUBIC Solutions.

(Imaging)

- 1) The RI-matched sample is immersed in a mixture of Mounting Solution 1 and Mounting Solution 2. Observe the CUBIC-cleared samples using fluorescence microscopy.

<Note>

The optimal ratio of Mounting Solution 1 and Mounting Solution 2 for transparency observation depending on the sample, so it is recommended to optimize the ratio according to the sample. Typically, 2:8 or 3:7 mixture of Mounting Solution 1 and 2 (RI ~ 1.49) can be used.


Figure 3. LSFM imaging of CAG-EGFP Tg Mouse (8 w Male) whole brain, whole heart and dissected liver.

## 【References】

- 1) E. A. Susaki. *et al.* : *Cell*, **157**, 726 (2014).
- 2) K. Tainaka. *et al.* : *Cell*, **159**, 911 (2014).
- 3) E. A. Susaki. *et al.* : *Nature Protocols*, **10**, 1709 (2015).

## 【Related Products】

Code No.	Description	Size
160-16061	Paraformaldehyde	100 g
164-25511	1 × PBS (-)	5 L

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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