[LBIS Mouse IFN- $\gamma$  ELISA Kit]

Cat # 630-44701

Please, read this instruction carefully before use.

This kit is manufactured by FUJIFILM Wako Pure Chemical Corporation. Use only the current version of the Instruction Manual enclosed with the kit!

**Intended use**

Interferon-gamma (IFN- $\gamma$ ) is the only type II interferon consisting of 143 amino acid residues (human), which is structurally different from type I interferons such as IFN- $\alpha/\beta$ . IFN- $\gamma$  is secreted predominantly by T cells and NK cells and functions as a homodimer. It is known to have various biological activities, including antiviral effect, cell growth inhibition, anti-tumor effect, alteration of macrophage activity, promotion of NK cell activity, immunomodulatory effects, modulation of the immune response, modulation of the differentiation process, and is attracting attention in various research fields.

This kit is a sandwich ELISA system for quantitative measurement of human IFN- $\gamma$ . It can measure IFN- $\gamma$  in mouse serum (plasma) specifically and with high sensitivity. This product is intended for research use only. Not for use in diagnostics procedures.

**Assay principle**

In FUJIFILM Wako Pure Chemical Corporation's LBIS Mouse IFN- $\gamma$  ELISA Kit, standards or samples are incubated in polyclonal antibody-coated wells to capture IFN- $\gamma$ . After 2 hours' incubation and washing, biotinylated anti-antibody is added and incubated further for 1 hour to bind with captured. After washing, HRP (horse radish peroxidase)-conjugated streptavidin is added, and incubated for 30 minutes. After washing, bound HRP-conjugated streptavidin is reacted with a chromogen (TMB) for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to concentration. The standard curve is prepared by plotting absorbance against standard IFN- $\gamma$  concentrations. concentrations in unknown samples are determined using this standard curve.

**Reagents supplied**

Components	State	Amount
(A) Anti-IFN- $\gamma$ antibody coated plate	Use after washing.	96 wells/1 plate
(B) Standard mouse IFN- $\gamma$	Freeze-dried. Use after reconstitution.	1 vial
(C) Buffer solution	Ready for use.	60 mL/1 bottle
(D) Biotinylated anti-IFN- $\gamma$ antibody	Freeze-dried. Use after reconstitution.	1 vial
(E) HRP-conjugated streptavidin	Concentrated. Use after dilution.	100 $\mu$ L/1 vial
(F) Chromogen (TMB)	Ready for use.	12 mL/1 bottle
(H) Stop solution <b>Be careful!</b>	Ready for use.	12 mL/1 bottle
( I ) Wash stock solution (10 $\times$ )	Concentrated. Use after dilution.	100 mL/1 bottle
Plate seal	—	4 sheets
Instruction Manual	—	1 copy

**Storage and expiration**

When the complete kit is stored at 2 °C - 8 °C (Do not freeze), the kit is stable until the expiration date shown on the label on the container. Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

**Equipments or supplies required but not supplied**  Use as a check box

Deionized water (or Distilled water)  Test tubes for preparation of standard solution series.  Glassware for dilution of Wash stock solution (10 $\times$ ) (a graduated cylinder, a bottle)  Pipettes (disposable tip type). One should be able to deliver 50  $\mu$ L precisely, and another for 100  $\mu$ L- 1000  $\mu$ L.  Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 100  $\mu$ L.  Paper towel to remove washing buffer remaining

LBIS Mouse IFN- $\gamma$  ELISA Kit

in wells. □A vortex-type mixer. □A shaker for 96 well-plate (600 rpm -1200 rpm) □An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. □A 96 well-plate reader (450 nm  $\pm$ 10 nm, 620 nm: 600 nm - 650 nm) □Software for data analysis.

**Preparation of reagents**

- ◆Bring all reagents of the kit to room temperature (20 °C - 25 °C) before use.
- ◆Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

**[Concentrated reagents]**

**[(B) Standard mouse IFN- $\gamma$  (Freeze-dried)]**


Reconstitute the (B) standard with volume\* of deionized water (or distilled water) described in a separate sheet to prepare the Original standard solution (5 ng/mL). And then prepare dilution series of standard solutions as shown below using the (C) Buffer solution brought to the room temperature.

\*Because the volume of deionized water (or distilled water) to be added to the freeze-dried standard differs depending on the lot, see the specified volume in a separate sheet.

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer solution	Concentration
Original standard solution 50 $\mu$ L	450 $\mu$ L	500 pg/mL
500 pg/mL solution 100 $\mu$ L	150 $\mu$ L	200 pg/mL
200 pg/mL solution 100 $\mu$ L	150 $\mu$ L	80.0 pg/mL
80.0 pg/mL solution 100 $\mu$ L	150 $\mu$ L	32.0 pg/mL
32.0 pg/mL solution 100 $\mu$ L	150 $\mu$ L	12.8 pg/mL
12.8 pg/mL solution 100 $\mu$ L	150 $\mu$ L	5.12 pg/mL
5.12 pg/mL solution 100 $\mu$ L	150 $\mu$ L	2.05 pg/mL
Blank	200 $\mu$ L	0

} Standard Curve

## LBIS Mouse IFN- $\gamma$ ELISA Kit

### [(D) Biotinylated anti-IFN- $\gamma$ antibody]

Add 100  $\mu$ L of deionized water (or distilled water) of room temperature to (D) Biotinylated anti-IFN- $\gamma$  antibody (Freeze-dried) ( Original (D) solution). Prepare working solution by dilution of original (D) solution with the (C) Buffer solution to **1:100**. 10 mL of the diluted solution is enough for 96 wells.

### [(E) HRP-conjugated streptavidin]

Prepare working solution by dilution of (E) with the (C) Buffer solution to **1:100**. 10 mL of the diluted solution is enough for 96 wells.

### [(I) Wash stock solution (10 $\times$ )]

Dilute 1 volume of the concentrated (I) Wash stock solution (10 $\times$ ) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example: 100 mL of concentrated (I) Washing buffer (10 $\times$ ) and 900 mL of deionized water (or distilled water).

### **【Storage and stability】**

#### [(A) Anti-IFN- $\gamma$ antibody coated plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used for well-plate container and store at 2 °C - 8 °C. The strip will be stable until expiration date.

#### [(B) Standard mouse IFN- $\gamma$ (Freeze-dried)]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted standard solutions after use.

#### [(C) Buffer solution] & [(F) Chromogen (TMB)]

Use only volume you need for your assay. Remaining reagents should be stored at 2 °C - 8 °C closing the cap tightly. It maintains stability until expiration date. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

#### [(D) Biotinylated anti-IFN- $\gamma$ antibody]

The reconstituted (D) solution should be used within 2 weeks if stored in a refrigerator. Dispose the remaining diluted working solution.

#### [(E) HRP-conjugated streptavidin]

Remaining working solution (already diluted) should be disposed. The rest of the undiluted solution (unused) : close the cap tightly and store at 2 °C - 8 °C. It is stable until expiration date. Once opened, we recommend using as soon as possible to avoid influence by environmental condition.

#### [(H) Stop solution]

Close the cap tightly and store at 2 °C - 8 °C. It maintains stability until expiration date.

#### [(I) Wash stock solution (10 $\times$ )]

The rest of undiluted solution (unused) : close the cap tightly and store at 2 °C - 8 °C, it is stable until expiration date. Dispose any remaining diluted buffer.

### **Preparation of samples**

- This kit is intended to measure IFN- $\gamma$  in mouse serum or plasma.
- We recommend to use EDTA-2Na, EDTA-2K or heparin as anticoagulant.
- Samples should be immediately assayed or stored below -35 °C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using buffer solution prior to adding them to wells

### **Assay procedure**

Remove the cover sheet of the antibody-coated plate after bringing up to room temperature.


- (1) Wash the anti-IFN- $\gamma$  antibody coated plate (A) by filling the wells with washing buffer and discard 4 times (\*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells. (\*⑤)
- (2) Pipette 50  $\mu$ L of standard solution to the wells designated for standards.
- (3) Pipette 25  $\mu$ L of (C) Buffer solution to the sample wells, and the add 25  $\mu$ L of each sample to the wells.
- (4) Shake the plate gently on a plate shaker (\*②).
- (5) Stick a plate seal (\*③) on the plate and incubate for 2 hours at 20 °C - 25 °C.
- (6) Discard the reaction mixture and rinse wells as step (1).

**LBIS Mouse IFN- $\gamma$  ELISA Kit**

- (7) Pipette 50  $\mu$ L of biotinylated anti-IFN- $\gamma$  antibody to all wells, and shake as step (4).
  - (8) Stick a plate seal (\*③) on the plate and incubate the plate for 1 hour at 20 °C - 25 °C.
  - (9) Discard the reaction mixture and rinse wells as step (1).
  - (10) Pipette 50  $\mu$ L of HRP-conjugated streptavidin to all wells, and shake as step (4).
  - (11) Stick a plate seal (\*③) on the plate and incubate the plate for 30 minutes at 20 °C - 25 °C.
  - (12) Discard the reaction mixture and rinse wells as step (1).
  - (13) Pipette 50  $\mu$ L of Chromogen (TMB) to wells, and shake as step (4).
  - (14) Stick a plate seal (\*③) on the plate and incubate the plate for 20 minutes at 20 °C - 25 °C.
  - (15) Add 50  $\mu$ L of the stop solution to all wells and shake as step (4).
  - (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (\*④)) immediately using a plate reader.
- \*Refer to the page 7-8 for notes of \*①, \*②, \*③, \*④, and \*⑤.

**Calculations**

- (1) Create a standard curve by reducing the data using computer software capable of generating 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse IFN- $\gamma$  concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the IFN- $\gamma$  concentration of samples at its absorbance, and multiply the assay value by dilution factor if the sample has been diluted. The standard operation method is 2-fold dilution.

**Calibration**

The measured value obtained by this kit can be converted to the unit concentration based on NIH standard IFN- $\gamma$  (Gg02-901-533) by multiplying by the conversion factor specified in a separate sheet according to the following calculation formula.

- NIH (Gg02-901-533) unit concentration (U/mL) = conversion factor (\*1)  $\times$  measured IFN- $\gamma$  value (pg/mL) (\*2)
- (\*1) The exchange rate is different depending on the lot. Please see the specified conversion factor in the annex document.
  - (\*2) Use the assay value that multiplied by dilution rate.

●Precision of assay (Within assay variation)  
2 serum samples, 5 replicates assay

n / ID	Sample1 (pg/mL)	Sample2 (pg/mL)

●Reproducibility (Between assay variation)  
3 serum samples, 4 days, 3 replicates assay

Day /ID	Sample1 (pg/mL)	Sample2 (pg/mL)	Sample3 (pg/mL)

## LBIS Mouse IFN- $\gamma$ ELISA Kit

1	288	31.6
2	281	30.8
3	276	30.2
4	281	31.2
5	275	29.6
mean	280	30.7
SD	5.2	0.79
CV(%)	1.9	2.6

Mean CV was within 15 %.

0	200	32.0	5.22
1	193	32.2	5.16
2	192	33.0	4.91
3	195	32.7	4.89
mean	195	32.5	5.04
SD	3.5	0.45	0.17
CV(%)	1.8	1.4	3.3

Mean CV was within 15 %.

- Recovery test (Standard Mouse IFN- $\gamma$  was added in 3 concentrations to serum or plasma sample and assayed.) The recoveries were 96.0 % – 109 %


Serum sample

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
–	8.70	–	–
20.0	27.9	19.2	96.0
80.0	95.5	86.8	109
165	174	165	100

Plasma sample (EDTA)

Added (pg/mL)	Found (pg/mL)	Recovered (pg/mL)	Recovery (%)
–	4.14	–	–
20.0	24.3	20.2	101
80.0	84.3	80.2	100
165	173	169	102

- Dilution test (Serum sample and plasma sample were serially diluted by 3 steps.) The dilution curves showed linearity with  $R^2=0.996$  –  $0.998$

### Trouble shooting

- Low absorbance in all wells

Possible explanations:

- The standard or samples might not be added.
- Reagents necessary for coloration such as Biotinylated anti-IFN- $\gamma$  antibody, HRP-conjugated streptavidin, or chromogen (TMB) might not be added.
- Wrong reagents related to coloration might have been added. Wrong dilution of biotinylated anti-IFN- $\gamma$  antibody or HRP-conjugated streptavidin.
- Contamination of enzyme inhibitor(s).
- Influence of the temperature under which the kits had been stored.

## LBIS Mouse IFN- $\gamma$ ELISA Kit

6) Excessive hard washing of the well plate.

7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD was higher than that of the lowest standard concentration (2.05 pg/mL).

Possible explanations:

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with HRP-conjugated streptavidin.)

- High coefficient of variation (CV)

Possible explanation:

1) Improper or inadequate washing.

2) Improper mixing of standard or samples.

3) Pipetting at irregular intervals.

- Q-1: Can I divide the plate to use it for the other testing?

A-1: Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2: I found 96 well-plate is empty when I opened the box.

A-2: As this kit is dried type, not preservation stabilizer is added.

### Technical tips / Precautions

- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20 °C to 25 °C. Avoid performing assay under the airstream velocity over 0.4 m/sec. and the humidity less than 30 %, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the antibody-coated plate is module type of 8 wells  $\times$  12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The chromogen (TMB) should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of stop solution. Greenish color means incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- The materials must not be pipetted by mouth.
- Do not use reagents with different lot numbers together.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Residual samples and used tips should be sterilized before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

## LBIS Mouse IFN- $\gamma$ ELISA Kit

### Summary of assay procedure : Use as a check box

\*First, read this instruction manual carefully and start your assay after confirmation of details.

- Bring the well-plate and all reagents back to **20 °C - 25 °C for 2 hours.**
- Concentrated washing buffer must be diluted to **10 times** by deionized water (or distilled water) that returned to 20 °C - 25 °C.
- Standard solution dilution example  
Reconstitute the (B) standard with volume\* of deionized water (or distilled water) described in a separate sheet to prepare the Original standard solution (5 ng/mL) . And then prepare dilution series of standard solutions as shown below using the (C) Buffer solution brought to the room temperature.  
\*Because the volume of deionized water (or distilled water) to be added to the freeze-dried standard differs depending on the lot, see the specified volume in a separate sheet.

Preparation of diluted IFN- $\gamma$  standard solutions:

Conc.(pg/mL)	500	200	80.0	32.0	12.8	5.12	2.05	0
Std. sol.( $\mu$ L)	50	100*	100*	100*	100*	100*	100	—
Buffer( $\mu$ L)	450	150	150	150	150	150	150	200

\*One rank higher standard.

<input type="checkbox"/> Anti-IFN- $\gamma$ antibody coated plate	
<input type="checkbox"/> ↓Washing 4 times (*①) , (*⑤)	
<input type="checkbox"/> (Diluted) samples, or Standards	50 $\mu$ L
<input type="checkbox"/> ↓Shaking (*②), Incubation for 2 hours at 20 °C - 25 °C. (Standing (*③)) Add deionized water (or distilled water) of room temperature to (D) Biotinylated anti-IFN- $\gamma$ antibody (Freeze-dried) ( Original (D) solution). Prepare working solution by dilution of original (D) solution to <b>1:100</b> with the (C) Buffer solution returned to 20 °C - 25 °C. (This should be prepared during incubation.)	
<input type="checkbox"/> ↓Washing 4 times (*①) , (*⑤)	
<input type="checkbox"/> Biotinylated anti-IFN- $\gamma$ antibody	50 $\mu$ L
<input type="checkbox"/> ↓Shaking (*②), Incubation for 1 hour at 20 °C - 25 °C. (Standing (*③)) Dilute HRP-conjugated streptavidin to <b>1:100</b> with the (C) Buffer solution returned to 20 °C - 25 °C. (This should be prepared during incubation.)	
<input type="checkbox"/> ↓Washing 4 times(*①) , (*⑤)	
<input type="checkbox"/> HRP-conjugated streptavidin	50 $\mu$ L
<input type="checkbox"/> ↓Shaking (*②), Incubation for 30 minutes at 20 °C - 25 °C. (Standing (*③))	
<input type="checkbox"/> ↓Washing 4 times(*①) , (*⑤)	
<input type="checkbox"/> Chromogen (TMB) After dispense, the color turns to blue depending on the concentration.	50 $\mu$ L
<input type="checkbox"/> ↓Shaking (*②), Incubation for 20 minutes at 20 °C - 25 °C. (Standing (*③))	
<input type="checkbox"/> Stop solution After dispense, the color turns to yellow depending on the concentration.	50 $\mu$ L
<input type="checkbox"/> ↓Shaking (*②) Immediately shake.	
<input type="checkbox"/> Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately. Ref. wave cancels the dirt in the back of plate.	

\*①After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec and remove the buffer.  
Guideline of washing volume: 300  $\mu$ L/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to

## LBIS Mouse IFN- $\gamma$ ELISA Kit

be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with HRP-conjugated streptavidin.

Standard of plate-washing pressure: 5 - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


\*② Guideline of shaking: 600 rpm - 1,200 rpm for 10 seconds  $\times$  3 times.

\*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.

\*④ 600 nm - 650 nm can be used as reference wavelength.

\*⑤ After removal of wash buffer, immediately dispense the next reagent.

### Assay worksheet

## LBIS Mouse IFN- $\gamma$ ELISA Kit

[Storage condition] Store the kit at 2 °C - 8 °C (Do not freeze)

[Term of validity] Expiration date is indicated on the container

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