**FUJIFILM****Wako**

Code No. 292-83801 (96 reactions)

**CD81-Capture Human Exosome  
ELISA Kit (Streptavidin HRP)****[Materials Supplied]**

Components	Amount	State	Preparation before use
(A) Anti-CD81 Antibody-immobilized 96 Well Plate	8 well × 12 strips / 1 plate	Use after washing	Wash 3 times with Washing Buffer (1×).
(B) Plate Seal	4 sheets	—	—
(C) Sample Reaction Buffer	50 mL × 1 vial	Ready-to-Use	—
(D) Antibody Reaction Buffer	50 mL × 1 vial	Ready-to-Use	—
(E) Washing Buffer (10×)	100 mL × 1 vial	Concentrated. Use after dilution.	Bring it to room temperature before opening and dispense the required volume. Dilute 10 times with purified water (distilled water). [Washing Buffer (1×)]
(F) Control Biotinylated Antibody Anti-CD81 (100×)	120 μL × 1 tube	Concentrated. Use after dilution.	Dilute 100 times with Antibody Reaction Buffer. [Biotinylated antibody]
(G) HRP-conjugated Streptavidin (100×)	240 μL × 1 tube	Concentrated. Use after dilution.	Dilute 100 times with Antibody Reaction Buffer. [HRP-conjugated Streptavidin (1×)]
(H) TMB Solution	12 mL × 1 vial	Ready-to-Use	—
(I) Stop Solution	12 mL × 1 vial	Ready-to-Use	—
(J) Instruction Manual	1 copy	—	—

**[Storage]**

Store at 2-10°C.

**[Materials Required, Not Supplied]**

Purified water (distilled water)
Biotinylated antibody against any extracellular vesicle surface marker protein (for detection except human CD81) If a biotinylated antibody is unavailable, antibodies can be labeled with biotin using Biotin Labeling Kit-SH (Code No. 348-90941) or Biotin Labeling Kit-NH <sub>2</sub> (Code No. 347-90891).
Reference Standard (for quantitative analysis) Prepare purified extracellular vesicles from appropriate culture supernatants of cell strain or body fluids using <b>MagCapture™ Exosome Isolation Kit PS</b> (Code No. 293-77601) or <b>PS Capture™ Exosome Isolation Resin Kit</b> (Code No. 290-80301). Measure protein concentration and particle number. Make sure that the purified extracellular vesicles have any surface marker protein and use it as the reference standard.

— 1/8 —

**[Other supplies Required]**

Test tubes for dilution of samples and antibodies
Glassware for preparation of Reaction/Washing Buffer (Graduated cylinder)
Micropipettes
Syringe-type repeating dispenser (If available) Such as Eppendorf Multipipette or ThermoFisher Finnpiipette
Disposable reservoir (As necessary)
Paper towels to remove the remained solution in the wells after washing steps
Vortex-type mixer
Shaker for microplate (500 rpm) Recommended products : Code No. 623-05671 MS 3 digital, Code No. 627-05691 MS 3 basic
Automatic washer for 96-well plate (If available) Recommended products : Code No. 510-22411 HydroFlex M8Ch2
Microplate reader capable of measuring at 450 nm, with the correction wavelength set at 600-650 nm Recommended products : Code No. 518-84231 InfiniteF50R
Software for data analysis Recommended products : Code No. 290-34631 PLATEmanager V5/I PC SET (for infinite)

**[Standard Curve Preparation for Quantitative Measurement]**  
(for quantitative analysis)

Dilute the reference standard of extracellular vesicles to an appropriate concentration with Sample Reaction Buffer in order to make its absorbance value at 450 nm (absorbance at 620 nm is subtracted) for ELISA measurement to be within the range of 0.1 to 3.0.

**[Sample Dilution]**

Dilute the sample appropriately with Sample Reaction Buffer using a tube, etc. so that the absorbance at 450 nm (value obtained by subtracting the absorbance at a secondary wavelength of 620 nm) measured by this ELISA is in the range of 0.2 to 2.5<sup>\*1, 2</sup>.

\*1 Dilute cell culture supernatant, serum, and heparin plasma samples twice or more. It is also recommended to dilute other samples at least 2 times.

\*2 Since the appropriate dilution rate of the sample varies greatly depending on the concentration of extracellular vesicles in the sample and the amount of surface protein marker to be detected, it is recommended to serially dilute the sample and preliminarily examine the appropriate dilution rate.

**[Assay Procedure]**

- 1) Wash each well of Anti-CD81 Antibody-immobilized 96 Well Plate (only strips to use) 3 times with 300 to 350 μL of Washing Buffer (1×)<sup>\*3</sup>.
- 2) Add 100 μL of sample dilutions, standard dilutions (when quantitative analysis is performed), and Sample Reaction Buffer as a blank to each well.
- 3) Cover with Plate Seal<sup>\*4</sup>. Incubate for 2 hours at room temperature (20 to 25°C) with shaking at about 500 rpm using a microplate shaker<sup>\*5</sup>.

— 2/8 —

- 4) Discard the solution and wash 3 times with 300 to 350  $\mu\text{L}$  of Washing Buffer (1 $\times$ )<sup>\*\*3</sup>.
- 5) Add 100  $\mu\text{L}$  of biotinylated antibody to each wells.
- 6) Cover with new Plate Seal<sup>\*\*4</sup>. Incubate for 1 hour at room temperature (20 to 25°C) with shaking at about 500 rpm using a microplate shaker<sup>\*\*5</sup>.
- 7) Discard the solution and wash 3 times with 300 to 350  $\mu\text{L}$  of Washing Buffer (1 $\times$ )<sup>\*\*3</sup>.
- 8) Add 100  $\mu\text{L}$  of HRP-conjugated Streptavidin (1 $\times$ ) to each wells.
- 9) Cover with new Plate Seal<sup>\*\*4</sup>. Incubate for 2 hours at room temperature (20 to 25°C) with shaking at about 500 rpm using a microplate shaker<sup>\*\*5</sup>.
- 10) Discard the solution and wash 5 times with 300 to 350  $\mu\text{L}$  of Washing Buffer (1 $\times$ )<sup>\*\*3</sup>.
- 11) Add 100  $\mu\text{L}$  of TMB Solution brought to room temperature to each wells. Shake for about 1 minute using microplate shaker.
- 12) Cover with new Plate Seal<sup>\*\*4</sup>. Incubate for 30 minutes at room temperature (20 to 25°C).
- 13) Add 100  $\mu\text{L}$  of Stop Solution brought to room temperature to each wells.
- 14) Shake for about 5 seconds by microplate shaker. Measure immediately the absorbance at 450 nm and 620 nm as reference wavelength (600-650 nm)<sup>\*\*6</sup> using a microplate reader.

<sup>\*\*3</sup> After the last wash, turn the plate upside down on a stack of paper towels and tap it lightly to remove the liquid remaining in the wells.

<sup>\*\*4</sup> When plate strips are separately used, cut a Plate Seal into a size of the strips for use.

<sup>\*\*5</sup> Incubation with shaking at about 500 rpm using a microplate shaker is recommended since standing condition may cause low detection sensitivity and large well to well variation.

<sup>\*\*6</sup> Reference wavelength can be used in the range of 600-650 nm.

#### **[Calculation of Results]**

In the following calculation, use the calculated value at 450 nm obtained by subtracting the absorbance readings at 620 nm (600-650 nm).

##### **<Qualitative measurement>**

Calculate the value by subtracting the absorbance readings of blank at 450 nm from the absorbance readings of diluted sample at 450 nm, and compare it between samples.

##### **<Quantitative measurement>**

- 1) Calculate the value by subtracting the absorbance readings of blank at 450 nm from the absorbance readings of diluted samples at 450 nm<sup>\*\*7</sup>.
- 2) Calculate the value by subtracting the absorbance readings of blank at 450 nm from the absorbance readings of dilution series of the reference standard at 450 nm.
- 3) Prepare a standard curve by plotting the absorbance readings (Y-axis) against values of standard concentration (X-axis)<sup>\*\*8</sup>.
- 4) Determine the unknown sample concentration from the Standard Curve<sup>\*\*9</sup> and multiply the value by the dilution factor.

<sup>\*\*7</sup> When the absorbance value of the diluted sample falls outside the absorbance range of standard curve, dilute it again to the appropriate concentration and perform the measurement.

<sup>\*\*8</sup> Prepare a new standard curve for each assay performed.

<sup>\*\*9</sup> The use of a 3rd order regression curve for log-log plot or 4 or 5 parameters method for log-normal plot in computer calculation is recommended.

#### **[Precaution]**

- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Handle any reagents carefully for preventing any transmissions during laboratory procedures. This kit contains reagents of animal origin.
- Used samples and tips should be rinsed in 1% formalin, 2% glutaraldehyde or more than 0.1% sodium hypochlorite solution for more than 1 hour or be treated by autoclaving before disposal. Dispose of waste in accordance to applicable national, regional or local regulations.
- When allowing the plate to incubate in each step, always affix a plate seal to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of the dispensed reagents.
- ELISA can be easily affected by the laboratory environment. Follow the assay procedure and the temperature in the laboratory should be strictly maintained at 20-25°C. Also, avoid measurement in an environment with wind current (including air conditioner wind) and low humidity.
- Avoid cross contamination of samples or reagents by changing tips between sample, standard and reagent additions.

#### **[Additional Information]**

Detailed information about this kit is available on our website. Please review this information of a better understanding of product features, performance, and operation.

Please access from the URL below and search by Product Code.

<https://labchem-wako.fujifilm.com>

## **FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### **FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

### **FUJIFILM Wako Chemicals Europe GmbH**

Fuggerstrasse 12  
 D-41468 Neuss  
 Germany  
 Telephone : + 49-2131-311-0  
 Facsimile : + 49-2131-311-100  
<http://www.wako-chemicals.de>

Code No. 292-83801 (96 reactions)

## CD81-Capture Human Exosome ELISA Kit (Streptavidin HRP)

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